

Product datasheet for **RC239029**

C13orf18 (RUBCNL) (NM_001286763) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C13orf18 (RUBCNL) (NM_001286763) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C13orf18
Synonyms:	C13orf18; KIAA0226L; PACER
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC239029 representing NM_001286763
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**

ATGGTTTCCAATCATTACTTTCTCTGTGTCAATCTGCCTCTTCGAGAAATTCATACTCTGGCCCTT
 CACCTTCGTGCCGCGGACTCCCTGGCAGAGACAACGTTGTCTGAGGATACCACAGACTCCGTTGGCAG
 CGCTTCTCCCATGGCTCGAGTGAAGAGTAGCAGCTTCTCTGTCTCAACAGAGGTACACATGGTC
 CGCCAGGATACTCTCATCGGTGTCTCTGCCACAAGCCCTGGGATTTTGGCCACCTCCCATATCCTG
 AGACTGACAGTGCTTTTTTGGCCTTCCCATCTGACATCTGCTGCTGATGAAGGTGCTGTTCAAGTCAG
 TAGAAGAACCATTCTTGAATTCCTTCTACCAGAGGTATTTGTGCTGCCTGTTGATGTAGAAAAGGAA
 AATGCCACTTTTATGTTGCAGATATGATTATATCAGCAATGGAGAAAATGAAGTGAACATTCTGAGTC
 AACAGCAGACAGAGAGCTGGAGTAAAGAAGTCAGTGGGTTACTTGGGAGTGATCAGCCTGACTCTGAAAT
 GACTTTTGATACCAACATAAAGCAAGAGTCTGGGTCTTCTACTTCTCATACAGTGGCTATGAAGTTGT
 GCTGTGTTACAGGTCAGCCAGTGAACACGTACTTACCATGATGTGAAAGAGATTGCAAAATGCG
 ATGTTGATGAATTTGTTATTTAGAGCTTGGAGATTTTATGATATCACAGAAACCTGTAGCTGTTCTCTG
 CAGCTCCTCTAAGAGTGTCACTTATGAGCCAGACTTCAATTCTGCAGAACTATTAGCCAAAGAGCTGTAC
 CGCGTGTTCAGAAAGTGTGGTACTGTGCTAGTAAATTCTCAGCTGGCAGGTTCCCTGAGTGCAGCTG
 GCTCGATAGTCGTAATGAAGAGTGTGTCCGAAAAGACTTTGAATCCAGTATGAATGTAGTACAGGAAAT
 TAAATTTAAGTCTAGGATCAGAGGGACTGAAGACTGGGCTCCTCTAGATTTCAAATCATATTTAATATT
 CATCCACCACTCAAGAGGGACCTTGTGGTGGCAGCCAGAATTTTCTGTGCCGCTGTGGAACCTCAG
 TAGAGCTAAGTTTGTGAAGCGCTCCGGTACTGCGAATACCTAGGGAAGTATTTCTGTGACTGCTGCCA
 CTCATATGCAGAGTCGTGCATCCCTGCCGAATCCTGATGATGTGGGACTTCAAGAAGTACTACGTGAGC
 AATTTCTCCAACAGCTGCTCGACAGCATATGGCACCAGCCATTTTCAATTTGCTGAGCATCGGCCAAA
 GCCTGTATGCGAAAGCCAAGGAGCTGGACAGAGTGAAGGAAATTCAGGAGCAGCTCTTCCATATCAAGAA
 GCTGTTGAAGACCTGTAGGTTTGTAAACAGTGCATTAAGGAGTTCGAGCAGGTGCCGGGACACTTGACT
 GATGAGCTCCACCTGTTCTCCCTTGAGGACCTGGTCAGGATCAAGAAAGGGCTGCTGGCACCCCTTACTCA
 AGGACATTCTGAAAGCTTCCCTGCACATGTGGCTGGCTGTGAGCTGTGTCAAGGAAAGGGCTTTATTTG
 TGAATTTTGCAGAATACGACTGTCATCTTCCATTTAGACAGCAACATGTAGAAGATGTTGAGCGTGC
 AGGGCTTGCTTTCAAAACAGTCTTCCAGTCTCCGAGTGCCCCGGTGTGCGAGGATCACAGCGAGGA
 GAAACTTCTGAAAGTGTGGCCTCTGCAGCAACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC239029 representing NM_001286763
 Red=Cloning site Green=Tags(s)

MVSNHYFLLCVNLPLREIHTPGSPSPCLGDSLAEITLSEDTTDSVGSASPHGSSEKSSSFSLSTEVMHV
 RPGYSHRVSLPTSPGILATSPYPETDSAFFEPShLTSAADEGAVQVSRRTISSNSFSPEVFLPVDVEKE
 NAHFYVADMIISAMEKMKCNILSQQTESWSKEVSGLLGSDQPDSEMTFDTNIKQESGSTSSYSGYEGC
 AVLQVSPVTETRTYHDVKEICKCDVDEFVILELGDNDITETCSCSCSSKSVTYEPDFNSAELLAKELY
 RVFQKCWILSVVNSQLAGSLAAGSIVVNEECVRKDFESSMNVVQEIKFSRIRGTEDWAPPRFQIIFNI
 HPPLKRDLVVAAQNFFCAGCGTPVEPKFVKRLRYCEYLKGYFDCCHSYAESCIIPARILMMWDFKKYYVS
 NFSKQLLDSIWHQPIFNLLSIGQSLYAKAKELDRVKEIQEQLFHIKLLKTCRFANSALKEFEQVPGHLT
 DELHLFSLLEDLVRIKKGLLAPLLKDILKASLAHVAGCELCQKGKGFICEFCQNTTVIFPFQTATCRRCSAC
 RACFHKQCFQSSECPRCARITARRKLLESVASAAT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001286763

ORF Size: 1785 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001286763.3](#)

RefSeq Size: 3592 bp

RefSeq ORF: 1788 bp

Locus ID: 80183


UniProt ID: [Q9H714](#)

Cytogenetics: 13q14.13

MW: 66.8 kDa

Gene Summary: This gene encodes a cysteine-rich protein that contains a putative zinc-RING and/or ribbon domain. The encoded protein is related to Run domain Beclin-1-interacting and cysteine-rich domain-containing protein, which plays a role in endocytic trafficking and autophagy. In cervical cancer cell lines, this gene is expressed at low levels and may function as a tumor suppressor. Promoter hypermethylation of this gene is observed in cervical cancer cell lines and tissue derived from human patients. [provided by RefSeq, Mar 2017]

Product images:



Circular map for RC239029