

Product datasheet for **RC239061**

DNA helicase HEL308 (HELQ) (NM_001297756) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNA helicase HEL308 (HELQ) (NM_001297756) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DNA helicase HEL308
Synonyms:	HEL308
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC239061 representing NM_001297756
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGTTGAAGACCTACAAAAGTTTCTTCAAGCAGAATATTATACCACTCAATTTAGACCACTATGTACCTA
 AGGTTGAGTTAAAAGAATATCTGAAAATAAATGATACAATATATGAAGTTGACAGCAAAGCTGAGAATGG
 CATGACTTTTTACGTCCTTCTTAATTATAAGTATTCTGATACCCTAAAAAAGATGGATCCTGATCACTTG
 GTAGCATTGGTGACAGAAGTTATCCCAATTATTCCTGCTTAGTTTTTTGTCTAGTAAGAAGAACTGTG
 AAAATGTAGCAGAAATGATATGCAAAATTTTAAGCAAGGAATATCTGAAACATAAGGAGAAAAGAAAATG
 TGAGGTGATTAAGAAGTTGAAGAATATTGGCAATGGCAACCTGTGCTGTTTTAAAGCGCACTATCCCA
 TTTGGAGTTGCCTATCACCACAGTGGCTTAACAAGTGATGAAAGGAACTCTGGAGGAGGCTACTCCA
 CAGGAGTGCTCTGTCTTTTACCTGCACATCTACCCTAGCGGCAGGTGTCAACCTACCAGCTCGAAGAGT
 TATTTTAAGAGCTCCCTATGTTGCTAAGGAATTTTAAAGAGGAATCAATATAACAGATGATTGCGAGA
 GCTGGTCGTGCTGGAATAGATACTATTGGGGAGAGTATCCTCATATTGCAAGAAAAAGACAAACAACAGG
 TATTGGAGTTAATAACTAAACCATTGAAAACTGTTACAGCCATCTTGTTTCAGGAATTCACCAAGGGAAT
 CCAACATTATTTCTCTCTTTGATTGGTTTGAAGATTGCAACGAATCTTGATGACATCTATCATTTTCATG
 AATGGTACATTTTTTGGTGTTCAGCAAAAGGTTTTATTGAAAGAAAAAGTCTCTGGGAAATAACTGTTG
 AATCACTTAGATACCTGACAGAAAAAGGACTCTACAAAAAGACACTATTTATAAGTCTGAAGAAGAGGT
 CCAATATAATTTTCATATTACAAAGTTGGGACGTGCTTCATTTAAGGGAATATAGATTTAGCTTATTGT
 GACATTCTGTACAGAGACTTGAAGAAAGGTCTTGAAGGACTTGTGCTTGAAGCCTTCTTCATCTAATCT
 ACCTAACAAACCCCTATGATCTGGTTTACAGTGTAACCTGATTGGATGATATACTTCAGGCAGTTTAG
 CCAACTCAGTCCAGCAGAACAAAATGTAGCTGCCATTCTTGGAGTCTCTGAAAGCTTTATTGGGAAGAAA
 GCATCAGGCCAAGCCATCGGAAAGAAGGTGGACAAGAACGTTGTCAACAGGCTATATCTGTCTTTTGTTT
 TTTATACCTTGCTCAAAGAGACCAACATTTGGACTGTATCTGAAAAATTTAATATGCCTCGAGGATATAT
 ACAAAATCTTCTCACTGGAAGTGCCTATTCTCATCTTGTGTGTACATTTCTGTGAGGAGCTTGAGGAG
 TTTTGGGTTTACAGAGCCCTTTTGGTAGAAGTTACCAAGAAGCTGACTTACTGTGTAAAGGCAGAAATTA
 TCCCTCTCATGGAAGTTACTGGAGTTTATAGAGGTCGAGCAAAACAGTTATACAGTGCAGGTTACAAAAG
 TCTAATGCACTTAGCTAATGCAATCCTGAAGTCTCGTAAGGACAATTGATCATTTATCAAGACGCCAA
 GCCAAGCAAATTGTTTCATCAGCAAAGATGCTGTTGCATGAAAAAGCAGAAGCCCTGCAAGAAGAGGTAG
 AAGAGTTACTAAGATTGCCTTCTGATTTCCCTGGTGTGTGGCTTCTTCCACTGACAAAGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC239061 representing NM_001297756
 Red=Cloning site Green=Tags(s)

MLKTYKSFFKQNIIPVNLQYVPKVELKEYLKINDTIYEVDSKAENGMTFSRLLNKYSDTLKKMDPDHL
 VALVTEVIPNYSCLVFCPSKNCENVAEMICKFLSKEYLKHKEKEKCEVIKNLKNIGNLCPVLKRTIP
 FGVAHHHSLTSDERKLEEAYSTGVLCFTCTSTLAAGVNLPRRVILRAPYVAKEFLKRNQYQKMIGR
 AGRAGIDTIGESILILQEKDQQLLELITKPLENCYSHLVQEFTKGIQTLFLSLIGLKIATNLDDIYHFM
 NGTFFGVQKQVLLKEKSLWEITVESLRYLTKGLLQKDTIYKSEEEVQYNFHITKLGRASFKGTIDLAYC
 DILYRDLKKGLEGLVLESLLHLIYLTPYDLVSQCNPDWMIYFRQFSQSPAEQNVAAILGVSEFSGKK
 ASGQAIGKKVDKNVVRNLYLSFVLYTLKETNIWTVSEKFNMPRGYIQNLLTGTASFSSCVLHFCEEELE
 FWYRALLVELTKKLYCVKAELIPLMEVTGVLEGRAKQLYSAGYKSLMHLANANPEVLVRTIDHLSRRQ
 AKQIVSSAKMLLHEKAEALQEEVEELLRLPSDFPGAVASSTDKA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001297756

ORF Size: 1812 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001297756.1](#), [NP_001284685.1](#)

RefSeq Size: 3636 bp

RefSeq ORF: 1815 bp


Locus ID: 113510

Cytogenetics: 4q21.23

MW: 69.2 kDa

Gene Summary: HEL308 is a single-stranded DNA-dependent ATPase and DNA helicase (Marini and Wood, 2002 [PubMed 11751861]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RC239061