

Product datasheet for **RC239082**

MCAK (KIF2C) (NM_001297657) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MCAK (KIF2C) (NM_001297657) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KIF2C
Synonyms:	CT139; KNSL6; MCAK
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239082 representing NM_001297657
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCCACTGTCTCAGAGCTTCGCATCACGGCTCAGGAGAATGACATGGAGGTGGAGCTGCCTGCAGCTG
 CAAACTCCCGCAAGCAGTTTTAGTTCCCTCCTGCCCCACTAGGCCTTCTGCCCTGCAGTGGCTGAAAT
 ACCATTGAGGATGGTCAGCGAGGAGATGGAAGAGCAAGTCCATTCATCCGAGGCAGCTCTTCTGCAAAC
 CCTGTGAACTCAGTTCGGAGGAAATCATGTCTTGTGAAGGAAGTGGAAAAATGAAGAACAAGCGAGAAG
 AGAAGAAGGCCAGAACTCTGAAATGAGAATGAAGAGAGCTCAGGAGTATGACAGTAGTTTTCAAACCTG
 GGAATTTGCCGAATGATTAAGAATTTCCGGCTACTTTGGAATGTCATCCACTACTATGACTGATCCT
 ATCGAAGAGCACAGAATATGTGTCTGTGTAGGAAACGCCCACTGAATAAGCAAGAATTGGCCAAGAAAG
 AATTGATGTGATTTCCATTCTAGCAAGTGTCTCCTCTTGGTACATGAACCAAGTTGAAAGTGGACTT
 AACAAAGTATCTGGAGAACCAAGCATTCTGCTTTGACTTTGCATTTGATGAAACAGCTTCGAATGAAGTT
 GTCTACAGGTTACAGCAAGGCCACTGGTACAGACAATCTTTGAAGGTGGAAAAGCAACTGTTTTGCAT
 ATGGCCAGACAGGAAGTGGCAAGACACATACTATGGGCGGAGACCTCTCTGGGAAAGCCAGAAATGCATC
 CAAAGGGATCTATGCCATGGCCTCCCGGAGCTTCTCCTCTGAAGAATCAACCTGTACCGGAAGTTG
 GGCCTGGAAAGTCTATGTGACATCTTCGAGATCTACAATGGGAAGCTGTTTGACCTGCTCAACAAGAAAG
 CCAAGCTGCGCGTGTGGAGGACGGCAAGCAACAGGTGCAAGTGGTGGGGCTGCAGGAGCATCTGGTTAA
 CTCTGCTGATGATGTCATCAAGATGATCGACATGGGCAGCGCTGCAGAACCTCTGGGCAGACATTTGCC
 AACTCCAATTCCTCCCGCTCCACGCGTGTCTTCAAATATTCTTCGAGCTAAAGGGAGAATGCATGGCA
 AGTTCTCTTTGGTAGATCTGGCAGGGAATGAGCGAGGCGCGACACTTCCAGTGTGACCGGCAGACCCG
 CATGGAGGGCGCAGAAAATCAACAAGAGTCTCTTAGCCCTGAAGGAGTGCATCAGGGCCCTGGGACAGAAC
 AAGGCTCACACCCGTTCCGTGAGAGCAAGCTGACACAGGTGCTGAGGGACTCCTTCATTGGGGAGAACT
 CTAGGACTTGATGATTGCCAGATCTCACCAGGCATAAGCTCCTGTGAATATACTTTAAACACCCTGAG
 ATATGCAGACAGGGTCAAGGAGCTGAGCCCCACAGTGGGCCAGTGGAGAGCAGTTGATTCAAATGGAA
 ACAGAAGAGATGGAAGCCTGCTTAACGGGGCGCTGATTCCAGGCAATTTATCCAAGGAAGAGGAGGAAC
 TGTCTTCCAGATGTCCAGCTTAAACGAAGCCATGACTCAGATCAGGGAGCTGGAGGAGAAGGCTATGGA
 AGAGCTCAAGGAGATCATACAGCAAGGACCAGACTGGCTTGGCTCTCTGAGATGACCGAGCAGCCAGAC
 TATGACCTGGAGACCTTTGTGAACAAAGCGGAATCTGCTCTGGCCAGCAAGCCAAGCATTCTCAGCCC
 TCGGAGATGTCATCAAGCCCTTGCCTGGCCATGCAGCTGGAAGAGCAGGCTAGCAGACAAATAAGCAG
 CAAGAAACGGCCCCAG

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC239082 representing NM_001297657
 Red=Cloning site Green=Tags(s)

MSTVSELRITAQENDMEVELPAAANSRKQFSVPPAPTRPSCPAVAEIPLRMVSEEMEEQVHSIRGSSSAN
 PVNSVRRKSCLVKEVEKMKNKREEKKAQNSEMRMKRAQEYDSSFNWEFARMIKEFRATLECHPLTMTDP
 IEEHRICVCRKRPLNKQELAKKEIDVISIPSKCLLLVHEPKLVDLTKYLENQAFCFDFAFDEASNEV
 VYRFTARPLVQTFEGGKATCFAYGQTGSGKTHTMGGDLGSKAQNASKGIYAMASRDVFLLNQPCYRKL
 GLEVYVTFEYIYNGKLFDLLNKKAKLRVLEDGKQVQVVLQEHVNSADDVIKIDMGSA CRTSGQTF
 NSNSSRSHACFQIILRAKGRMHGKFSVLDLAGNERGADTSSADRQTRMEGAEINKSLLALKECIRALGQN
 KAHTPFRESKLTQVLRDSFIGENSRTCMIAISPGISSEYTLNTRYADRVKELSPHSGPSGEQLIQME
 TEEMEACSNALIPGNLSKEEELSSQMSSFNEAMTQIRELEEKAMEELKEIIQQGPDWLELSEMTEQPD
 YDLETFVNKAESALAQAKHFSALRDVIKALRLAMQLEEQASRQISSKKRPQ


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297657

ORF Size: 1836 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001297657.2
RefSeq Size:	2701 bp
RefSeq ORF:	1839 bp
Locus ID:	11004
UniProt ID:	Q99661
Cytogenetics:	1p34.1
Protein Families:	Druggable Genome
MW:	69.3 kDa
Gene Summary:	This gene encodes a kinesin-like protein that functions as a microtubule-dependent molecular motor. The encoded protein can depolymerize microtubules at the plus end, thereby promoting mitotic chromosome segregation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]