

Product datasheet for **RC239084**

CLASRP (NM_001278439) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLASRP (NM_001278439) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CLASRP
Synonyms:	CLASP; SFRS16; SWAP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

```
>RC239084 representing NM_001278439
Red=Cloning site Blue=ORF Green=Tags(s)
```

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

[illegible]

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC239084 representing NM_001278439
Red=Cloning site Green=Tags(s)

MWHEARKHERKL R G M M V D Y K K R A E R R E Y Y E K I K D P A Q F L Q E S D E R K C N Y E R Y R G L V Q N D F A G I S E E Q C
 L Y Q I Y I D E L Y G G L Q R P S E D E K K L A E K K A S I G Y T Y E D S T V A E V E K A A E K P E E E E S A A E E E S N S D E D E V I P
 D I D V E V D V D E L N Q E Q V A D L N K Q A T T Y G M A D G D F V R M L R K D K E E A E A I K H A K A L E E E K A M Y S G R R S R R Q R R
 E F R E K R L R G R K I S P P S Y A R R D S P T Y D P Y K R S P S E S S S E S R S R S R S P T P G R E E K I T F I T S F G G S D E E A A A A
 A A A A A S G V T T G K P P A P P Q P G G P A P G R N A S A R R R S S S S S S S S A S R T S S S R S S S R S S R R G G G Y Y R S G
 R H A R S R S R S W S R S R S R R Y S R S R S R G R R H S G G S R D G H R Y S R S P A R R G G Y P R R R S R S R S H S G D R Y R R G
 G R G L R H H S S S R S R S S W S L S P S R S R S L T R S R S H S P S P S Q S R S R S R S R S Q S P S P A R E K L T R P A A S P A V G E
 K L K K T E P A A G K E T G A A K P K L T P Q E K L K L R M Q K A L N R Q F A D K K A A Q E K M I Q Q E H E R Q E R E D E L R A M A R K I
 R M K E R E R R E K E R E E W E R O Y S R O S R S P S P R Y S R E Y S S R R R S R S R S R S P H Y R H

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001278439

ORF Size: 1836 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001278439.2](#)

RefSeq Size: 2064 bp

RefSeq ORF: 1839 bp

Locus ID:	11129
UniProt ID:	Q8N2M8
Cytogenetics:	19q13.32
MW:	70.7 kDa
Gene Summary:	Probably functions as an alternative splicing regulator. May regulate the mRNA splicing of genes such as CLK1. May act by regulating members of the CLK kinase family (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239084