

Product datasheet for **RC239116**

ZNF221 (NM_001297589) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF221 (NM_001297589) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF221
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC239116 representing NM_001297589
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGATTTACCTTCACCTTGAAGTCTTCATTAGGACTCTGCAAATTCCTGAAGTAGAAGGAAAAATGA
 CCACATTCAAAGAGGCAGTGACATTCAAGGATGTGGCTGTGGTCTTCACTGAGGAGGAGCTGGGGCTGCT
 GGACCCTGCCAGAGGAAGCTGTACCGAGATGTGATGCTAGAGAAGTTCAGGAACCTGCTCTCAGTAGGG
 AATCAACCATTCCACCAAGATACTTTCACTTCTTAGGGAAGGAAAAAGTTTGAAGATGAAGACAACAA
 GCCAAAGAGAAGGAATTCAGGAGGCAAGATCCAAATTGAGATGGAGACTGTTCCAGAAGCAGGACCACA
 TGAAGAGTGGTCTGTACGAAATATGGGAACAAATTGCAAGTGACCTAACCAAGTCTCAAACTCCATA
 AGGAACAGCTCTCAGTTCTTCAAGAAGGTGATGTCCCTGCCAGATTGAGGCAAGACTATCTATAAGTC
 ACGTGCAACAGAAACCTTACCGTTGTATGAATGTAACAGTCTTCAGTGATGTTTCTGTCTTTGATCT
 TCATCAACAATCACACTCAGGAGAGAAATCTCATACATGTGGTGAGTGTGAAAAAGCTTCTGTTACAGC
 CCAGCCCTTCATATTCATCAGAGAGTCCATATGGGAGAAAAATGCTATAAGTGTGATGTGTGTGGTAAGG
 AATTTAATCAGAGCTCACATCTGCAAACTCATCAGAGAGTCCATACTGGAGAGAAACCATTCAAATGTGG
 GCAATGTGGGAAAGGCTTCCATAGTAGATCAGCACTTAATGTTTCATTGCAAAATGCACACAGGAGAGAAA
 CCTTATAATTGTGAGGAATGTGGGAAAGCCTTCATTATGATTACAGCTTCAAGAACATCAGAGAATCC
 ATACTGGGAGAGCCATTCAAATGTGATATATGTGGTAAGAGCTTCCGTGTAGATCAAGACTTAATAG
 GCATTCATGGTTACACAGGAGAAAAACCATTCAGATGTGATACATGTGGCAAGAACTTTCGTGAGAGA
 TCAGCACTTAATAGTCATTCCATGGTCCACATAGAAGAGAAGCCATACAAATGTGAGCAATGTGAAAAAG
 GCTTCATTTGTAGGCGAGATTTTGTAGCATCAGATGGTCCACACAGGAGAGAAACCATATAATTGTAA
 AGAATGTGGGAAGACCTTCAGATGGTCCCTCATGTCTTTTGAACCATCAGCAAGTCCACAGTGGACAAAA
 TCCTTCAAATGTGAAGAATGTGGGAAGGATTTTATACAAATTCACGACGATCTTCCATCAGAGATCCC
 ACAATGGAGAAAAAGCCATATACTGTGAGGAGTGTGGTAAGGACTATAAAGGAGGTTGGATCTTGAGTT
 TCACCAGAGGGTCCACACGGGTGAGAGACCTATAATTGTAAGGAATGTGGCAAGAGCTTTGGCTGGGCC
 TCCTGTCTTTGAAACATCAGAGACTCCACAGTGGGAAAAACCTTTCAAATGTGAAGAGTGTGGAAGA
 GATTTACTCAGAGTTCACAACTTCATCCCATCAGACATGCCATACTGGAGAAAAGCTATACAAATGTGA
 GCAGTGTGAGAAGGGGTACAACAGTAAATTTAATCTTGACATGCACCAGAGGGTCCACGGGGGAGAGCGA
 CCCTATAATTGTAAGGAATGTGGAAGAGCTTTGGCTGGGCTTCATGTCTTTGAAACATCAGAGACTCC
 ACAGTGGAGAAAAGCCATTGAAATCTGGAGTGTGGGAAGAGATCTACTCAGAATTCACAGCTTCATTTAC
 ATCAGTAAGTCTATGTGGGAGAAAAGCCATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC239116 representing NM_001297589
 Red=Cloning site Green=Tags(s)


MISPSLELLHSLGCKFPEVEGKMTTFKEAVTFKDVAVVFTEEELGLLDPAQRKLYRDVMLENFRNLLSVG
 NQPFHQDTFHFLLGKEKFWKMKTTTSQREGNSGGKIQIEMETVPEAGPHEEWSCQIWEQIASDLTRSQNSI
 RNSSQFFKEGDVPCQIEARLSISHVQKPYRCNECKQSFSDVSVFDLHQQSHSGEKSHTCGECGKSFCYS
 PALHHIQRVHMGEKCYKCDVCGKEFNQSSHLQTHQRVHTGEKPFKCGQCGKFHRSRSLNVHCKLHTGEK
 PYNCEECGKAFIHDSQLQEHQRIHTGEKPFKCDICGKSFRVRSRLNRHSMVHTGEKPFRCDCGKNFRQR
 SALNSHSMVHIEEKPYKCEQCGKGFICRRDFCKHQMVTGEKPYNCKEKGKTFRWSSCLLNHQVHSGQK
 SFKCEECGKGFYNSRRSSHQRSHNGEKPYNCEECGKDYKRRDLFEHQRVHTGERPYNCKEKGKSGFWA
 SCLLKHQRLHSGEKPFCCEECGKRFTQSSQLHSHQTCHTGEKLYKCEQCEKGYNSKFNLDMMHQRVHGGER
 PYNCKEKGKSGFWASCLLKHQRLHSGEKLKSGVWEEIYSEFTASFTSVSLCGRKAI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001297589

ORF Size: 1851 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001297589.1](#), [NP_001284518.1](#)


RefSeq Size: 2288 bp

RefSeq ORF: 1854 bp

Locus ID: 7638

UniProt ID: Q9UK13
 Cytogenetics: 19q13.31
 MW: 71.6 kDa
 Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239116