

## Product datasheet for **RC239119**

### ZNF539 (ZNF254) (NM\_001278661) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF539 (ZNF254) (NM_001278661) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF254
Synonyms:	BMZF-5; HD-ZNF1; ZNF91L; ZNF539
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC239119 representing NM\_001278661  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTTAGAGAACTACAGAAACCTGGCCTTCTGGGTATTGCTGTCTCTAAGCCAGACCTGATCACCTGTC  
 TGGAAACAAGGAAAGAGCCCTGGAATATGAAGCGACATGAGATGGTGGATGAACCCAGGTATGTGTCC  
 TCATTTTCTCAAGACCTTTGGCCAGAGCAGGGCATGGAAGATCTTTTCAAAAAGCAACTACTGAGAAGA  
 TATGGAAAATATGGACATGAGAATTTACAGTTAAGAAAAGGCTGTAAAAGTGTGGATGAGTAAAGTGA  
 ACAAAGAAGGTTATAATGGACTTAACAGTGTTTCACTGCCCAGAGCAAAGTATTTCAATGTGATAA  
 ATATTTGAAAGTCTTCTATAAATTTTAAATTCAAACAGACCTAAGATAAGACATACTGAAAAGAAATCT  
 TTCAAATGTA AAAACGTGTCAAATTTTTCATGCTTTACATAAAAACCCACACAAAAGCATTATC  
 ATAGAGAGAAGTCTACAAATGTAAGAATGTGAAAAACCTTTAATTGGTCTCAACCTTACTAATCA  
 TAGGAAAAATTTACTGAAGAGAAACCTTACAAATGTGAAGAATAACAAATCTCTAAGCAACTCTCA  
 ACCCTTACTACACATGAAATAATTCATGCTGGAGAGAACTCTACAAATGTGAAGAATGTGGTGAAGCTT  
 TTAATCGATCCTCAAATCTTACTACACATAAGATAATTCATACTGGAGAGAAACCTTACAAGTGAAGA  
 ATGTGGCAAAGCATTTATCTGGTCTCAACCTTACTGAGCATAAGAAAATTCATACTAGAAAAGAAACCC  
 TACAAGTGTGAAGAATGTGGCAAAGCATTTATATGGTCTCAACCTAAGTACATAAGAGGATGCACA  
 CTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCTTTTAGCCAGTCTCAACCTTACTACACA  
 TAAGATAATTCATACTGGAGAGAAACGCTACAAATGCTTAGAATGTGGCAAAGCTTTAAGCAACTCTCA  
 ACTCTTACTACACATAAAAATAATTCATGTTGGAGAGAACTCTACAAATGTGAAGAATGCGGCAAAGGTT  
 TTAATCGATCTTCAAATCTTACTACACATAAGATAATTCATACTGGAGAGAAACCTTACAAGTGTGAAGA  
 ATGTGGCAAAGCATTTATCTGGTCTCAACCTTACTAAACATAAGAGAATTCATACTAGAGAGAAACCC  
 TACAAATGTGAAGAATGTGGCAAAGCATTTATATGGTCTCAACCTAAGTACATAAGAGGATGCACA  
 CTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAATCTTTTAGCCAATCCTCAACCTTACTACACA  
 TAAGATAATTCATACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCTTTAAGTGGTCTCA  
 ACTCTTACTAAACATAAGATAATTCATACTGAAGAGAAACCTACAAATGTGAAAATGTGGCAAAGCTT  
 TTAAGCAGTCTTCAATCCTTACTAACCATAAGAGAATTCATACTGGAGAGAAACCTATAAATGTGAAGA  
 ATGTGGCAAATCTTTAACCAGTCTTCACTTTTACTAAACATAAGGTAATTCATACTGGAGTAAAACCC  
 TACAAATGTGAAGAATGTGGCAAAGCATTTTCTGGTCTCAACCTAAGTAAACATAAGAGAATTCATA  
 CTGGAGAGCAACCTACAAATGGAAAAATTTGGCAAAGCTTTAATCGTCTCGCACCTACCCACAGA  
 TAAGATAACTCATTGGAGAGAAATCTTACAAGTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>RC239119 representing NM\_001278661  
 Red=Cloning site Green=Tags(s)

MLENYRNLAFLGI AVSKPDLITCLEQKPEPWNMKRHEMVDPEPPGMCPHF AQDLWPEQGMEDSFQKAILRR  
 YGKYGHENLQLRKGCKSVDEYKVNKEGYNGLNQCFTTAQSKVFQCDKYLKVFYKFLNSNRPKIRHTEKKS  
 FKCKKRVKLF CMLSHKTQH KSIYHREKSYKCKEKGKTFNWSSTLTNHRKIYTEEKPYKCEEYKNSPKQLS  
 TLTTHEIIHAGEKLYKCEECEGAFNRSSNL TTHKIIHTGEKPYKCEEKGAFIWSSTL TEHKKIHTRKKP  
 YKCEEKGAFIWSSTL TRHKRMHTGEKPYKCEEKGAFS QSSTL TTHKIIHTGEKRYKCECGKAFKQLS  
 TLTTHKIIHVGEKLYKCEECEGKFNRSSNL TTHKIIHTGEKPYKCEEKGAFIWSSTL TKHKRIHTREK  
 YKCEEKGAFIWSSTL TRHKRMHTGEKPYKCEEKGSFSQSSTL TTHKIIHTGEKPYKCEEKGAFNWSS  
 TLTTHKIIHTEEKPYKCEKCGKAFKQSIL TNHKRIHTGEKPYKCEEKGSFNRSSTFTKHKVIHTGVK  
 YKCEEKGAFIWSSTLTKHKRIHTGEQPYKWEKFGKAFNRSSHL TTDKITHWREILQV


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001278661

ORF Size: 1854 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001278661.1</a> , <a href="#">NP_001265590.1</a>
<b>RefSeq Size:</b>	4262 bp
<b>RefSeq ORF:</b>	1857 bp
<b>Locus ID:</b>	9534
<b>Cytogenetics:</b>	19p12
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	72.8 kDa
<b>Gene Summary:</b>	Zinc finger proteins have been shown to interact with nucleic acids and to have diverse functions. The zinc finger domain is a conserved amino acid sequence motif containing 2 specifically positioned cysteines and 2 histidines that are involved in coordinating zinc. Kruppel-related proteins form 1 family of zinc finger proteins. See ZFP93 (MIM 604749) for additional information on zinc finger proteins.[supplied by OMIM, Jul 2002]