

Product datasheet for **RC239140**

BMAL1 (ARNTL) (NM_001297722) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BMAL1 (ARNTL) (NM_001297722) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ARNTL
Synonyms:	bHLHe5; BMAL1; BMAL1c; JAP3; MOP3; PASD3; TIC
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239140 representing NM_001297722
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCAGACCAGAGAATGGACATTTCTCAACCATCAGTGATTTTCATGTCCCCGGGCCCCACCGACCTGC
 TTTCCAGCTCTCTTGGTACCAGTGGTGGATTGCAACCGCAAACGGAAAGGCAGCTCCACTGACTACCA
 AGAAAGCATGGACACAGACAAAGATGACCTCATGGAAGGTTAGAATATACAGAACCAAGGAAGGATA
 AAAAATGCAAGGGAAGCTCACAGTCAGATTGAAAAGCGGCGTCGGGATAAAAATGAACAGTTTTATAGATG
 AATTGGCTTCTTTGGTACCAACATGCAACGCAATGTCCAGGAAATTAGATAAACTTACTGTGCTAAGGAT
 GGCTGTTCCAGCACATGAAAACATTAAGAGGTGCCACCAATCCATACACAGAAGCAAACCTACAAACCACT
 TTTCTATCAGACGATGAATTGAAACACCTATTCTCAGGGCAGCAGATGGATTTTTGTTTGTGCTAGGAT
 GTGACCGAGGGAAGATACTCTTTGTCTCAGAGTCTGTCTTCAAGATCCTCAACTACAGCCAGAATGATCT
 GATTGGTCAGAGTTTGTGGTACTACCTGCATCCTAAAGATATTGCCAAAGTCAAGGAGCAGCTCTCTCC
 TCTGACACCGCACCCCGGAGCGGCTCATAGATGCAAAAACCTGGACTTCCAGTAAAACAGATATAACCC
 CTGGGCCATCTCGATTATGTTCTGGAGCAGCAGCTTCTTTCTTGTAGGATGAAGTGTAAACAGGCCTTC
 AGTAAAGGTTGAAGACAAGGACTTCCCCTCTACCTGCTCAAAGAAAAAGCAGATCGAAAAAGCTTCTGC
 ACAATCCACAGCAGGCTATTTGAAAAGCTGGCCACCCACAAAGATGGGGCTGGATGAAGACAACGAAC
 CAGACAATGAGGGGTGTAACCTCAGTGCCTCGTCGCAATTGGACGACTGCATTCATGTAGTTCCACA
 ACCAGTGAACGGGAAATCAGGGTGAATCTATGGAATATGTTTCTCGGCACGCGATAGATGGAAAGTTT
 GTTTTTGTAGACCAGAGGCAACAGCTATTTTGGCATATTTACCACAAGAACTTCTAGGCACATCGTGT
 ATGAATATTTTCAACAAGATGACATAGGACATCTTGCAGAATGTCATAGGCAAGTTTTACAGACGAGAGA
 AAAAATTACAACATAATTGCTATAAATTTAAAATCAAAGATGGTTCTTTTATCACACTACGGAGTCGATGG
 TTCAGTTTCATGAACCTTGGACCAAGGAAGTAGAATATATTGTCTCAACTAACACTGTTGTTTTAGCCA
 ACGTCTTGAAGGCGGGGACCAACCTTCCCACAGCTCACAGCATCCCCCACAGCATGGACAGCATGCT
 GCCCTCTGGAGAAGTGGCCAAAGAGGACCCACCCACTGTTCCAGGGATTCCAGGGGAAACCCGGCT
 GGGGCAGGAAAAATAGGCCGAATGATTGCTGAGGAAATCATGGAATCCACAGGATAAGAGGGTCATCGC
 CTTCTAGCTGGCTCCAGCCATTGAACATCACGAGTACGCCTCCCCTGATGCCTCTTCTCCAGGAGG
 CAAGAAGATTTAAATGGAGGACTCCAGACATTCCTTCCAGTGGCTACTATCAGGCCAGGCTCAGGAG
 AACCCAGTTATCCATATCTGATAGTTCTTCTATTCTTGGTGAACCCCCACATAGGTATAGACATGA
 TTGACAACGACCAAGGATCAAGTAGTCCAGTAATGATGAGGCAGCAATGGCTGTATCATGAGCCTCTT
 GGAAGCAGATGCTGGACTGGGTGGCCCTGTTGACTTTAGTGACTTGCCATGGCCGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC239140 representing NM_001297722
 Red=Cloning site Green=Tags(s)

MADQRMDISSTISDFMSPGPTDLLSSSLGTSGVDCNRKRKGSSTDYQESMDTDKDDPHGRLEYTEHQGRI
 KNAREAHSQIEKRRRDKMNSFIDELASLVPTCNAMSRKLDKLTVLRMAVQHMKTLRGATNPYTEANYKPT
 FLSDELKHLILRAADGFLFVVGCDRGKILFVSESVFKILNYSQNDLIGQSLFDYLHPKDIKVKQLSS
 SDTAPRERLIDAKTGLPVKTDITPGPSRLCSGARRSFFCRMKNRPSVKVEDKDFPSTCSKKKADRKSFC
 TIHSTGYLKSWPPTKMGLDEDNEPDNEGCNL SCLVAIGRLHSHVVPQVNGEIRVKSMEYVSRHAIDGKF
 VFVDQRATAILAYLPQELLGTSCYEYFHQDDIGHLAECHRQVLQTREKITTNCYKFKIKDGSFITLRSRW
 FFSFMPWTKEVEYIVSTNTVVLANVLEGGDPTFPQLTASPHSMDSMLPSGEGGPKRTHPTVPGIPGGTRA
 GAGKIGRMIAEEIMEIHRIRGSSPSSCGSSPLNITSTPPPDASSPGGKILNGGTPDIPSSGLLSGQAQE
 NPGYPYDSSSILGENPHIGIDMIDNDQSSSPSNDEAAMAVIMSLLEADAGLGGPVDFSDLPWPL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297722

ORF Size: 1878 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001297722.1 , NP_001284651.1
RefSeq Size:	2811 bp
RefSeq ORF:	1881 bp
Locus ID:	406
UniProt ID:	O00327
Cytogenetics:	11p15.3
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Circadian rhythm - mammal
MW:	69.2 kDa
Gene Summary:	The protein encoded by this gene is a basic helix-loop-helix protein that forms a heterodimer with CLOCK. This heterodimer binds E-box enhancer elements upstream of Period (PER1, PER2, PER3) and Cryptochrome (CRY1, CRY2) genes and activates transcription of these genes. PER and CRY proteins heterodimerize and repress their own transcription by interacting in a feedback loop with CLOCK/ARNTL complexes. Defects in this gene have been linked to infertility, problems with gluconeogenesis and lipogenesis, and altered sleep patterns. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]