

Product datasheet for **RC239153**

GOLGA8J (NM_001282472) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA8J (NM_001282472) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GOLGA8J
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239153 representing NM_001282472
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGAAGAACTCAACACAACAAATTGGCTGCAGCCAAGAAAAAGTTAAAAGAATATTGGCAGAAAA
 ACAGCCCTAGAGTTCCAGCAGGAGCGAACAGGAACAGGAAAAACAAATGGCAGTATCCCTGAGAAAGCCAC
 TTCTGGTGGTTGCCAGCCACCTAGGGATTAGCAACAGGTTTCCACAGGGAAGGCCCTACATCATCTGCT
 ACCCTGAAAGATCTGGAGAGCCCGTGCCAGAAGCAGCAGTAGTCCTGGATTCAAGGTCCGTAGAAATCA
 GTCAACTGAAGAACACCATCAAATCTCTGAAACAACAGAAGAAACAAGTGAACATCAGCTGGAAGAAGA
 AAAGAAAGCAAACAACAAGAAACAGAAAGCCAAAAGGGTGCTAGAGGTTCAAATCCAGACATTGAACATA
 CAGAAAGAGGAACTAAATACGGACCTGTACCACATGAAACGTTCTCTCAGATACTTTGAAGAAAAGTCCA
 AGGATCTGGCTGTCCGCTGCAACATTCATTGCAGCGTAAAGGAGAGTTAGAGAGTGTCTCTCTAATGT
 CATGGCCACACAGAAGAAGAAGGCAAACAGTTGTCCAGCCGAGCAAGCACGTACGGAGTGGAAGTTA
 GAGCAGTCCATGCGGGAGGAGGCACTACTGAAAGTGCAGCTGACACAGTTCAAGGAGTCATTTCAACAAG
 TCCAATTAGAAAGAGATGAGTATTCTGAACATCTAAAAGGAGAGAGGGCCCGGTGGCAGCAGAGGATGAG
 AAAATGTGCGAGGAGATTTGCACATTAAGAAAGAGAAGCAGCAAGATATGCGTCGGGTAGAGAAGCTG
 GAGAGGAGCTTGTCCAACTCAAAAACAGATGGCTGAACCCCTTGCCCCCGAGCCCCCAGCAGTGCCCT
 CTGAGGTGGAGCTGCAGCACCTGAGGAAGGAAGTAGAGAGAGTGGCAGGAGAGCTCCAGGCCCAGGTCAA
 AAACAATCAGCGCATAGTCTCTGAACACGCGACAAGAAGAGAGGATTTCGGGAGCAGGAAGAGAGGCTT
 CGGAAGCAGGAGGAGAGGATTGAGGAGCAGCACAAGAGCCTTCAGCAGCTGGCCAAGCCACAGAGCGTCT
 TCAAGGAGCCGAACAATGAGAACAAGAACGCACTGCAGTTGGAGCAGCAAGTAAAGGAGCTACAGGAGAA
 GCTTGGCGAGGAGCACCTGGAAGCTGCCAGCCAGCAGCAAGCAGCAGCTAACGGGCCAGCTGAGCCTCATG
 GCTCTCCCTGGGAAGGACACGGAGGAGAACATCTGGACAGTGAGGGGGAGGAGGCACCTCGGCCCATGC
 CGAGTGTCCAGAGGACCCGGAGAGCAGGGAGGCCATGAGCAGCTTTATGGACCACCTGGAGGAGAAGGC
 AGACCTGAGTGAGCTGGTGAAGAAAAAGAACTCTGCTTCATCCACCACTGGCGAGAGAGATGCCATCAG
 AAAACCCATCACCTTTTATCAGAACAGGGGGCCGTGCCAAAGATGCGGCACTGGGAGGAGGACACCATC
 AGGCTGGAGCTCAGGGAGGAGATGAAGGTGAAGCTGCTGGAGCTGCAGCAGATGGTATTGCGGCTTACAG
 CAACTACAACAATGGGCACAGAAAATTCCTGGCCGCTGCCCAACTCTGCTGATGAGCCCGGTCCAGGA
 GCCCCAGCCCCCAGGAGCTTGGGGCTGCAGACAAGCATGGTCATCTTTGTGAGGTGAGCCTCACCTCCT
 CTGCCCCAAGGAGAGGCCAGGGAGGATCCTCTCCTTGACAAGCCTACTGCACAGCCGATCGTGCAGGACCA
 CCAGGAGCACCCAGGCTTGGGCAGCAACTGCTGTGTGCCATTCTGTGTCTGGGCTTGGCTGCCAAGAAGA
 AGGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239153 representing NM_001282472
 Red=Cloning site Green=Tags(s)


MAEETQHNKLA AAKKLKEYWQKNSPRVPAGANRNKTNNGSIPEKATSGGCQPPRDSATGFHREGPTSSA
 TLKDLESPCQERAVVLD SRSVEISQLKNTIKSLKQKKQVEHQLEEEKANNKKQKAKRVLEVQIQT LNI
 QKEELNTDLYHMKRSLRYFEESKDLAVRLQHSLQRKGELESVLSNVMATQKKKANQLSSRSKARTEWKL
 EQSMREEALLKVQLTQFKESFQQVQLERDEYSEHLKGERARWQQRMRKMSQEICTLKKEKQDMRRVEKL
 ERSLSKLNQMAEPLPPEPPAVPSEVELQHLRKE LERVAGELQAQVKNQRI SLLNQRQEERIREQEERL
 RKQEERIQEQHKSLQQLAKPQSVFKEPNNENKNALQLEQQVKELQEKLGEEHLEAASQQNQQLTAQLSLM
 ALPGEGHGGHELDSEGEAPRPMPSPVEDPESREAMSSFMDHLEEKADLSELVKKKELCFIHHWRERCHQ
 KTHHLLSEPGGRAKDAALGGGHHQAGAQQGDEGEAAGAAADGIAAYSNNYNGHRKFLAAAHNSADEPGPG
 APAPQELGAADKHGHLCEVSLTSSAQGEAREDP LLDKPTAQPIVQDQHEHPGLGSNCCVPFLCWAWLPRR
 RR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001282472

ORF Size: 1896 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001282472.1](#), [NP_001269401.1](#)

RefSeq Size: 5199 bp

RefSeq ORF: 1899 bp


Locus ID: 653073

UniProt ID: [A6NMD2](#)

Cytogenetics: 15q13.2

MW: 72.2 kDa

Product images:



Circular map for RC239153