

Product datasheet for **RC239156**

GOLGA8H (NM_001282490) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA8H (NM_001282490) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GOLGA8H
Synonyms:	GOLGA6L11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239156 representing NM_001282490
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGAAGAACTCAACACAACAAATTGGCTGCAGCCAAGAAAAAGTTAAAAGAATATTGGCAGAAAA
 ACAGCCCTAGAGTTCCAGCAGGAGCGAACAGGAACAGGAAAAACAAATGGCAGTGTCCCTGAGAAAGCCAC
 TTCTGGTGGTTGCCAGCCACCTGGGGATTAGCAACAGGTTTCCACAGGGAAGGCCCTACATCATCTGCT
 ACCCTGAAAGATCTGGAGAGCCCGTGCCAGAAGCAGCAGTAGTCCTGGATTCAAGGTCCGTAGAAATCA
 GTCAACTGAAGAACACCATCAAATCTCTGAAACAACAGAAGAAACAAGTGAACATCAGCTGGAAGAAGA
 AAAGAAAGCAAACAACAAGAAACAGAAAGCCAAAAGGGTGCTAGAGGTTCAAATCCAGACATTGAACATA
 CAGAAAGGGAACTAAATACGGACCTGTACCACATGAAACGTTCTCTCAGATACTTTGAAGAAAAGTCCA
 AGGATCTGGCTGTCTGCCTGCAACATTCATTGCAGCGTAAAGGAGAGTTAGAGAGTGTCTCTCTAATGT
 CATGGCCACACAGAAGAAGAAGGCAAACAGTTGTCCAGCCGAGCAAGCACGTACGGAGTGGAAGTTA
 GAGCAGTCCATGCGGGAGGAGGCACTACTGAAAGTGCAGCTGACACAGTTGAAGGAGTCATTTCAACAAG
 TCCAATTAGAAAGAGATGAGTGTGCTGAACATCTAAAAGGAGAGAGGGCCCGGTGGCAGCAGAGGATGAG
 AAAATGTGCGAGGAGATTTGCACATTAAGAAAGAGAAGCAGCAAGATATGCGTCGGGTAGAGAAGCTG
 GAGAGGAGCTTGTCCAACTCAAAAACAGATGGCTGAACCCCTTGCCCCCGGAGCCCCCAGCAGTGCCCT
 CTGAGGTGGAGCTGCAGCACCTGAGGAAGGAAGTAGAGAGAGTGGCAGGAGAGCTCCAGGCCCAGGTCAA
 AAAGAATCAGCGCATAGTCTCCTGAACCGCACAAGAAGAGAGGATTCAGGAGCAGGAAGAGAGGCTT
 CGGAAGCAGGAGGAGAGGATTGAGGAGCAGCACAAGAGCCTTCAGCAGCTGGCCAAGCCACAGAGCGTCT
 TCGAGGAGCCGAACAATGAGAACAAGAACGCACTGCAGTTGGAGCAGCAAGTAAAGGAGCTACAGGAGAA
 GCTTGGCGAGGAGCACCTGGAAGCTGCCAGCCAGCAGAACCAGCAGCTAACGGGCCAGCTGAGCCTCATG
 GCTCTCCCTGGGGAAGGACACGGAGGAGAACATCTGGACAGTGAGGGGGAGGAGGCACCTCGGCCCATGC
 CGAGTGTCCAGAGGACCCGGAGAGCAGGGAGGCCATGAGCAGCTTTATGGACCACCTGGAGGAGAAGGC
 AGACCTGAGTGAGCTGGTGAAGAAAAAGAACTCTGCTTCATCCACCACTGGCGAGACAGATGCCATCAG
 AAAACCCATCACCTTTTATCAGAACAGGGGGCCGTGCCAAAGATGCGGCACTGGGAGGAGGACACCATC
 AGGCTGGAGCTCAGGGAGGAGATGAAGGTGAAGCTGCTGGAGCTGCAGCAGATGGTATTGCGGCTTACAG
 CAACTACAACAATGGGCACAGAAAATTCCTGGCCGCTGCCACAACCTCTGCTGATGAGCCCGGTCCAGGA
 GCCCCAGCCCCCAGGAGCTTGGGGCTGCAGACAAGCATGGTGATCTTTGTGAGGTGAGCCTCACCTCCT
 CTGCCCCAAGGAGAGGCCAGGGAGGATCCTCTCCTTGACAAGCCTACTGCACAGCCGATCGTGCAGGACCA
 CCAGGAGCACCCAGGCTTGGGCAGCAACTGCTGTGTGCCATTGTTGTGCTGGGCTTGGCTGCCAAGAAGA
 AGGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239156 representing NM_001282490
 Red=Cloning site Green=Tags(s)


MAEETQHNKLA AAKKKLKEYWQKNSPRVPAGANRNKRTNGSVPEKATSGGCQPPGDSATGFHREGPTSSA
 TLKDLESPCQERAVVLD SRSVEISQLKNTIKSLKQKKQVEHQLEEEKANNKKQKAKRVLEVQIQT LNI
 QKGK LNTDLYHMKRSLRYFEESKDLAVCLQHSLQRKGELESVLSNVMATQKKKANQLSSRSKARTEWKL
 EQSMREEALLKVQLTQLKESFQQVQLERDECAEHLKGERARWQQRMRKMSQEICTLKKEKQDMRRVEKL
 ERSLSKLNQMAEPLPPEPPAVPSEVELQHLRKE LERVAGELQAQVKNQRI SLLNQRQEERI QEQEERL
 RKQEERI QE QHKS LQQLAKPQSVFEEPNNENKNALQLEQQVKELQEKLGEHLEAASQQNQQLTAQLSLM
 ALPGEGHGGHLDSEGE EAPRPMPSVPEDPESREAMSSFMDHLEEKADLSELVKKKELCFIHHWRDRCHQ
 KTHHLLSEPGGRAKDAALGGGHHQAGA QGGDEGEAAGAAADGIAAYSNNNGHRKFLAAAHNSADEPGPG
 APAPQELGAADKHGDLCEVSLTSSAQGEAREDP LLDKPTAQPIVQDQHEHPGLGSNCCVPLLCWAWLPRR
 RR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001282490

ORF Size: 1896 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001282490.1](#), [NP_001269419.1](#)

RefSeq Size: 5261 bp

RefSeq ORF: 1899 bp


Locus ID: 728498

UniProt ID: [P0CJ92](#)

Cytogenetics: 15q13.2

MW: 71.8 kDa

Product images:



Circular map for RC239156