

Product datasheet for **RC239173**

MTMR12 (NM_001294344) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MTMR12 (NM_001294344) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MTMR12
Synonyms:	3-PAP; PIP3AP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239173 representing NM_001294344
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGCTGGGAAAGGAGTAGTCGGCGGTGGCGGCGGCACCAAGGCCCCCAAGCCCTCCTTCGTGTCGTACG
 TACGCCCTGAGGAAATTCACACAAACGAAAAGGAAGTAACAGAGAAGGAAGTAACCTTCACTTGTGGC
 AGGTGAACAGCTGCTTTGTGAAGCCAGCACAGTACTGAAGTATGTCCAGGAAGATTCTGTGAGCATGGG
 GTCTATGGGAGGCTTGTCTGCACAGACTTCAAGATTGCCTTCTTGGGTGATGATGAATCTGCATTGGATA
 ATGATGAACTCAATTTAAGAATAAGGTTATAGGAGAAAATGACATTACACTCCACTGTGTGATCAGAT
 TTATGGAGTGTGATGAGAAAAAGAACTCTCTTGGACAACCTGAAGAAATACCCTGAGAAGCTCATC
 ATCCACTGCAAAGACCTTCGAGTGTCCAGTTTGTCTGAGGTACACAAAGGAAGAGGAAGTCAAAAGGA
 TTGTGAGTGGCATAATTCATCATACCCAGGCTCCTAACTGCTTAAACGATTATTTCTGTTTCTATGC
 GACTGCTGCACAAAACAATACAGTCACTGATCCCAAGAACCATAACCGTAATGTTTGACACACTTAAGGAC
 TGGTGTGGGAACTGGAACGGACCAAGGCAACATGAAGTACAAAGCAGTGAGTGTCAACGAAGGCTATA
 AAGTCTGTGAGAGATTGCCAGCATACTTGTGTGCTCCCACTCTTCTCTGAAGAGAATGTGCAGCGCTT
 TCAGGGTCATGGCATACCAATATGGTGTGGTCTGCCACAATGGAAGTGCCCTTTTGAAAATGTGACGA
 CTGCCCCAAGAACAGGATGACGGCATTTTACAAATCCAAAGAGCTTCTTAGATGGAATTTACAAGACCA
 TCCACAGGCCACCCTATGAAATTGTTAAACCGAAGACCTGTCAAGCAACTTCTGTCCCTGCAGGAAAT
 CCAGACTGCATACTCTAAATTTAAACAGCTATTTCTGATAGATAACAGTACTGAATTTGGGACACAGAT
 ATAAAATGGTTTTCTCTGTTGAAAGTAGCAGCTGGCTTGACATAATCAGACGTTGCCTGAAAAAGCAA
 TAGAGATTACAGAATGTATGGAAGCACAAAACATGAATGTTCTTCTTTAGAGGAGAATGCATCCGACCT
 CTGCTGTCTCATTTCTCTCTGGTGCAACTGATGATGGACCCCACTGCAGAACCAAGATTGGTTTTCCAG
 AGCCTCATCCAAAGGAGTGGGTGATGGGTGGCCACTGTTTCTTGGATCGCTGCAACCATCTCCGCCAGA
 ACGACAAAGAGGAGCATCAACGACAACCTTCTTGGCACTTACACAATCTAAGTCATCTCCCAAGAGG
 ATTTTTCAGGGAAGAACAGATCATTTAATTAACCTTCTGGGCAAGAGAATTAGCAACTTATTAAT
 TCTTCTGATGAGCTCCAAGACAACCTTCGAGAGTTCTATGACAGCTGGCATAGCAAGTCCACTGACTACC
 ATGTTTTGTTGTTACCGCATATCGAGGGGCCGAAATCAAAGTCTGGGCCAGCGCTACCTACGTTGGAT
 TCCAGAAGCCCAATCTAGGTGGTGGCAAGTGGCCACTCTGAGCAAACTCTTGAAATGATGGAGGAA
 GTCCAGAGTTTACAAGAGAAGATCGACGAGAGACACCACAGCCAGCAGGCCCCAGGCTGAGGCCCCCT
 GCCTGCTGAGGAACTCTGCCCGCTCTCGTCTTTGTTTCTTTCTGCTCTGCTCCAGCGACATTCTCTAA
 GCCCGTCTTACCCACAGTGGCTGGAAAGCTTTGGGAGATGAAGACGATTGGCCAAACGAGAAGATGAG
 TTCGTGGACCTAGGGGATGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239173 representing NM_001294344
Red=Cloning site Green=Tags(s)

MLLKGKGVGGGGGTAKAPKPSFVSYPREEIHTNEKEVTEKEVTLHLLPGEQLLCEASTVLKYVQEDSCQH
VYGRVLCTDFKIAFLGDDESDLNDETQFKNKVIGENDITLHCVQDIYGVFDEKKKTLFGQLKKYPEKL
IHCKDLRVFQFLRYTKEEEVKRIVSGIIHHTQAPKLLKRLFLFSYATAAQNNVTDPKNHTVMFDTLKD
WCWELERTKGNMKYKAVSVNEGKVCERLPAYFVPTPLPEENVQRFQGHGIPWCWSCHNGSALLKMSA
LPKEQDDGILQIQKSFLDGIYKTIHRPPYEIVKTEDLSSNFLSLQEIQTAYSKFKQLFLINDSTEFWDT
IKWFSLESSWLDIRRCLKKAIEITECMEAQNMMNVLLLEENASDLCCLISSVLQMLMDPHCRTRIGFQ
SLIQKEWVMGGHCFLDRCNHLRQNDKEEHQRQLSLPLTQSKSSPKRGFFREETDHLIKNLLGKRISKLIN
SSDELQDNFREFYDSWHSKSTDYHGLLLPHIEGPEIKVWAQRYLRWIPEAQILGGGQVATLSKLEEMEE
VQSLQEIKIDERHHSQAPQAEAPCLLRNSARLSSLFPFALLQRHSSKPVLP TSGWKALGDEDDLAKREDE
FVDI.GDV

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001294344

ORF Size: 1911 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001294344.2</u>
RefSeq Size:	4868 bp
RefSeq ORF:	1914 bp
Locus ID:	54545
UniProt ID:	<u>Q9C0I1</u>
Cytogenetics:	5p13.3
Protein Families:	Druggable Genome
MW:	73.7 kDa
Gene Summary:	Phosphatidylinositide 3-kinase-derived membrane-anchored phosphatidylinositides, such as phosphatidylinositol 3-phosphate (PtdIns(3)P), regulate diverse cellular processes. The protein encoded by this gene functions as an adaptor subunit in a complex with an active PtdIns(3)P 3-phosphatase. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2014]

Product images:



Circular map for RC239173