

Product datasheet for **RC239176**

ACSBG2 (NM_001289179) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACSBG2 (NM_001289179) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ACSBG2
Synonyms:	BGR; BRGL; PRTD-NY3; PRTDNY3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239176 representing NM_001289179
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**

ATGACTGGAACCCCAAAGACTCAAGAAGGAGCTAAAGATCTTGAAGTAGACATGAATAAAACAGAAGTTA
 CTTCCAGGCTGTGGACCACCTGTCGAGATGGAGAAGTCCTTCTGAGGCTATCCAAACACGGACCAGGCCA
 TGAGACCCCGATGACCATCCCTGAATTTTTCGAGAGTCAAGTCAACCGATTGGAACTTATCCAGCCCTC
 GCATCCAAGAATGGCAAAAAGTGGGAAATCTGAATTTCAACCAGTACTATGAGGCTTGTGCGAAGGCTG
 CAAAATCCTTGATCAAGCTGGGTTTGGAGCGTTTCCACGGAGTTGGTATCCTGGGGTTAACTCTGCAGA
 GTGGTTTATCACTGCTGTTGGTGCCATCTAGCCGGGGTCTTTGTGTTGGTATTTATGCCACCAACTCT
 GCCGAGGTTTGTCAATATGTCATCACTCATGCCAAAGTGAACATCTTGCTGGTTGAGAATGATCAACAGT
 TACAGAAAATCCTTCGATTCCACAGAGCAGCCTAGAGCCCTAAAAGCGATCATCCAGTACAGACTGCC
 AATGAAGAAGAACAACACTTGTACTCTTGGGATGATTTTATGGAAGTGGCAGAAGTATCCTGACACC
 CAACTGGAGCAGGTCATCGAGAGCCAGAAGGCGAATCAATGCGCAGTGTCTATCACTTCAGGGACCA
 CAGGCATACCCAAGGGAGTGATGCTCAGTCATGACAACATCACGTGGATTGCAGGAGCAGTGACAAAGGA
 CTTTAACTGACAGACAAGCATGAGACGGTGGTTAGCTACCTCCCACTCAGCCATATTGCAGCACAGATG
 ATGGACATCTGGGTACCCATAAAGATTGGGGCGCTCACATACTTTGCTCAAGCAGATGCTCTCAAGGGCA
 CCTTGGTAAGTACTCTAAAGGAGGTAACCTACTGTCTTCATTGGAGTGCCTCAAATTTGGGAGAAGAT
 ACATGAGATGGTGAAGAAAAATAGTGCCAAGTCCATGGGCTTGAAGAAGAAGGCATTTCGTGTGGGCAAGA
 AACATTGGCTTCAAGGTCAACTCAAAAAAGATGTTGGGGAAATATAATACTCCCGTGAGTACCGCATGG
 CTAAGACTCTCGTGTTCAGCAAAGTCAAGACATCCCTTGGCTTGGATCACTGTCACTCTTTATCAGTGG
 GACTGCGCCCTCAACCAAGAGACTGCCGAGTTCTTTCTAAGCTTGGACATACCTATAGCGGAGTTGTAT
 GGGTTGAGTGAGAGCTCGGGACCCACACGATATCCAACCAGAATAACTACAGGCTTCTAAGCTGTGGCA
 AGATCTTGACTGGGTGTAAGAATATGCTGTTCCAGCAGACAAGGATGGCATTGGGGAGATCTGCCTCTG
 GGGTAGGCACATCTTCATGGGCTATCTGAAAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAG
 TGGCTACACTCTGGGATCTGGGCCAGCTGGACGGTCTGGGTTTCTCTATGTCACCGGCCACATCAAAG
 AAATCCTTATCACTGCTGGTGGTGAAGTGTGCCCCCATTCTGTTGAGACCTTGGTTAAGAAGAAGAT
 CCCCATCATCAGTAACGCCATGTTAGTAGGAGATAAACTGAAGTTTCTGAGCATGTTGCTGACGCTGAAG
 GCATCCACCGTACTGAGATTGTGAAGCAGCAAGACCCCTGGTCTACAAGGCCATCCAGCAAGGCATCA
 ATGCTGTGAACCAGGAAGCCATGAACAATGCACAGAGGATTGAAAAGTGGGTATCTTGGAGAAGGACTT
 TTCCATCTATGGTGGAGAGCTAGGTCCAATGATGAACTTAAGAGACATTTTGTAGCCGAGAAATACAAA
 AAACAAATTGATCACATGTACCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239176 representing NM_001289179
Red=Cloning site Green=Tags(s)


MTGTPKTQEGAKDLEVDMNKTEVTPRLWTTCTRDGEVLLRLSKHGPGETPMTIPEFFRESVNRFGTYPAL
ASKNGKKWEILNFNQYYEACRKAASLKLGLERFHGVGILGFNSAEWFITAVGAILAGGLCVGIYATNS
AEVCQYVITHAKVNILLVENDQQLQKILSIPQSSLEPLKAIQYRLPMKNNNLYSWDDFMELGRSIPDT
QLEQVIESQKANQCAVLIYTS GTTGIPKGVMLSHDNITWIAGAVTKDFKLT DKHETVVSYLPLSHIAAQM
MDIWVPIKIGALTYFAQADALKGTLVSTLKEVKPTVF IGV PQIWEKIHVMKNSAKSMGLKKKAFVWAR
NIGFKVNSKKMLGKYNTPVSYRMAKTLVFSKVKTSGLDHCFSISGTAPLNQETAEFFLSLDIPIGELY
GLSESSGPHTISNQNNYRLSCGKILTGCKNMLFQQNKDGI GEICLWGRHIFMGYLESETETTEAIDDEG
WLHSGDLGQLDGLGFLYVTGHIKEILITAGGENVPPIPVETLVKKKIPIISNAMLVGDCLKFLSMLLT LK
ASTVTEIVKQQDPLVYKAIQQGINAVNQEAMNNAQRIEKWVILEKDFSIYGGELGPMMLKRHFVAQKYK
KQIDHMYH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001289179

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001289179.2](#)

RefSeq Size: 2705 bp

RefSeq ORF: 1917 bp

Locus ID: 81616


UniProt ID: [Q5FVE4](#)

Cytogenetics: 19p13.3

MW: 71.7 kDa

Gene Summary: Mediates activation of long-chain fatty acids for both synthesis of cellular lipids, and degradation via beta-oxidation. Able to activate long-chain fatty acids. Also able to activate very long-chain fatty acids; however, the relevance of such activity is unclear in vivo. Has increased ability to activate oleic and linoleic acid. May play a role in spermatogenesis. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239176