

Product datasheet for **RC239190**

EYA4 (NM_001301013) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EYA4 (NM_001301013) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	EYA4
Synonyms:	CMD1J; DFNA10
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239190 representing NM_001301013
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAAGACTCCCAGGATTTAAATGAACAATCAGTAAAGAAAACGTGCACAGAATCAGATGTTTCACAAT
CTCAGAATTCAGGTCTATGGAAATGCAGGACCTAGCAAGTCTCATACTTGTGGAGGTGGTGATAC
TCCAGGTAGTCCAAACTGGAAAAATCTAATCTCAGCAGCACATCAGTTACTACAATGGGACAGGAGGG
GAAAAATGACTGTTTTAAACACAGCAGACTGGTTGCTGAGTTGCAACACCCCTCTTCTGCAACAATGT
CTCTTCTGCACTCAAACAGAGCCCTTGAACAGCAGTGAACCACAGCCAGACTGGAGATGGAGCGCT
TGACACTTTTACTGGGTCACTAATTACAAGTAGTGGCTACAGCCCAGATCAGCACATCAGTATCCCA
CAGCTGTATCCTTCCAAGCCCTATCCACACATTCTTCTACACCAGCAGCTCAAACAATGTCTGCATG
CAGGCCAGACTCAGTATTCGGGGATGCAGCAGCCAGCCGTCTACACAGCCTACTCACAGACAGGACGCC
CTACAGCTTGCCACTTACGATTTGGGTGTGATGTTGCCAGCCATCAAGACAGAGAGTGGACTTTCCAA
ACTCAGTCCCCATTACAGAGTGGCTGCCTCAGTTACAGCCCAGGGTTCTTACCCACAGCCAGGCCAGA
CACCTTATTCTTACCAAATGCCAGGTTCTAGTTTTGCACCATCATCTACTATTTATGCAAAATAATCAGT
TTCCAATTCACGAATTTCAAGTGGTTCACAACAGGATTATCCATCCTATACAGCCTTTGGCCAAAACAG
TATGCACAGTATTATTCAGCATCAACGTATGGAGCGTATATGACATCGAATAACACAGCCGATGGCACAC
CCTCTTCAACCTCTACTATCAGTTGCAGGAATCTCTCCAGGACTGACTAACCAACCAGGTACAGATCT
TCACCCAGGAGAGTTCGATACCATGCAGAGTCCCTCCACCCATCAAAGATCTTGATGAGAGAACCTGT
AGGAGTTCTGGTCAAAGTCCAGAGGAAGAGGCCGGAAAAATAATCCCTCCCGCCTCCTGATAGTGACC
TGGAGCGTGTGTTGTCTGGGATTTGGATGAAACCATCATTGTTTTTCACTCACTGCTCACCGGCTTAA
TGCACAGAAGTATGGCAAGGATCCCCCATGGCTGTAACCCTGGACTCCGCATGGAAGAAATGATTTTT
AATCTTGCTGATACTCATTGTTTTTAATGATTTAGAGGAGTGTGATCAAGTTCATATAGATGATGTTT
CCTCTGATGATAATGGGCAGGACTTAAGTACCTACAGTTTTGCAACTGATGGCTTCCATGCAGCTGCAAG
TAGTGCAAACCTTTGTTTCCCAACAGGTGAAGAGGAGGGTTGACTGGATGAGGAAGTTGGCTTTTCGT
TACAGAAGAGTAAAAGAATTATATAACACCTACAAGAACAACGTTGGAGGACTCCTTGGCCCTGCCAAGA
GGGATGCCTGGCTACAGTTAAGGGCAGAGATTGAAGGTCTGACAGATTCTGGCTAACAAATGCACTTAA
GTCTTTATCAATTATTAGCACTAGGAGTAAGTGCATAAATGTCTTGGTAACGACAACCTCACTGATCCCA
GCACTTGGCAAGTTCTACTCTATAGTTTAGGAGGTGCTTCCCCATTGAGAATATTACAGTGCACCTA
AAATAGGAAAAGAAAGTTGCTTTGAACGAATAATGCAAAGGTTTGGCAGAAAAGTAGTGTATGTTGTAAT
TGGGGATGGTGTAGAAGAAGAACAGGCAGCAAAAAGCACAACATGCCCTTCTGGAGGATATCCAGTCAC
TCAGACCTCCTGGCTCTCCACCAAGCACTGGAATTAGAGTATTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239190 representing NM_001301013
 Red=Cloning site Green=Tags(s)


MEDSQDLNEQSVKKTCTESDVSQSQNSRSMEMQDLASPHLTVGGGDTPGSSKLEKSNLSSTSSTVTTNGTGG
 ENMTVLNTADWLLSCNTPSSATMSLLAVKTEPLNSSETTATTGDGALDTFTGSVITSSGYSRPSAHQYSP
 QLYPSKPYPHILSTPAAQTMSAYAGQTQYSGMQQPAVYTAYSQTGQPYSLPTYDLGVMLPAIKTESGLSQ
 TQSPQLQSGCLSYSPGFSTPQPGQTPYSYQMPGSSFAPSSTIYANNSVSNSTNFGSQDQDYPSTAFGQNQ
 YAQYYSASTYGAYMTSNNTADGTPSSTSTYQLQEESLPLTNQPGTDLHPGEFDTMQSPSTPIKDLDERTC
 RSSGSKSRGRGRKNPSPPPDSLERVFWDLDETIIVFHSLLTGSYAQYKGDPPMAVTLGLRMEEMIF
 NLADTHLFFNDLEECDQVHIDDVSSDDNGQDLSTYSFATDGFHAAASSANLCLPTGVRGGVDWMRKLAFR
 YRRVKELYNTYKNNVGLLGPAKRDAWLQLRAEIEGLTDSWLTNALKSLSIISTRNCINLVTTTQLIP
 ALAKVLLYSLGGAFPIENIYSATKIGKESCFERIMQRFGRKVVVVVIGDGVVEEQAAKKNMPFWRISSH
 SDLLALHQALELEYL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:


ACCN: NM_001301013

ORF Size: 1935 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001301013.1](#), [NP_001287942.1](#)

RefSeq Size: 5573 bp

RefSeq ORF: 1938 bp

Locus ID: 2070

UniProt ID: [O95677](#)

Cytogenetics: 6q23.2

Protein Families: Druggable Genome, Phosphatase, Transcription Factors

MW: 70.6 kDa

Gene Summary: This gene encodes a member of the eyes absent (EYA) family of proteins. The encoded protein may act as a transcriptional activator through its protein phosphatase activity, and it may be important for eye development, and for continued function of the mature organ of Corti. Mutations in this gene are associated with postlingual, progressive, autosomal dominant hearing loss at the deafness, autosomal dominant non-syndromic sensorineural 10 locus. The encoded protein is also a putative oncogene that mediates DNA repair, apoptosis, and innate immunity following DNA damage, cellular damage, and viral attack. Defects in this gene are also associated with dilated cardiomyopathy 1J. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2014]