

Product datasheet for **RC239193**

RIPOR2 (NM_001286446) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RIPOR2 (NM_001286446) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RIPOR2
Synonyms:	C6orf32; DFNB104; DIFF40; DIFF48; FAM65B; MYONAP; PL48
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239193 representing NM_001286446
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGTGTACCTTCGCGATTACACTGGATCAGAAACAGGCCATGGAGGGACAGGATCCGGAGGCGCCAGC
 TCAACCGACTACCGACCAGACTCCCGAAATCATGTTGGTAGGATCCAGTCTTTTTCGCCTGGAGGGCC
 CAATGGGATCATTAGAAGCCAGTCCTTTGCGGGTTTCAGCGGCCTCCAGGAAAGCGATCCAGGTGAAC
 TCCTTCATTGAAAATTCCTCGCTCTCAAGAAGCCTCAGGCCAACTGAAGAAAATGCACAATTTAGGCC
 ACAAAAACAACATCCCCCAAAGAGCCTCAGCCTAAAAGGGTGAAGAAGTCTACAGGGCCTTGAAAAA
 TGGACTTGATGAATATCTGGAGGTTACCAGACGGAGCTGGACAAGTTGACAGCTCAGTTAAAGATATG
 AAAAGAACTCTCGCCTGGGTGTACTGTATGACCTAGACAAGCAAATTAACAATTGAAAGATACATGA
 GACGCCTGGAGTTTCATATAAGTAAGGTAGATGAACCTATGAAGCTTATTGTATCCAGCGACGCCTCCA
 GGATGGTGCCAGCAAAATGAAGCAAGCCTTCGCAACATCCCTGCCAGCAAAGTGCCCGGAGAGTCTG
 ACAGAGATCAATCGGAGCTTCAAGGAGTACACAGAGAATATGTGACCATTGAAGTGGAGCTAGAGAATC
 TGCTGGGAGAATTCTCCATCAAGATGAAAGGTCTGGCTGGCTTTGCACGCCTCTGTCTGGAGATCAATA
 TGAAATTTTATGAAGTATGGCCGGCAGCGGTGAAACTGAAAGGCAAAATAGAAGTAAATGGCAAGCAG
 AGCTGGGATGGAGAAGAAACAGTTTTCTGCCCTGATAGTTGGGTTCATCTCCATCAAGGTCACGGAGC
 TCAAAGGGCTAGCAACTCACATCCTGGTAGGTAGCGTACCTGTGAGACCAAAGAGCTGTTTGCAGCCCG
 ACCTCAGGTAGTGGCTGTGACATCAATGACCTTGGTACCATCAAACCTGAAATCACCTGGTAT
 CCATTTGACGTGGAGGACATGACCGCATCCTCAGGCGCTGGGAACAAGGCAGCAGCCCTTCAGAGGAGAA
 TGTCCATGTACAGCCAGGGTACCCCGAAACGCCACCTTCAAAGACCACTCCTCTTTTCAAATCTACC
 TGATGACATCTTTGAAAATGGAAGGCAGCCGAGGAGAAAATGCCACTGTGCTCAGCTTCAGTGACCTG
 CCCAACGGGGACTGCGCCTCACCTCCCACTCAACAGGCTCCCCTTCAAACCTCAACAAATCCAGAAATTA
 CCATCACCCCTGCGGAGTTTAACCTCAGCAGCTTGGCCTCCAGAATGAGGGTATGGATGACACCAGCTC
 AGCATCTTCCAGGAACCTCCCTGGGAGAAGGCCAAGAGCCAAAGTCACACCTGAAGGAGGAAGACCCAGAG
 GAGCCAGAAAACCTGCCTCGGCCCATCTGAGGCTTGCCGCCGACAGTCTCAGGTGCTGGGGCTGAGC
 ACCTGTTCTTGAGAATGATGTTGCAGAAGCACTTCTGCAAGAGTCTGAGGAGGCTCTGAGCTCAAGCC
 TGTGGAACCTGGACACTTCGGAAGGAAACATCACAAAGCAGCTGGTCAAGAGGCTCACATCTGCAGAGGTG
 CCAATGGCCACAGACAGGCTGCTCTCTGAGGGTTCTGTTGGTGGAGAATCTGAAGGCTGCAGATCCTTTC
 TAGATGGAAGCTTAGAGGATGCTTTTAATGGGCTTTTACTTGCATTAGAACCACATAAAGAGCAGTATAA
 AGAGTTTCAGGATCTGAACCAAGAAGTCATGAATTTGGATGATATTCTAAAAGATGCTAAACACCTTGAC
 GATCAGAAACTAAATGGTGCTGCCTCTTGGACGGAGATCACTGAGGGTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239193 representing NM_001286446
 Red=Cloning site Green=Tags(s)


MVLPSRLHWIRNRPWRDRIRRRQLNRLPTRLPEIMLVGSQSFSPGGPNGIIRSQSFAGFSGLQERRSRCN
 SFIENSALKKPQAKLKKMHNHGHKNNPPKEPQPKRVEEVYRALKNGLDEYLEVHQTELDKLTAKLKDM
 KRNSRLGVL YLDKQIKTIERYMRRLEFHISKVDELYEAYCIQRRQLQDGASKMKQAFATSPASKAARESL
 TEINRSFKEYTENMCTIEVELENLLGEFSIKMKGLAGFARLCPGDQYEIFMKYGRQRWKLKGKIEVNGKQ
 SWDGEETVFLPLIVGFIISKVTELKGLATHILVGSVTCETKELFAARPQVVAVDINDLGTIKLNLEITWY
 PFDVEDMTASSGAGNKAALQRRMSMYSQGTPTPTFKDHSFFSNLPDDIFENGKAAEEKMPLSLSFSDL
 PNGDCALTSHTGSPSNSTNPEITITPAEFNLSSLASQNEGMDTSSASSRNSLGEGQEPKSHLKEEDPE
 EPRKPASAPSEACRRQSSGAGAEHLFLENDVAEALLQESEEASELKPVELDTSEGNITKQLVKRLTSAEV
 PMATDRLLSEGSVGGESEGRSFLDGSLEDAFNGLLLALEPHKEQYKEFDLNQEVMLNDDILKDAKHL
 DQKLNGAASWTEITEGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001286446

ORF Size: 1941 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001286446.2](#)

RefSeq Size: 3541 bp

RefSeq ORF: 1944 bp

Locus ID: 9750


UniProt ID: [Q9Y4F9](#)

Cytogenetics: 6p22.3

MW: 72.9 kDa

Gene Summary: This gene encodes an atypical inhibitor of the small G protein RhoA. Inhibition of RhoA activity by the encoded protein mediates myoblast fusion and polarization of T cells and neutrophils. The encoded protein is a component of hair cell stereocilia that is essential for hearing. A splice site mutation in this gene results in hearing loss in human patients. [provided by RefSeq, Sep 2016]

Product images:



Circular map for RC239193