

Product datasheet for **RC239196**

ZNF202 (NM_001301779) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ZNF202
Synonyms:	ZKSCAN10; ZSCAN42
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC239196 representing NM_001301779
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCTACAGCCGTGGAACAGAGGACCAGGATCTTTGGGAAGAAGAGGAATTCTGATGGTGAACCTGG
 AAGATGATTTACCTGTGCGCCAGAGTCTGTCTTACAGAGGGATGACCCGGTGTGGAACCTCCACCA
 GAATTCGACGCTTCCGCTACCAGGAGGAGCAAGCCCTAGAGAAGCTCTCATCAGACTCCGAGAACTT
 TGTCACCACTGGTGTAGAGGAGGCGGACAAAGGAGCAGATCTAGAGTGTCTGTGCTGGAACAAT
 TTCTTACCGTCTACCTGGAGAATACAGAGCTGGGTGCGGGGCCAACGGCCAGAAAGTGGCAGGAGGC
 AGTGACGCTGGTGGAGGGTTTGCAGAAACAACCCAGGAGACCAAGCGGTGGGTGACTGTCCATGTTTAC
 GGCCAGGAAGTCTGTGAGAGGAGCGGTGCATTTAGGAGCGGAGCCTGAGTCACTAATGAGCTGCAGG
 ATCCTGTGCAAGCTCGACCCCGAGCAGTCTCCTGAGGAAACACACAGAGCCAGATCTGGGGCACC
 GGCAGAGCAGCGTCCACACCAGGAAGAGGAGCTCCAGACCTGCAGGAGAGCGAGGTCCAGTCCCGAG
 GACCCAGACCTTCTGCAGAGAGGAGCTCTGGAGACTCAGAGATGGTTGCTCTTCTTACTGCTCTGTAC
 AGGGACTGGTAACGTTCAAGGATGTGGCCGTATGCTTTTCCAGGACCAGTGGAGTGTCTGGACCAAC
 ACAGAAAGAGTTCTATGGAGAATATGTCTTGAAGAAGACTGTGGAATTGTTGTCTCTGTCAATTTCCA
 ATCCCCAGACCTGATGAGATCTCCAGGTTAGAGAGGAAGACCTTGGGTCCCAGATATCCAAGAGCCTC
 AGGAGACTCAAGAGCCAGAAATCCTGAGTTTACCTACACAGGAGATAGGAGTAAAGATGAGGAAGAGTG
 TCTGGAGCAGGAAGATCTGAGTTTGGAGGATATACACAGGCCTGTTTGGGAGAACCAGAAATTCACCAG
 ACTCCAGATTGGGAAATAGTCTTTGAGGACAATCCAGGTAGACTTAATGAAAGAAGATTGGTACTAATA
 TTTCTCAAGTGAATAGTTTGTGAACCTTCGGGAAACTACACCCGTCCACCCCTGTTAGGGAGGCATCA
 TGACTGTTCTGTGTGTGGAAGAGCTTCACTTGTAACTCCACCTTGTAGACACCTGAGGACTCACACA
 GGAGAGAAACCTATAATGTATGGAATGTGGAAGGTTACACACGAAGCTCACATCTTGCAGGCACC
 AAAAGGTTCAAGATGAACGCGCCTTACAAATATCCCTAAACCGGAAGAATTTGGAAGAGACCTCCCC
 TGTGACACAGGCTGAGAGAATCCATCAGTGGAGAAACCTATAGATGTGATGATTGCGGAAGCACTTC
 CGCTGGACTTCAGACCTTGTGAGACATCAGAGGACACATACTGGAGAAAAACCTTCTTTTGTACTATT
 GTGGCAAAAGCTTCAGCCAGAAATCTGTGTTAAACAACACACCAAGAATCCACCTGGGAGGCAAAACCTA
 CTTGTGTGGAGAGTGTGGTGGAGACTTCAGTGAACACAGCGGTACCTGGCGCACCGGAAGACGCACGCT
 GCTGAGGAATCTACCTCTGCAGCGAGTGGGGCGCTGCTTACCCACAGCGCAGCGTTCGCCAAGCACT
 TGAGAGGACACGCCTCAGTGGGCGCTGCCGATGCAACGAATGTGGGAAGAGCTTCACTCGCAGGGACCA
 CCTCGTCAGGCATCAGAGAACACACTGGGGAGAAACCATTCACGTGCCCTACCTGTGAAAAAGCTTC
 AGCAGAGGATATCACTTAATTAGGCATCAGAGGACCCACTCAGAAAAGACCTCC

ACGCGTACGCGGCCGTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239196 representing NM_001301779
 Red=Cloning site Green=Tags(s)


MATAVEPEDQDLWEEGILMVKLEDDFTCRPESVLQRDDPVLETSHQNFRFRFYQEASPREALIRLREL
 CHQWLRPERRTKEQILELLVLEQFLTVLPGELQSWVRGQRPESGEEAVTLVEGLQKQPRRPRRWTVHVVH
 GQEVLSSEETVHLGAEPESPNELODPVQSSTPEQSPETTQSPDLGAPAEQRPHQEELQTLQESEVPVPE
 DPDLPAERSSGDSEMVALLTALSQGLVTFKDVAVCFSDQWSDLDPTQKEFYGEYVLEEDCGIVVLSFP
 IPRPDEISQVREEEPWVPDIQEPQETQEPEILSFTYTGDRSKDEEECLEQEDLSLEDIHRPVLGEPEIHQ
 TPDWEIVFEDNPGRLNERRFGTNIISQVNSFVNLRETTVPVHLLGRHHDCSVCGKSFTCNHSLVRHLRTH
 GEKPYKMECGKSYTRSSHLARHQVHKMNAPYKYPLNRKNLEETSPVTQAERTPSVEKPYRCDDCGKHF
 RWTSDLVRHQRTHTGEKPFCTICGKSFSQKSVLTHQRIHLGGKPYLCEGCEGDFSEHRRYLAKRTHA
 AEELYLCSECGRCFTHSAAFAKHLRGHASVPRCRCNECGKSFRRDHLVRHQRTHTGEKPFCTCPTCGKSF
 SRGYHLIRHQRTHTSEKTS

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6153_c03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001301779

ORF Size: 1944 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001301779.1, NP_001288708.1

RefSeq Size: 4393 bp

RefSeq ORF: 1947 bp

Locus ID: 7753

UniProt ID: O95125


Cytogenetics: 11q24.1

Protein Families: Transcription Factors

MW: 74.7 kDa

Gene Summary: Transcriptional repressor that binds to elements found predominantly in genes that participate in lipid metabolism. Among its targets are structural components of lipoprotein particles (apolipoproteins AIV, CIII, and E), enzymes involved in lipid processing (lipoprotein lipase, lecithin cholesteryl ester transferase), transporters involved in lipid homeostasis (ABCA1, ABCG1), and several genes involved in processes related to energy metabolism and vascular disease.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239196