

Product datasheet for **RC239198**

ALS2CR8 (CARF) (NM_001282910) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALS2CR8 (CARF) (NM_001282910) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ALS2CR8
Synonyms:	ALS2CR8; NYD-SP24
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239198 representing NM_001282910
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGAACAATCTAATGATTCTTAAGAGTCAACCATAATGACGGTGAAGAGTCAAAAACCAAGTGCTCAAG
 TATTTGAGATGATCGTTGCCAGCCCAACAGAAAATGGACAGGTACTTCGTGTAATTCACCTACCCAGAC
 AGGAATGGCACAAGTGATTATACCTCAGGGGCAACTTGTGGATGTGAATAGTCCTCGGGATGTCCCTGAA
 GAGAAACCCAGTAACAGAACTTACCAACTGTAAGAGTGGATACTCTAGCAGACAATACCAGCAATTACA
 TTCTTCATCTCAACATCCTTCCATTGCCCAAAAAGTCAGTGACCGGAATGCTGGAAGAACCCCTTCT
 GGGGCTCTTCAGCCACTTCTTCTAATACACCTATATGGGCTGCCGTCTTAGGAGCTGTGAGAAAATT
 GGAGATTCATACCGTGGCTACTGTGAAGTGAGACTGAATTAGAAAGTGCTAACATTTACAAGCAGC
 AAACACAGAGTGTTGGGGACCCGTAGTCTCCAAGCCAGCCAAGCCTGTACACGCTTGATGTGGAA
 ATCCCAGTATGTTCCATATGATGGAATCCATTTGTTAATGCAGGGAGTAGAGCTGTGGTAATGGAGTGT
 CAGTATGGGCCAAGAAGAAAAGGTTTCCAGTTAAAAAAGTCAGTGAGCAGGAAAGCAGGTCTTGTGAGC
 TCTACAAAGCCACTTGTCCAGCTCGGATTACATTA AAAAGGTACAGAAGTTCTGAAATATAGATTCC
 TACAGACCCCAAAATTGACAAGAAAATTATCAGAAATGGAGCAGGAGAAAGCTTTTAACATGCTAAAGAAG
 AACTTGGTAGATGCTGGTGGTGTCTTAGGTGGTATGTACAGTTACCTACACAGCAAGCTCATCAGTATC
 ATGAATTAGAGACTCCCTGCCTCACTTTGTACCTTCTCTTTTCTGTGCTTCTCTTGAAGAAGAGGA
 AACTGCAGTTAGAGATGAGAATTGTGCATTACCCTCACGTTTACATCCTCAAGTAGCACATAAGATTCAA
 GAATTAGTATCACAGGAATAGAACAAGTGTATGCAGTAAGGAAACAGCTAAGAAAATTTGTGGAAGGG
 AACTGTTCAAACCCGATGAGGTACCTGAAAGACATAATTTATCTTTTTTCCAAGTGAATGATATAAA
 AAATCACATCCATGAGGTACAGAAATCCTTGAGAAATGGAGATACGGTATATAACTCAGAGATTATTCCA
 GCAACGCTTCAATGGAAGTACAGACAGTGGGAATATTCTCAAAGAGACCATGACAGTTACATTTGCAGAAG
 GAAATTCACCAGGAGAATCAATTACCACCAAAGTGAAACAAATCAGACCAGGGGTTCTTTGTCTCCTGA
 GCCAACCACTTGTCTCCTCACTCTCCTCATTTCAGCCCAAAATATTTACAACTACAGGGTTTGCAG
 TTACAACCAAGGTACACCTCTCCTGATGAATCACCAGCTGTGGTATCAGTAAATAACCAGCCGCTCTCTA
 GTCCTTCAGGACTTCTGGATACAATAGGAAGTGTGTAATGAATAATAATTCTCTACTGCTTGGTCAAAG
 TCATAGCCTTCAAAGAGATACATGCTTAACCCAAAACAATAGTACTGCCTCCACCATGGGTAACCTTCCA
 GAACCAGATCAAAATCTAGTTGCAATGGACGAGCTGGTAGAAGTTGGAGATGTTGAGGATACAGGGAATC
 TGGAAGGAACTGTTTCATCGGATTCTGTTGGGAGATGTGCAGACTATTCCAATACAGATTATAGACAACCA
 CTCAGCTCTTATTGAAGAAAATCCAGAAAGTACCATTTCTGTGAGCCAAGTTAAACAAGAACCAAGAA
 CCAGCATTGTCTATGGAAGCAAAAAAAGTGTGGACTATAAGAAATTATCTGCTACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239198 representing NM_001282910
 Red=Cloning site Green=Tags(s)


MEQSNDSLRVNHNDGEESKTSAQVFEMIVASPTENGQVLRVIPPTQTGMAQVVIIPQGQLVDVNSPRDVPE
 EKPSNRNLPTRVDTLADNTSNYLHPQTSFPLPKKSVTGMLEELLGPLQLSSNTPIWACRLRSCEKI
 GDSYRGYCVSETELESVLTFFHKQQTQSVWGTRQSPSPAKPATRLMWKSQYVPYDGIPFVNAGSRVVMCE
 QYGPRRGFQLKKVSEQESRSCQLYKATCPARIYIKKVQKFPEYRVPTDPKIDKKIIRMEQEKAFNMLKK
 NLVDAGGVLRWYVQLPTQQAHQYHELETPCLTLSPSPFPVSSLEEEETAVRDENCALPSRLHPQVAHKIQ
 ELVSQGIEQVYAVRKQLRKFKVERELFKPDEVPERHNLFFPTVNDIKNHIHEVQKSLRNGDTVYNSIIP
 ATLQWTTDSGNILKETMTVTFAEGNSPGESITTKVETNQTRGSLSPEPTHLLSSLSFQPKIFTQLQGLQ
 LQPRYTSPDESPAVSVNNQPSPPSGLDITGSAMNNNSLLLGQSHSLQRDTCLTQNNSTASTMGNL
 EPDQNLVAMDELVEVGDEDTGNLEGTVHRILLGDVQTIPIQIIDNHSALIEENPESTISVSQVKQEPKE
 PALSMEAKKTVDYKKLSAT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001282910

ORF Size: 1947 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001282910.3](#)

RefSeq Size: 5128 bp

RefSeq ORF: 1950 bp

Locus ID: 79800


UniProt ID: [Q8N187](#)

Cytogenetics: 2q33.2

MW: 72.8 kDa

Gene Summary: Acts as a transcriptional activator that mediates the calcium- and neuron-selective induction of BDNF exon III transcription. Binds to the consensus calcium-response element CaRE1 5'-CTATTTCGAG-3' sequence.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239198