

Product datasheet for **RC239244**

LIN54 (NM_001288996) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIN54 (NM_001288996) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LIN54
Synonyms:	CXCDC1; JC8.6; MIP120; TCX1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239244 representing NM_001288996
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGAGGTGGTGCCAGCTGAGGTGAATAGTTTGCTTCCAGAGGAAATAATGGACACTGGTATAACTTTAG
 TGGATGATGATAGTATTGAGGCTGTTATTGTTTCATCCCAATTCCCATGGAGACAGAAGTGAAGAAAT
 TGCAACATAAATTCTACTGGTACTCTACAGCCACGCCATTTCCACGGAACCAATCACAGTGACAGT
 AACCACTAACCAAGTTGCAGTGAATACCACAATTACTAAAGCAGATTCTAATACCACAGTGAAACAG
 CTTTTCAGTGGCCTTCAAAAACCTGGTGCTCAGACTCCTGTGACTATATCAGCCAATCAGATTATTTT
 AAACAAAGTATCACAGACATCTGATCTTAACTTGGAATCAGACCTTAAACCAGATGGACAGAAGTTA
 ATTTTAACTTTGGGCAAGTCTGGTTCACCAATTGTTTAGCACTACCCATAGCCAATACCCAGG
 CTCAGAAAGTTACAACCTAGGCCAGTCAGGAGATGCTAAGTTACCACCGCAGCAAAATTAAGTAGTTAC
 CATTGGAGGGAGGCCAGAGGTGAAACCTGTCATTGGTGTCTCAGCATTGACCCAGGAAGTCAACTGATT
 AATACTACAACCTAGCCCTCTGTGTTACAGACCAACAGTTAAAAACAGTACAGGCAGTGAATCAACTG
 TGCAGACCATCACTGTTGGAGGAGTGAGCACATCACAGTTAAGACAATTATTCCTCTGGCAACTGCTCC
 CAATGTCCAGCAGATTCAAGTGCCTGGAAGCAAGTTTCATTATGTCCGACTTGTTACTGCCACATCAGCC
 AGTAGCTCAACCCAGCCAGTTAGTCAGAATCCAGTACAAACACTCAGCCTCTTCAGCAAGCAAAGCCAG
 TGGTTGTTAATACAACCCAGTGCAGTGTGAGTTCCAATTGTCTCAGCTCAGGCTGTCAAACAAGTTGT
 TCCAAAACCAATCAATCCAATTCACAAATAGTAACTACTAGCCAGCCACAGCAACGGCTTATCATGCCCT
 GCCACACCACTGCCACAGATCCAGCCCAACCTCACTAACCTGCCACCAGGCACTGCTCTGGCACCAGCTC
 CGGGAACAGGGAATGTGGGTTATGCAGTGCTTCCAGCTCAGTATGTTACTCAGCTACAGCAGTCTTCATA
 TGTATCAATAGCAAGCAACTCTACCTTTACTGGAACATCTGGTATCCAGACCCAGGCACGGCTTCCATT
 AATGGCATAATCCCATCAGAGTCGGCCAGTCGGCCCCGAAAGCCCTGTAATTGTACAAAATCACTGTGTT
 TGAAATTGTATTGTGATTGCTTTGCAAATGGTGAATTTTGCAACAAGTGAATTTGTTACTAATTGTTACAA
 CAATTTGGAACATGAAAATGAAAGGCAAAAAGCAATAAAGGCATGCCTTGACAGAAATCCAGAAGCCTTT
 AAGCCTAAGATAGGGAAGGAAAGGAGGAGAGTCTGATCGACGTCATAGCAAAGGTGTAATTGCAAAAC
 GATCAGGATGTCTTAAAACTACTGTGAATGCTATGAGGCAAAAATAATGTGTTCTCAATATGCAAAATG
 TATTGGCTGTAAGAATTTTGAAGAAAGCCGGAAGGAAGACATTGATGCATTTGGCAGATGCAGCTGAA
 GTAAGGGTACAGCAACAAACAGCAGCCAAGACGAAGTTATCCTCTCAAATTTCAGACTTGCTTACTAGGC
 CAACACCAGCTTTAAATAGTGGAGGCGGAAAATTGCCATTTACATTTGTAAGTAAGGAAGTAGCTGAAGC
 CACATGTAATTGCTCCTTGCCAGGCAGAGCAGGCAGACAAGAAGGGAATCAAGGCAGCAGCGGAA
 CGGATGATACTTGAGGAATTCGGACGATGTTTGATGAGTGTCACTCACTCTGCAGGAAAGGCAAAAAGTG
 ACCCTTGTGCCATGAATTGC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239244 representing NM_001288996
Red=Cloning site Green=Tags(s)

MEVVPAEVNSLLPEEIMDTGITLVDDDSIEAVIVSSPIMETELEEEIVNINSTGDSTATPISTEPITVYS
NHTNQVAVNTTITKADSNTTVKPAFPSGLQKLGAQTPVTISANQIILNKVSQTSDLKLGNTLKPDPGQKL
ILTTLGKSGSPIVLALPHSQLPQAQKVTTQAQSGDAKLPPQIQIKVVTIGGRPEVKPVIGVSALTPGSQLI
NTTTQPSVLQTQQLKTVQAVKSTVQTITVGGVSTSQFKTIIPLATAPNVQQIQVPGSKFHYVRLVTATSA
SSSTQPSVQNPNSTNTQPLQQAQPVVVNTTPVRMSVPIVSAQAVKQVVPKPINPTSQIVTTSQPQQRLIMP
ATPLPQIQPNLTNLPPGTVLAPAPGTGNVGYAVLPAQYVTLQSSYVSIASNSTFTGTSIGIQTQARLPF
NGIIPSESASRPKPCNCTKSLCLKLYCDFANGFENNCNCTNCTNLEHENERQKAIAKLDRNPEAF
KPKIGKGESESDRRHSGKCNCKRSGCLKNYCECYEAKIMCSSICKIGCKNFESPERKTLMLADAAE
VRVQQQTAAKTKLSSQISDLLTRPTPALNSGGGKLPTFTVTKVEAEATCNCLLAQAEQADKKGSKAAAE
RMILEEFGRCLMSVINSAGKAKSDPCAMNC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001288996

ORF Size: 1980 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001288996.2](#)

RefSeq Size: 5876 bp

RefSeq ORF: 1983 bp

Locus ID: 132660


UniProt ID: [Q6MZP7](#)

Cytogenetics: 4q21.22

MW: 70.8 kDa

Gene Summary: LIN54 is a component of the LIN, or DREAM, complex, an essential regulator of cell cycle genes (Schmit et al., 2009 [PubMed 19725879]).[supplied by OMIM, Dec 2010]

Product images:



Circular map for RC239244