

## Product datasheet for **RC239253**

### **HJURP (NM\_001282963) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	HJURP (NM_001282963) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	HJURP
Synonyms:	FAKTS; hFLEG1; URLC9
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC239253 representing NM\_001282963  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCTGGGTACGCTGCGCGCCATGGAGGGCGAGGACGTGGAAGACGACCAGCTGCTGCAGAAGCTCAGGG  
CCAGTCGCGCGCCCTTCCAGAGGCGCATGCAGCGGCTGATAGAGAAGTACAACCAGCCCTTCGAGGACAC  
CCCGGTGGTGCAAATGGCCACGCTGACCTACGAGACGCCACAGGGATTGAGAATTTGGGGTGAAGACTA  
ATAAAGGAAAGAAACGAAGGAGAGATCCAGTGTGCAGGTAACAGAGCTGGAAGGGATGTACGTGTGACTC  
CGCTGCCCTTCACTGGCCTCACCTGCCGTGCCTGCCCGGATACTGCAGTCGTATCTCCAGAAAGAGTCC  
TGGTGACCCAGCGAAACCAGCTTCATCTCCAGAGAATGGGATCCTTTGCATCCTTCCACAGACATG  
GCCTTAGTACCTAGAAATGACAGCCTCTCCCTACAAGAGACCAGTAGCAGCAGCTTCTTAAGCAGCCAGC  
CCTTTGAAGATGATGACATTTGCAATGTGACCATCAGTGACCTGTACGCAGGGATGCTGCACTCCATGAG  
CCGGCTGTTGAGCACAAAGCCATCAAGCATCATCTCCACAAAACGTTTCATCATGAAAACCTGGAACCTC  
AGGAGGAGGCACAGATATAAGAGCAGGATGAACAAAACATATTGCAAAGGAGCCAGACGTTCTCAGAGGA  
GCTCCAAGGAGAACTTCATACCTGCTCTGAGCCTGTGAAAGGGACAGGGGCATTAAGAGATTGCAAGAA  
CGTATTAGATGTTTTCTTGCCGTAAGACAGGTTTAAAAATTGGAAAAGCTTTTCTTGAAGTCAACAGACCC  
CAAATCCATAAGTTAGATCCAAGTTGGAAGGAGCGCAAAGTGACACCCTCGAAGTATTCTTCTTGATTT  
ACTTCGACTCCAGTGCAACATATAATCTTGATGAGGAAAATAGATTTAGGACATTAATAATGGTTAATTTT  
TCCTGTAAAAATAGTTTCCAGACCAACAATACGACAGGGCCATGGAGAGAACCCTCAGAGGGAGATTGAA  
ATCCGATTTGATCAGCTTCATCGGGAATATTGCCTGAGTCCAGGAACCAGCCTCGCCGGATGTGCCTCC  
CGGACTCCTGGCCATGAACATGTACAGAGGGGGTCTGCGAGTCTGGTGGCCTCAGGGCTTAGAAAC  
CCGCAGGCTGAGTTTACCTTCCAGCAAAGCAAAGCAAAGTTTAAAGTGAAGCTTTTGAACCTAGGC  
AAAAGATCTCTGGAAGCAGGTAGGTGCCTGCCAAGAGCGATTTCATCTTCACTTCCAAAGACCAACC  
CCACACACAGCGCAACTCGCCCGCAGCAGACATCTGACCTTACGTTTCAGGGAAATAGTTCTGGAATATT  
TAGAAAGTCAGTGTACCCAGCAAAACTCTTTCAGTCCCAGATAAAGAAGTGCCAGGCCACGGAAGGAAT  
CGTTACGATGAAATTAAGAAGAATTTGACAAGCTTCATCAAAGTATTGCCTCAAATCTCCTGGGCAGA  
TGACAGTGCCTTTATGTATTGGAGTGTCTACAGATAAAGCAAGTATGGAAGTTCGATATCAAACAGAAGG  
CTTCTTAGGAAAATTAATCCAGACCCTCACTTCCAGGGTTTCCAGAAGTTGCCATCATACCCCTGGGG  
TGCAGAAAAAGTCTACTGGGCTCAACTGCAATTGAGGCTCCTTACATCTACATGTGTGTCTGCGTCCATCA  
CGAGGGATGGCAGGAGGACCATCAGTTCCTGCAAAAAGACCCAGGCTATCAGAACCCAGGGCTCCGG  
ACGCCAGGGCAATTCCTGGGTGCCTCAGATGGGTGGACAACCCGTCAGACCCGGAGACCAGGGCAGC  
TCTTACAGCCCAACTCAGAAGAGAGAGGAGAGAACACGCTTACAGGATGGAAGAGAAAAGTGATTTCA  
TGCTAGAAAAATTGGAACCTAAAAGTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC239253 representing NM\_001282963  
 Red=Cloning site Green=Tags(s)


MLGTLRAMEGEDVEDDQLLQKLASRRRFQRRMQRLIEKYNQPFEDTPVVQMATLTYETPQGLRIWGGRL  
 IKERNEGEIQCAGNRAGRDRVVTPLPSLASPAPVAPGYCSRI SRKSPGDPAPASSPREWDPLHPSSTDM  
 ALVPRNDSLQETSSSSFLSSQPFEDDDICNVTISDL YAGMLHMSRLLSTKPSIIISTKTFIMQNWNS  
 RRRHRYKSRMNKTYCKGARRSQRSSKENFIPCSEPVKGTGALRDCKNVL DVSCRKTGLKLEKAFLEVNRP  
 QIHKLDPSWKERKVTPSKYSSL IYFDSSATYNLDEENRRTLKWLISPVKIVSRPTIRQGHGENRQREIE  
 IRFDQLHREYCLSPRNQPRRMCLPDSWAMNMYRGGPASPGLQGLETRRLSLPSSKAKAKSLSEAFENLG  
 KRSLEAGRCLPKSDSSSLPKTNPTHSA TRPQQTSDLHVQGNSSGIFRKS VSPSKTLSVPDKEVPGHGRN  
 RYDEIKEEFDKHLHQYCLKSPGQMTVPLCIGVSTDKASMEVRYQTEGFLGKLNPDPHFQGFQKLPSSPLG  
 CRKSLGSAIEAPSSTCVARAITRDGTRDHQFPAKRPRLSEPQGSGRQNSLGASDGDVNTVRPGDQGS  
 SSQPNSEERGENTS YRMEEKSDFMLEKLETKSV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001282963

**ORF Size:** 1989 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001282963.1](#), [NP\\_001269892.1](#)

RefSeq Size: 2934 bp

RefSeq ORF: 1992 bp

Locus ID: 55355

UniProt ID: [Q8NCD3](#)

Cytogenetics: 2q37.1

Protein Families: Druggable Genome

MW: 74.8 kDa

**Gene Summary:** Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres. Acts as a specific chaperone for CENPA and is required for the incorporation of newly synthesized CENPA molecules into nucleosomes at replicated centromeres. Prevents CENPA-H4 tetramerization and prevents premature DNA binding by the CENPA-H4 tetramer. Directly binds Holliday junctions.[UniProtKB/Swiss-Prot Function]