

Product datasheet for **RC239258**

ARMC9 (NM_001291656) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARMC9 (NM_001291656) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ARMC9
Synonyms:	ARM; JBTS30; KU-MEL-1; NS21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239258 representing NM_001291656
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGGACATTCTGGCTCATGAATCTGAATTACTTGGACTAGTGAAAGAGTATTTAGATTTTGCTGAAT
 TTGAAGACACCTTGAAAACATTTTCAAAAGAATGCAAAATAAAAGGAAAACCATTTGTGTAAACAGTAGG
 CGGATCTTTCAGAGACTCCAATCATTGACAATTCAGAAGGATCTTGTGCTGCAATTTGACAACGGAGAC
 CAGAAGGTGTTCTTCGATCTGTGGGAGGAGCACATTTCAAGTTCCATCCGAGATGGGGACTCCTTTGCC
 AGAAGCTGGAATTCTATCTCCACATCCATTTGCCATCTATCTTTGAAGTACTCTGTGGGAGACCGGA
 CAAAGAGGAGCTGGATGAAAAGATTTCTACTTCAAAACCTACCTGGAGACCAAAGGGGAGCCTTGAGC
 CAGACCACAGAGTTTCTCCTTTCTATGCCCTTCTTTTGTCCCAACCTATGGTGCACCCCTCATTTA
 AAGAACTCTTCCAGGATTCCTGGACTCCAGAGTTAAAGTTGAAGTTGAAAAGTTTCTAGCTTTAATATC
 TAAAGCCAGCAACACGCCAAAGCTTTTAACAATATATAAGGAGAATGGACAAAGTAACAAAGAAATCTTG
 CAGCAGCTCCACCAGCAGCTGGTTGAAGCTGAACGTAGGTCACTGACATACCTCAACCGGTACAATAAGA
 TCCAGGCCGACTACCACAATCTCATTGGAGTCACAGCAGAGCTGGTGGATTCTCTAGAGGCCACAGTCAG
 CGGCAAGATGATACCCCTGAGTACCTCCAGAGCGTCTGTGTCCGCTGTTCAAGTAAACAGATGCGGCAG
 AGCCTGGCGCATAGTGTGGACTTCACGAGGCTGGGACGGCATCCACCATGTTACGAGCCTCCTTGGCAC
 CCGTGAATTTGAAGGATGTCCTTACTGCCCTCCTTGGATTATGAGAACTGAAGAAGGATTTGATTTT
 GGGGAGTGACCGCTTGAAGCCTTCTTGTGACGGCTCTGCGCTGGCGCTTGACCACATCCCATCCTGGA
 GAGCAGAGGGAGACCGTTCTGCAAGCCTACATCAGCAATGACCTCTTGGACTGTTATAGCCACAACCAGA
 GGAGTGTGCTTCAGTTGCTGCACTCCACGAGCGACGTGGTGGCAGTACATGGCCAGGCTCATCAATGC
 TTTTGCCTCACTGGCAGAAGGTCGCTCTACCTTGCCCAAGAACAAAGGTGCTGCAGATGCTGGAGGGA
 AGGCTGAAGGAGGAGGACAAGGATATCATCACCAGGGAGAATGTTCTTGGGGCCCTGCAGAAGTTAGTC
 TCAGGCGCCCGCTGCAGACAGCGATGATTCAAGACGGCCTCATCTTCTGGCTGGTTGATGTTCTGAAGGA
 CCCTGACTGCCTGTCTGACTACACGCTGGAGTACTCGGTGGCTTTGCTCATGAACCTCTGCCTCCGACG
 ACAGGGAAGAATGTGTGCCAAGGTGGCAGGCTCGTGCTCAAAGTCTTTGGATCTTCTTGGCCATG
 AAAACCATGAGATACAGCCGTATGTGAATGGAGCTCTGTACAGCATCCTTTCTGTTCCATCCATTCTGTA
 GGAAGCAAGAGCAATGGGAATGGAAGACATCCTACGCTGCTTCATCAAAGAAGGCAATGCTGAAATGATC
 CGCCAGATAGAATTCATCATCAAGCAGCTAAATCCGAAGAGCTACCAGATGGTGTCTTGAATCTGATG
 ATGATGAAGATGAAGATGATGAAGAGGACCATGACATCATGGAAGCCGATCTGGACAAAGACGAATGAT
 CCAGCCCCAGCTCGGAGAACTCTCAGGAGAGAAGCTTCTGACCACGGAGTACCTGGGGATCATGACCAAC
 ACGGGGAAGACAAGGCGGAAGGGGCTGGCTAATGTGCAGTGGAGCGGGATGAGCCCCTGCAAAGGCCCG
 TCACCCCGCGGCCACAGAAACGGGTACCCAGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239258 representing NM_001291656
Red=Cloning site Green=Tags(s)

MGDILAHASELLGLVKEYLDFAEFEDTLKTSKECKIKGKPLCKTVGGSFRDSKSLTIQKDLVAAFNDNG
QKVFFDLWEEHISSIRDGDSFAQKLEFYLIHFAIYLLKYSVGRPDKEELDEKISYFKTYLETKGAALS
QTTEFLPFYALPFVPNPMVHPSFKELFQDSWTPELKLKLEKFLALISKASNTPKLLTIYKENGQSNKEIL
QQLHQQLVEARRSVTYLKRYNKIQADYHNLIGVTAELVDSLEATVSGKMITPEYLQSVCVRLFNSQMRQ
SLAHSVDFTRPGTASTMLRASLAPVKLKDVPLPSLDYELKKDLILGSDRLKAFLLQALRWRLTTSHPG
EQRETVLQAYISNDLLDCYSHNQSVLQLLHSTDVVRQYMARLINAFAASLAEGRLYLAQNTKVLQMLEG
RLKEEDKDIITRENVLQALQKFSRLRPQTAMIQDGLIFWLVDVLKDPDCLSDYTL EYSVALLMNLCLRS
TGKNMCAKVAGLVKVLSDLLGHENHEIQPYVNGALYSILSVPSIREEARAMGMEDIILRCFIKEGNAEMI
RQIEFIIKQLNSEELPDGVLSEDDDEDEDDEEDHIMEADLKDDELIQPQLGELSGEKLTTTEYLGIMTN
TGKTRRKGLANVOWSGDEPLORPVTGGHNRGYPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001291656

ORF Size: 1995 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001291656.1](#), [NP_001278585.1](#)

RefSeq Size: 2315 bp

RefSeq ORF: 1998 bp

Locus ID: 80210


UniProt ID: [Q7Z3E5](#)

Cytogenetics: 2q37.1

MW: 76 kDa

Gene Summary: Required for ciliogenesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239258