

Product datasheet for **RC239263**

Calpastatin (CAST) (NM_001284212) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calpastatin (CAST) (NM_001284212) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CAST
Synonyms:	BS-17; PLACK
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide
Sequence:

>RC239263 representing NM_001284212
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAATCCCACAGAAACCAAGGCTGTA AAAACAGAACCTGAGAAGAAGTCACAGTCAACCAAGCCAAAA
GCCTACCCAAGCAGGCATCAGATACAGGAAGTAACGATGCTCACAATAAAAAAGCAGTTTCCAGATCAGC
TGAACAGCAGCCATCAGAGAAATCAACAGAACCAAGACTAAACCAAGACATGATTTCTGCTGGTGGG
GAGAGTGTGCTGGTATCACTGCAATATCTGGCAAGCCGGGTGACAAGAAAAAGAAAAGAAATCATTAA
CCCCAGCTGTGCCAGTTGAATCTAAACCGGATAAACCATCGGGAAAGTCAGGCATGGATGCTGCTTTGGA
TGACTTAATAGATACTTTAGGAGGACCTGAAGAACTGAAGAAGAAAATACAACGTATACTGGACCAGAA
GTTTCAGATCCAATGAGTTCACCTACATAGAGGAATTGGGTAAAAGAGAAGTACAATTCCTCCAAAAT
ATAGGGAATATTGGCTAAAAGGAAGGGATCACAGGGCCTCCTGCAGACTCTTCGAAACCCATAGGGCC
AGATGATGCTATAGACGCTTGTCTGACTTCACCTGTGGGTGCCTACAGCTGCTGGAAAGAAAAGT
GAAAAAGAGGAATCTACAGAAGTTTTAAAAGCTCAGTCAGCAGGGACAGTCAGAAGTGTCTCCACCCC
AAGAGAAGAAAAGAAAGGTGGAGAAGGATACAATGAGTGTCAAGCACTCGAGGCTCTGTGCGCTTCACT
GGGCACCCGGCAAGCAGAACCTGAGCTCGACCTCCGCTCAATTAAGGAAGTCGATGAGGCAAAAGCTAAA
GAAGAAAAGTAGAGAAGTGTGGTGAGGATGATGAAACAATCCCATCTGAGTACAGATTAACCAGCCA
CGGATAAAGATGGAAAACCACTATTGCCAGAGCCTGAAGAAAACCAAGCCTCGGAGTGAATCAGAACT
CATTGATGAACTTTCAGAAGATTTGACCGGTCTGAATGTAAGAGAAAACCATCTAAGCCAACTGAAAAG
ACAGAAGAATCTAAGGCCGCTGCTCCAGCTCCTGTGTGCGAGGCTGTGTGTCGGACCTCCATGTGTAGTA
TACAGTCAGCACCCCTGAGCCGCTACCTTGAAGGGCACAGTGCCAGATGATGCTGTAGAAGCCTTGGC
TGATAGCCTGGGAAAAAGGAAGCAGATCCAGAAGATGGAAAACCTGTGATGGATAAAGTCAAGGAGAAG
GCCAAAGAAGAAGACCGTAAAAGCTTGGTAAAAAGAAGAAAACAATTCCTCCTGATTATAGATTAGAAG
AGGTCAAGGATAAAGATGGAAAGCCACTCTGCCAAAAGAGTCTAAGGAACAGCTTCCACCCATGAGTGA
AGACTTCTTCTGGATGCTTTGTCTGAGGACTTCTCTGGTCCACAAAATGCTTCTCTTAAATTTGAA
GATGCTAACTTGCTGCTGCCATCTCTGAAGTGGTTTCCAAACCCAGCTTCAACGACCCAAAGCTGGAG
CCCCACCCGTGATACCTCGAGTGACAAAGACCTCGATGATGCCTTGGATAAACTCTCTGACAGTCTAGG
ACAAAGGCAGCCTGACCCAGATGAGAACAACCAATGGAAGATAAAGTAAAGGAAAAAGCTAAAGCTGAA
CATAGAGACAAGCTTGGAGAAAGAGATGACACTATCCACCTGAATACAGACATCTCCTGGATGATAATG
GACAGGACAAACCAAGTGAAGCCACTACAAGAAATCAGAGGATTCAAAGAAACCTGCAGATGACCAAGA
CCCCATTGATGCTCTCTCAGGAGATCTGGACAGCTGTCCCTCCACTACAGAAACCTCACAGAACACAGCA
AAGGATAAGTGAAGAAGGCTGCTTCCAGCTCCAAGCACCTAAGAATGGAGGTAAGGCGAAGGATTCAG
CAAAGACAACAGAGGAACTTCCAAGCCAAAAGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001284212

ORF Size: 1998 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001284212.3](#)

RefSeq Size: 4195 bp

RefSeq ORF: 2001 bp

Locus ID: 831

UniProt ID: [P20810](#)

Cytogenetics: 5q15

MW: 72.2 kDa

Gene Summary: The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]