

Product datasheet for **RC239272**

WDR48 (NM_001303403) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR48 (NM_001303403) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	WDR48
Synonyms:	Bun62; P80; SPG60; UAF1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239272 representing NM_001303403
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCCCATCACCGCAGAACACAGCAGGGCGGAGGAAAGTGCAGGTTTCCTATGTTATTTCGAGATG
AAGTGGAGAAGTACAACCGAAATGGAGTCAATGCTCTGCAGCTGGATCCAGCACTAAATAGACTTTTCAC
AGCCGGTCGAGACTCTATCATAAGAATATGGAGTGTCAATCAGCACAAGCAAGATCCATATATAGCATCT
ATGGAACACCATACTGATTGGGTAACGACATTGACTCTGTTGTAATGGGAAAACATTAAGATGGA
ATGCACACAAGGGATTTGCATGTCAACATTAAGGACACATAAGGATTACGTAAGGCCTTAGCATATGC
CAAGGATAAAGAAGTACTAGCATCAGCTGGGTTGGACAGACAAATATTCCTTTGGGATGTGAATACTCTA
ACAGCATTGACTGCCTCAAATAACTGTCACTTCTCTTTAAGTGGAAACAAGATTCATTATA
GCCTGGCCATGAATCACTGGGAACAATCATTGTATCAGGGTCCACTGAAAAGGTGTACGGGTATGGGA
TCCAAGAACATGTGCAAACTAATGAAGCTTAAGGGCACACGGATAATGTGAAGGCATTGCTATTAAC
AGAGATGGCACGCAATGCCTGTACAGGCAGTTCTGATGGGACAATTCGCCTTTGGTCCCTTGCCAGCAGA
GATGTATAGCAACATACCGAGTCCATGATGAAGGTGTTGGGCGCTGCAAGTCAATGATGCCTTCACACA
TGTGATTCTGGTGAAGGGACAGGAAGATTTATTGTACAGACCTAAGAAACCCTGACATTCGGGTGCTA
ATTTGTGAAGAAAAGCACCAGTTCTCAAGATGGAGCTTGATAGATCAGCTGATCCTCCTCTGCAATTT
GGGTTGCAACAATAAGTCTACAGTAAATAAATGGACTTTGAAAGGAATTCATAATTTAGAGCCTCTGG
AGATTATGACAATGACTGTACAAATCCTATAACACCTCTTTGTACACAACCTGACCAGGTTATTAAGGG
GGTGTAGTATTATTCAGTGCCACATTTCTAATGATAAGAGACATATTAACCAAGATACCAATAATA
ATGTGGCATATTGGGATGATTGAAGGCATGTAAGTTGAAGATCTGGGCAAAGTGGATTTGAAGATGA
AATTAAGAAAAGATTTAAAATGGTGTATGTGCCAAATTTGGTTCTCAGTAGACTTAAAAACAGGGATGTTA
ACTATTACTTTGGATGAAAGTGATTGTTTTGCTGCCTGGGTTTCTGCAAAAAGATGCTGGTTTCAGCAGCC
CTGATGGGTGAGATCCAAAATTGAATTTAGGAGGACTTTTACTCCAAGCACTCCTGGAATATTGGCCTAG
AACACATGTGAATCCAATGGATGAAGAGGAAAATGAAGTAAACCATGTAATGGGGAGCAGGAGAACCGA
GTGCAGAAGGGAAATGGATATTTCAAGTGCCCCACATACACCCGTGATCTTTGGTGAAGCTGGAGGTC
GCACACTGTTCAAGCTGCTCTGCCGAGATTCGGGGGTGAGACTGAGTCTATGCTTCTTAATGAAACAGT
GCCACAATGGGTAAATGACATCACTGTGGATAAAAAATATGCCCAAATTCACAAAATTCCTTTCTACCTC
CAACCTCATGCATCTTCAGGAGCAAAAACCTTAAAAAAGATAGACTCTCTGCTAGTGACATGCTCCAAG
TCCGAAAAGTTATGGAACATGTTTATGAAAAAATTATCAACTGGATAATGAGTCTCAAACCACTAGCTC
TTCTAATAATGAAAAACCAGGAGAACAGGAAAAAGAAAGATATTGCTGTGTTGGCAGAGGAGAAAATT
GAACCTTTGTGCCAGGACCAGGTTTTGGATCCAAATATGGACCTTCGAACAGTGAACACTTCATATGGA
AGAGCGGTGGAGACCTCACCTCCATTACCGTCAGAAGTCCACG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239272 representing NM_001303403
 Red=Cloning site Green=Tags(s)


MAAHRQNTAGRRKVVQSYVIRDEVEKYNRNGVNALQLDPALNRLFTAGRDSIIRIWSVNHKQDPYIAS
 MEHHTDWVNDIVLCCNGKTLKVVNAHKGFCEMSTLRTHKDYVKALAYAKDELVASAGLDRQIFLWDVNTL
 TAL TASNNTVTTSSLGNKDSIYSLAMNQLGTII VSGSTEKVLRVWDPRTCAKLMKLGHTDNVKALLN
 RDGTQCLSGSSDGTIRLWSL GQQRCIATYRVHDEGVWALQVNDAFTHVYSGGRDRKIYCTDLRNPDIRVL
 ICEEKAPVLKMELDRSADPPPAIWWATTSTVNKWTLKG IHNFRASGDYDNDCTNPITPLCTQPDQVIK
 GASIIQCHILNDRHILTKDTNNNVAYWVWLKACKVEDLGKVD FEDEIKKRFKMYYVVPNWFVSDLKGTGML
 TITLDESDCFAA WSAKDAGFSSPDGSDPKNLGGLLLQALLEYWRTHVNPMEDEEENEVNHVNGEQENR
 VQKNGYFQVPPHTPVIFGEAGRTLFRLLCRDSGGETESMLLNETVPQWVIDITVDKNMPKFNKIPFYL
 QPHASSGAKTLKKDRLSASDMLQVRKVMHVYEKIINLDNESQTTSSSNNEKPGEQEKEEDIAVLAEEKI
 ELLCQDQVLDPNMDLRTVKHFIWKS GGD LTLHYRQKST

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001303403

ORF Size: 2004 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001303403.2](#)

RefSeq Size: 4000 bp

RefSeq ORF: 2007 bp

Locus ID: 57599

UniProt ID: [Q8TAF3](#)

Cytogenetics: 3p22.2

MW: 75.8 kDa

Gene Summary: The protein encoded by this gene has been shown to interact with ubiquitin specific peptidase 1 (USP1), activating the deubiquitinating activity of USP1 and allowing it to remove the ubiquitin moiety from monoubiquitinated FANCD2. FANCD2 is ubiquitinated in response to DNA damage. [provided by RefSeq, Sep 2016]