

Product datasheet for **RC239277**

MCAK (KIF2C) (NM_001297656) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MCAK (KIF2C) (NM_001297656) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KIF2C
Synonyms:	CT139; KNSL6; MCAK
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239277 representing NM_001297656
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGATTGATTTTGTATGATGTGGCTGCAATAAACCCAGAAGCTTTACAGCTTCTCCCTTACATCCGAAGG
ACAATCTGCCCTTGCAGGAAAATGTAACAATCCAGAAAACAAAACGGAGATCCGTCAACTCCAAAATTCC
TGCTCCAAAAGAAAGTCTTGAAGCCGCTCCACTCGCATGTCCACTGTCTCAGAGCTTCGCATCACGGCT
CAGGAGAATGACATGGAGGTGGAGCTGCCTGCAGCTGCAAACTCCCGCAAGCAGTTTTTCAGTTCTCTCTG
CCCCACTAGGCCTTCTGCCCTGCAGTGGCTGAAATACCATTGAGGATGGTCAGCGAGGAGATGGAAGA
GCAAGTCCATTCCATCCGAGGCAGCTCTTCTGCAAACCCTGTGAACTCAGTTCGGAGGAAATCATGTCTT
GTGAAGGAAGTGGAAAAATGAAGAACAAGCGAGAAGAGAAGAAGGCCCAAGCTCTGAAATGAGAATGA
AGAGAGCTCAGGAGTATGACAGTAGTTTTCCAAGTGGGAATTTGCCCGAATGATTAAGAATTTCCGGC
TACTTTGGAATGTCATCCACTTACTATGACTGATCCTATCGAAGAGCACAGAATATGTGTCTGTGTAGG
AAACGCCCACTGAATAAGCAAGAATTGGCCAAGAAAGAAATTGATGTGATTTCCATTCTAGCAAGTGTC
TCCTCTTGGTACATGAACCCAAGTTGAAAGTGGACTTAACAAAGTATCTGGAGAACCAAGCATTCTGCTT
TGACTTTGCATTTGATGAAACAGCTTCGAATGAAGTTGTCTACAGGTTACAGCAAGGCCACTGGTACAG
ACAATCTTTGAAGGTGGAAAAGCAACTTGTTCATATGGCCAGACAGGAAGTGGCAAGACACATACTA
TGGGCGGAGACCTCTCTGGGAAAGCCAGAATGCATCCAAAGGGATCTATGCCATGGCCTCCCGGACGT
CTTCTCTGAAGAATCAACCCTGCTACCGGAAGTTGGCCTGGAAGTCTATGTGACATTTCTCGAGATC
TACAATGGGAAGCTGTTTGACCTGCTCAACAAGAAGCCCAAGCTGCGCGTCTGGAGGACGGCAAGCAAC
AGGTGCAAGTGGTGGGCTGCAGGAGCATCTGGTTAACTCTGCTGATGATGTCATCAAGATGATCGACAT
GGGCAGCGCCTGCAGAACCTCTGGGCAGACATTTGCCAACTCCAATTCTCCCGCTCCCACGCGTGCTTC
CAAATTATTCTCGAGCTAAAGGGAGAATGCATGGCAAGTTCTCTTTGGTAGATCTGGCAGGGAATGAGC
GAGGCGCGGACACTTCCAGTGTGACCGGCAGACCCGCATGGAGGGCGCAGAAATCAACAAGAGTCTCTT
AGCCCTGAAGGAGTGCATCAGGGCCCTGGGACAGAACAAGGCTCACACCCGTTCCGTGAGAGCAAGCTG
ACACAGGTGCTGAGGGACTCCTTATTGGGGAAGTCTAGGACTTGCATGATTGCCACGATCTCACCAG
GCATAAGCTCCTGTAATATACTTTAACACCCTGAGATATGCAGACAGGGTCAAGGAGCTGAGCCCCA
CAGTGGGCCAGTGGAGAGCAGTTGATTCAAATGGAACAGAAGAGATGGAAGCCTGCTCTAACGGGGCG
CTGATTCCAGGAATTTATCCAAGGAAGAGGAGGAAGTCTTCCAGATGTCCAGCTTTAACGAAGCCA
TGACTCAGATCAGGGAGCTGGAGGAGAAGGCTATGGAAGAGCTCAAGGAGATCATACAGCAAGGACCAGA
CTGGCTTGAGCTCTCTGAGATGACCGAGCAGCCAGACTATGACCTGGAGACCTTTGTGAACAAAGCGGAA
TCTGCTCTGGCCAGCAAGCCAAGCATTCTCAGCCCTGCGAGATGTCATCAAGGCCTTGCCTGGCCA
TGCAGCTGGAAGAGCAGGCTAGCAGACAATAAGCAGCAAGAAACGGCCCCAG

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239277 representing NM_001297656
 Red=Cloning site Green=Tags(s)


MIDFDDVAAINPELLQLLPLHPKDNLPLQENVTIQKQKRRSVNSKIPAPKESLSRSTRMSTVSELRITA
 QENDMEVELPAAANSRKQFSVPPAPTRPSCPVAEIPLRMVSEEMEEQVHSIRGSSSANPVNSVRRKSCS
 VKEVEKMKNKREEKKAQNSEMRMKRAQEYDSSFNWEFARMIKEFRATLECHPLTMTDPIEEHRICVCVR
 KRPLNKQELAKKEIDVISIPSKCLLLVHEPKLVDLTKYLENQAFCFDFAFDETASNEVVYRFTARPLVQ
 TIFEGGKATCFAYGQTSGKTHTMGGDL SGKAQNASKGIYAMASRDVFLLNQPCYRKLGLEVYVTFE
 YNGKLFDLLNKKAKLRVLEDGKQVQVYVGLQEHLVNSADDVIKIMDMGSACRTSGQTFANSNSSRSHAC
 QIILRAKGRMHGKFSVLDLAGNERGADTSSADRQTRMEGAEINKSLLALKECIRALGQNKHAHTPFRESK
 LQVLRDSFIGENSRTCMIATISPGISSCEYTLNTRLRYADRVKELSPHSGPSGEQLIQMETEEMEACNGA
 LIPGNLSKEEEEELSSQMSFNEAMTQIRELEEKAMEELKEIIQQGPDWLELSEMTEQPDYDLETFVNKAE
 SALAQQAKHFSALRDVIKALRLAMQLEEQASRQISSKKRPQ

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297656

ORF Size: 2013 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297656.2](#)

RefSeq Size: 3014 bp

RefSeq ORF: 2016 bp

Locus ID: 11004

UniProt ID: [Q99661](#)

Cytogenetics: 1p34.1

Protein Families: Druggable Genome

MW: 76 kDa

Gene Summary: This gene encodes a kinesin-like protein that functions as a microtubule-dependent molecular motor. The encoded protein can depolymerize microtubules at the plus end, thereby promoting mitotic chromosome segregation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]