

Product datasheet for **RC239279**

CWH43 (NM_001286791) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CWH43 (NM_001286791) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CWH43
Synonyms:	CWH43-C; PGAP2IP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239279 representing NM_001286791
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGATCTATTACTTTCTTTGCAAACTAGAACTCACTGGGCTTGAAGTTTTAGTATAGCATTTCCTT
 CTCCAATATTCTAACAATTACTCTTTCTGGAAATTGGTTAACAAGAAGTGGATGCTAACCTGCTGAG
 GATAATCACTATTGGCAGCATAGCCTCCTTCCAGGCTCCAAATGCCAACTTCGACTGATGGTTCTTGCG
 CTTGGGGTGCTTCTCACTGATAGTGCAAGCTGTGACTTGGTGGTCAGGAAGTCATTTGCAAAGGTACC
 TCAGAATTTGGGATTCACTTTAGGACAGATTGTTCTTGTGTTCTACGCATATGGTATACTTCACTAAA
 CCCAATCTGGAGTTATCAGATGTCCAACAAAGTGATACTGACATTAAGTGCCATAGCCACACTTGATCGT
 ATTGGCACAGATGGTACTGCAGTAAACCTGAAGAAAAGAAGACTGGTGAGGTAGCCACGGGGATGGCCT
 CTAGACCCAACTGGCTGCTGGCAGGGGCTGCTTTTGGTAGCCTTGTGTTCTCACCCTGGGTTTTTGG
 AGAAGTCTCTCTGTTTCCAGATGGGCAGTGAGTGGGCATCCACATCCAGGGCCAGATCCTAACCCATTT
 GGAGGTGCAGTACTGCTGTGCTTGGCAAGTGGATTGATGCTTCCATCTTGTGTTGGTTTTCTGGTACTG
 GTTTGATCTGGTGGGTACAGGAACAGCTTCAGCTGCGGGGCTCCTTTACCTGCACACATGGGCAGCTGC
 TGTGCTGGCTGTGCTTTCGCCATCTTTACTGCATCCATGTGGCCCCAAACACTTGGACACCTTATTAAC
 TCAGGGACAAACCTGGGAAAACCATGACCATTGCCATGATATTTATCTTCTAGAAATATTTTTCTGTG
 CCTGGTGACAGCTTTTAAAGTTTGTCCCAGGAGGTGTCTACGCTAGAGAAAGATCAGATGTGCTTTTGGG
 GACAATGATGTTAATTATCGGGCTGAATATGCTATTTGGTCTAAGAAAAACCTTGACTTGCTTCTTCAA
 ACAAAAAACAGTTCTAAAGTGCTTTTCAGAAAGAGTGAAAAATACATGAACTTTTTCTGTGGCTGCTTG
 TTGGTGTGGGATTGTTGGGATTAGGACTACGGCATAAAGCCTATGAGAGAAAACCTGGGCAAAGTGGCACC
 AACCAAGAGGTCTCTGCTGCCATCTGGCCTTTACAGTTTGGATATGACAATGAAGGGTGGTCTAGTCTA
 GAAAGATCAGCTCACCTGCTCAATGAAACAGGTGCAGATTTTATAACAATTTTGGAGAGTGATGCTTCTA
 AGCCCTATATGGGGAACAATGACTTAACCATGTGGCTAGGGGAAAAGTTGGGTTTCTATACAGACTTTGG
 TCCAAGCACAAAGGTATCACACTTGGGGATTATGGCTTTGTCAAGATACCCAATTGTGAAATCTGAGCAT
 CACCTTCTCCGTCACCAGAGGGCGAGATCGCACCAGCCATCACATTGACCGTTAACATTTTCGGGCAAGC
 TGGTGGATTTTGTGCTGACACACTTTGGGAACCACGAAGATGACCTCGACAGGAACTGCAGGCTATTGC
 TGTTTCAAACTACTGAAAAGTAGCTCTAATCAAGTGATATTTCTGGGATATCACTTCAGCACCTGGC
 TCCAGAGATTATCTACAGCTCACTGAACATGGCAATGTGAAGGATATCGACAGCACTGATCATGACAGAT
 GGTGTGAATACATTATGTATCGAGGGCTGATCAGGTTGGGTTATGCAAGAATCTCCCATGCTGAAGTGA
 TGATTCAGAAATTCAGATGGCAAAATTTAGGATCCCTGATGACCCCACTAATTATAGAGACAACCAGAAA
 GTGGTCATAGACCACAGAGAAGTTTCTGAGAAAATTCATTTTAATCCCAGATTTGGATCCTACAAAGAAG
 GACACAATTATGAAAACAACCATCATTTTCATATGAATACTCCCAATACTTTTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239279 representing NM_001286791
 Red=Cloning site Green=Tags(s)


MIYYFPLQTELTGLEGFSIAFLSPIFLTITPFWKLVNKKWMLTLLRIITIGSIAFQAPNAKLRLMVLALGVSSSLIVQAVTWWSGSHLQRYLRIWGFILGQIVLVVLRWYTSNPIWSYQMSNKVILTL SAIATLDRIGTDGDCSKPEEKKTGEVATGMSRPNWLLAGAAFVSLVFLTHWVFGVSLVSRWAVSGHPHPGDPNPFGGAVLLCLASGLMLPSCWFRGTGLIWWVTGTASAAGLLYLHTWAAVSGCVFAIFTASMWPQTLGHLINSGTNPGKMTMIAMIFYLLEIFFCAWCTAFKFVPGGVYARERSDVLGTMMLIIGLNMLFGPKKNL D L L L QTKNSSKVLFRKSEKYMKLFLLVGVGLLGLRLHAYERKLGKVAPTKEVSAAIWPFRFGYDNEGWSLERSAHLNETGADFITILESDASKPYMGNDLTMWLGEKLGFTDFGPSTRYHTWGIMALSRYPIVKSEHLLPSPEGEIAPAITLTVNISGKLVDFVTHFGNHEDDLDRKLQAIKAVSKLLKSSSNQVIFLGYITSAPGSRDYLQLTEHGNVKDIDSTDHWRCEYIMYRGLIRLGYARISHAELSDSEIQMAKFRIPDDPTNYRDNQKVVIDHREVSEKIHFNPRFGSYKEGHNYENNHHFHMNTPKYFL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001286791

ORF Size: 2016 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001286791.1](#), [NP_001273720.1](#)

RefSeq Size: 2405 bp

RefSeq ORF: 2019 bp

Locus ID: 80157

UniProt ID: [Q9H720](#)


Cytogenetics: 4p11

Protein Families: Transmembrane

MW: 75.9 kDa

Gene Summary: Involved in lipid remodeling during GPI-anchor maturation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239279