

Product datasheet for **RC239289**

Estrogen induced gene 121 protein (KIAA1324) (NM_001284353) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Estrogen induced gene 121 protein (KIAA1324) (NM_001284353) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ELAPOR1
Synonyms:	EIG121; KIAA1324
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239289 representing NM_001284353
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTTCATTTCATCTCTGAGACACCAGACACAACACTCATGTACAAATGGGCCAAGCCGAAAATCTGTAGCG
AGGACCTTGAGGGGGCAGTGAAGCTGCCTGCTCTGGTGTGAAGACCCACTGCCACCCTGCAACCCAGG
CTTCTTCAAACCAACAACAGCACCTGCCAGCCCTGCCCATATGGTTCCTACTCCAATGGCTCAGACTGT
ACCCGCTGCCCTGCAGGGACTGAACCTGCTGTGGGATTTGAATACAAATGGTGAACACGCTGCCACAA
ACATGGAACGACCGTTCTCAGTGGGATCAACTTCGAGTACAAGGGCATGACAGGCTGGGAGGTGGCTGG
TGATCACATTTACACAGCTGCTGGAGCCTCAGACAATGACTTCATGATTCTCACTCTGGTTGTGCCAGGA
TTTAGACCTCCGAGTCGGTGTGGCAGACACAGAGAATAAAGAGGTGCCAGAATCACATTTGTCTTTG
AGACCCTCTGTTCTGTAACTGTGAGCTCTACTTCATGGTGGGTGAATTCTAGGACCAACACTCCTGT
GGAGACGTGGAAGGTTCCAAGGCAAACAGTCTATACCTACATCATTGAGGAGAACACTACCACGAGC
TTCACCTGGGCTTCCAGAGGACCACTTTTCATGAGGCAAGCAGGAAGTACACCAATGACGTTGCCAAGA
TCTACTCCATCAATGTCACCAATGTTATGAATGGTGTGGCCTCCTACTGCCGTCCCTGTGCCCTAGAAGC
CTCTGATGTGGGCTCCTCCTGCACCTCTTGTCTGCTGTTACTATATTGACCGAGATTCAGGAACCTGC
CACTCCTGCCCACTAACACAATTTGAAAGCCCACCAGCCTTATGGTGTCCAGGCTGTGTGCCCTGTG
GTCCAGGGACCAAGAACAACAAGATCCACTCTCTGTGCTACAACGATTGCACCTTCTCACGCAACACTCC
GACCAGGACTTTCACTACAACCTCTCCGCTTTGGCAAACACTGTCACTTTGCTGGAGGGCCAAGCTTC
ACTTCAAAGGGCTGAAATACTTCCATCACTTTACCCTCAGTCTCTGTGGAACCCAGGGTAGGAAAATGT
CTGTGTGCACCGACAATGTCACCTGACCTCCGGATTCCTGAGGGTGAGTCAGGGTTCTCCAAATCTATCAC
AGCCTACGTCTGCCAGGCAGTCATCATCCCCCAGAGGTGACAGGCTACAAGGCCGGGTTTCTCACAG
CCTGTCAGCCTTCTGATCGACTTATTGGGGTGACAACAGATATGACTCTGGATGGAATCACCTCCCCAG
CTGAACCTTTCCACCTGGAGTCTTGGGAATACCGGACGTGATCTTCTTTATAGGTCCAATGATGTGAC
CCAGTCTGCAGTTCTGGGAGATCAACCACCATCCGCGTCAGGTGCAGTCCACAGAAAATGTCCCTGGA
AGTTTGTGCTGCCAGGAACGTGCTCGGATGGGACCTGTGATGGCTGCAACTTCCACTTCTGTGGGAGA
GCGCGGCTGCTTGCCGCTCTGCTCAGTGGCTGACTACCATGCTATCGTCAGCAGCTGTGTGGCTGGGAT
CCAGAAGACTACTTACGTGTGGCGAGAACCAAGCTATGCTCTGGTGGCATTCTCTGCCTGAGCAGAGA
GTCACCATCTGCAAACCATAGATTTCTGGCTGAAAGTGGGCATCTCTGCAGGCACCTGTACTGCCATCC
TGCTCACCGTCTTGACCTGCTACTTTTGGAAAAAGAATCAAAAATAGAGTACAAGTACTCAAGCTGGT
GATGAATGCTACTCTCAAGGACTGTGACCTGCCAGCAGCTGACAGCTGCGCCATCATGGAAGGCGAGGAT
GTAGAGGACGACCTCATCTTTACCAGCAAGAAGTCACTCTTTGGGAAGATCAAATCATTTACCTCCAAGA
GGACTCCTGATGGATTTGACTCAGTGCCGCTGAAGACATCCTCAGGAGGCTAGACATGGACCTG

ACGCGTACGCGGCCGCTCGAGCAGAAAACACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239289 representing NM_001284353
 Red=Cloning site Green=Tags(s)

MFHSSLRHQQLMYKWKPKICSEDLGAVKLPASGVKTHCPPCNPFFKTNNSTCQPCPYGSYSNGSDC
 TRCPAGTEPAVGFEYKWWNTLPTNMETTTLSGINFEYKGMTGWEVAGDHIYTAAGASDNDFMILTLVVP
 FRPPQSVMAADTENKEVARITTFVETLCSVNCEL YFMVGVNSRTNTPVETWKGSKGKQSYTYIEENTTS
 FTWAFQRTTFHEASRKYTNDVAKIYSINVTNMVNGVASYCRPCALEASDVGSSTSCPAGYYIDRDSGTC
 HSCPTNTILKAHQPYGVQACVPCPGPTKNNKIHSLCYNDCTFSRNTPTRTFNYNFSALANTVTLAGGPSF
 TSKGLKYFHHFTLSL CGNQGRKMSVCTDNVTDLRIPEGESGFSKISITAYVCQAVIIPPEVTGYKAGVSSQ
 PVSADRLIGVTTDMTLDGITSPAELFHLES LGIPDVIFFYRNDVTQSCSSGRSTTIRVRCSPQKTVP
 SLLLPGTCDGTCDCGNFHLWESAAACPLCSVADYHAI VSSCVAGIQKTTYVWREPKLCSGGISLPEQR
 VTICKTIDFWLKVGISAGTCTAILLTVLTCYFWKKNQKLEYKSKLVMNATLKDCDLPAADSCAIMEGED
 VEDDLIFTSKSLFGKIKSFTSKRTPDGFDSVPLKTS SGGLDMDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001284353

ORF Size: 2025 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001284353.1](#), [NP_001271282.1](#)

RefSeq Size: 6914 bp

RefSeq ORF: 2028 bp

Locus ID: 57535

Cytogenetics: 1p13.3

Protein Families: Transmembrane

MW: 74.2 kDa

Gene Summary: Expression of this gene is induced by estrogen and the encoded protein has been characterized as a transmembrane protein. The encoded protein has been found in to correlate with survival in certain carcinomas (PMID: 21102415) and may be important for cellular response to stress (PMID: 21072319). Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2012]