

## Product datasheet for **RC239298**

### ITGB6 (NM\_001282355) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ITGB6 (NM_001282355) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ITGB6
Synonyms:	A11H
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC239298 representing NM\_001282355  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGGGATTGAACTGCTTTGCCCTGTTCTTTCTATTCTAGGAAGGAATGATCACGTACAAGTGGCTGTG  
CCCTGGGAGGTGCAGAAACCTGTGAAGACTGCCTGCTTATTGGACCTCAGTGTGCCTGGTGTCTCAGGA  
GAATTTTACTCATCCATCTGGAGTTGGCGAAAAGGTGTGATACCCAGCAAACCTTTTAGCTAAAGGATGT  
CAATTAACCTTCAACGAAAACCTGTCTCCCAAGTAGAAATACTTAAAAATAAGCCTCTCAGTGTAGGCA  
GACAGAAAAATAGTTCTGACATTGTTTCTGATTGCGCCTCAAAGCTTATCCTTAAGTTGAGACCAGGTGG  
TGCGCAGACTCTGCAGGTGCATGTCCGCCAGACTGAGGACTACCCGGTGGATTTGTATTACCTCATGGAC  
CTCTCCGCTCCATGGATGACGACCTCAACACAATAAAGGAGCTGGGCTCCCGCTTTCAAAGAGATGT  
CTAAATTAACCAGCAACTTTAGACTGGGCTTCGGATCTTTTGTGAAAAACCTGTATCCCCTTTTGTGAA  
AACACACCAGAAGAAATGGCAACCTTGCAGTAGTATCCATACTTCTGTTTACCTACATTTGGATTC  
AAGCACATTTTGGCATTGACAAAATGATGCTGAAAGATTCAATGAAATTGTGAAGAATCAGAAAAATTTCTG  
CTAATATTGACACACCCGAAGGTGGATTTGATGCAATTATGCAAGCTGCTGTGTGAAGGAAAAATTTGG  
CTGGCGGAATGACTCCCTCCACCTCCTGGTCTTTGTGAGTGATGCTGATTCTCATTTTGGAAATGGACAGC  
AACTAGCAGGCATCGTCATTCTAATGACGGGCTCTGTCACTTGGACAGCAAGAATGAATACTCCATGT  
CAACTGTCTTGGAAATATCCAACAATTGGACAACCTATTGATAAACTGGTACAAAACAACGTGTTATTGAT  
CTTCGCTGAACCAAGAACAAGTTCATTTATATGAGAATTACGCAAACTTATTCCTGGAGCTACAGTA  
GGTCTACTTCAGAAGGACTCCGGAACATTCTCCAGCTGATCATCTCAGCTTATGAAGAACTGCGGTCTG  
AGGTGGAACGGAAAGTATTAGGAGACACTGAAGGACTCAACTTGTCAATTTACAGCCATCTGTAACAACGG  
TACCCTCTTCCAACACCAAAAAGAAATGCTCTCACATGAAAGTGGGAGACACAGCTTCTTTCAGCGTGACT  
GTGAATATCCCACTGCGAGAGAAGAAGCAGGCACATTATCATAAAGCCTGTGGGGCTGGGGGATGCC  
TGGAACTACTTGTGAGCCAGAATGCAACTGCGACTGTCAGAAAGAAGTGGAAAGTGAACAGCTCCAATG  
TCACCACGGGAACGGCTCTTTCCAGTGTGGGTGTGTGCCTGCCACCCTGGCCACATGGGGCTCGCTGT  
GAGTGTGGCGAGGACATGCTGAGCACAGATTCTGCAAGGAGGCCAGATCATCCCTCTGCAGCGGAA  
GGGGTACTGCTACTGTGGCAGTGTATCTGCCACTTGTCTCCCTATGGAACATTTATGGGCTTATTG  
CCAGTGTGACAATTTCTCTGCGTGAGACAAAAGGGCTGCTCTGCGGAGATTTCTCAAAGGATGGTTCT  
GTTTCTGCTCTCTGCAAGGAGAAAATGAATGTCTTATTACATTCCTAATAACTACAGATAATGAGGGGA  
AAACCATCATTACAGCATCAATGAAAAAGATTGTCCGAAGCCTCAAACATTCCCATGATCATGTTAGG  
GGTTTCCCTGGCTATTCTTCTCATCGGGTTGTCTACTGTGCATCTGGAAGCTACTGGTGTCAATTCAT  
GATCGTAAAGAAGTTGCCAAATTTGAAGCAGAACGATCAAAGCCAAGTGGCAAACGGGAACCAATCCAC  
TCTACAGAGGATCCACAAGTACTTTAAAAATGTAACCTATAAACACAGGGAAAAACAAAAGGTAGACCT  
TTCCACAGATTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC239298 representing NM\_001282355  
 Red=Cloning site Green=Tags(s)


MGIELLCLFFLFLGRNDHVQGGCALGGAETCEDCLLIGPQCAWCAQENFTHPSGVGERCDTPANLLAKGC  
 QLNFIENPVSQVEILKNKPLSVGRQKNSSDIVQIAPQSLILKLRPQGAQTLQVHVRQTEDYPVDLYLMD  
 LSASMDDDLNTIKELGSRLSKEMSKLTSNFRLLGFGSFEKPVSPFVKTTPEEIANPCSSIPYFCLPTFGF  
 KHILPLTNDARFNEIVKNQKISANIDTPEGGFDAIMQAAVCKEKGWRNDSLHLLVFDADSHFGMS  
 KLAGIVIPNDGLCHLDSKNEYSMTVLEYPTIGQLIDKLVQNNVLLIFAVTQEQVHLYENYAKLIPGATV  
 GLLQKDSGNILQLIISAYEELRSEVELEVLGDTEGLNLSFTAICNNGTLFQHQKCKSHMKVGDTSFVS  
 VNIPHCERRSRHIIKPVGLDALELLVSPENCDCQKEVEVNSSKCHHGNGSFQCGVCACHPGHMGPRC  
 ECGEDMLSTDSCKEAPDHPSCSGRGDCYCGQCICHLSPYGNIIYGPYCQDNFSCVRHKGLLGGDFSKDGS  
 VCSLQGENECLITFLITDNEGKTIHISINEKDCPKPPNIPMIMLVSLAILLIGVLLCIWKLVSFH  
 DRKEVAKFEAERSKAKWQTGNPLYRGSTSTFKNVTYKHREKQKVDLSTDC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**


SgfI-MluI

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

## Plasmid Map:



ACCN: NM\_001282355

ORF Size: 2043 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001282355.2](#)

RefSeq Size: 4372 bp

RefSeq ORF:	2046 bp
Locus ID:	3694
UniProt ID:	<a href="#">P18564</a>
Cytogenetics:	2q24.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Regulation of actin cytoskeleton
MW:	75.4 kDa
Gene Summary:	This gene encodes a protein that is a member of the integrin superfamily. Members of this family are adhesion receptors that function in signaling from the extracellular matrix to the cell. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. The encoded protein forms a dimer with an alpha v chain and this heterodimer can bind to ligands like fibronectin and transforming growth factor beta 1. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]