

## Product datasheet for **RC239300**

### **MCAK (KIF2C) (NM\_001297655) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	MCAK (KIF2C) (NM_001297655) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KIF2C
Synonyms:	CT139; KNSL6; MCAK
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC239300 representing NM\_001297655  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCCATGGACTCGTCGCTTCAGGCCCGCTGTTCCCGGTCTCGCTATCAAGATCCAACGCGTAATG  
GTTTAATTCACAGTGCCAATGTAAGGACTGTGAACTTGGAGAAATCCTGTGTTTCAGTGGAAATGGGCAGA  
AGGAGGTGCCACAAAGGGCAAAGAGATTGATTTTGATGATGTGGCTGCAATAAACCCAGAACTCTTACAG  
CTTCTTCCCTTACATCCGAAGGACAATCTGCCCTTGCAGGAAAATGTAACAATCCAGAAACAAAAACGGA  
GATCCGTCAACTCCAAAATTCCTGCTCCAAAAGAAAGTCTTGAAGCCGCTCCACTCGCATGTCCACTGT  
CTCAGAGCTTCGCATCACGGCTCAGGAGAATGACATGGAGGTGGAGCTGCCTGCAGCTGCAAACCTCCCGC  
AAGCAGTTTTTCAGTTCCTCTTCGGAGGAAATCATGTCTTGTGAAGGAAGTGAAAAAATGAAGAACAAGC  
GAGAAGAGAAGAAGGCCAGAACTCTGAAATGAGAATGAAGAGAGCTCAGGAGTATGACAGTAGTTTTCC  
AAACTGGGAATTTGCCGAATGATTAAGAATTCGGGCTACTTTGGAATGTCATCCACTTACTATGACT  
GATCCTATCGAAGAGCACAGAATATGTGTCTGTGTTAGGAAACGCCACTGAATAAGCAAGAATTGGCCA  
AGAAAGAAATTGATGTGATTTCCATTCCTAGCAAGTGTCTCCTCTTGGTACATGAACCCAAGTTGAAAGT  
GGACTTAAACAAGTATCTGGAGAACCAAGCATTCTGCTTTGACTTTGCATTTGATGAAACAGCTTCGAAT  
GAAGTTGTCTACAGGTTACAGCAAGGCCACTGGTACAGACAATCTTTGAAGGTGAAAAAGCAACTTGTT  
TTGCATATGGCCAGACAGGAAGTGGCAAGACACATACTATGGGCGGAGACCTCTCTGGGAAAGCCAGAA  
TGCATCCAAAGGGATCTATGCCATGGCTCCCGGGACGCTTCTCCTCTGAAGAATCAACCCTGTACCGG  
AAGTTGGGCTGGAAGTCTATGTGACATTCTTCGAGATCTACAATGGGAAGCTGTTTGACCTGCTCAACA  
AGAAGGCCAAGCTGCGCGTGTGGAGGACGGCAAGCAACAGGTGCAAGTGGTGGGGCTGCAGGAGCATCT  
GGTTAACTCTGCTGATGATGTCATCAAGATGATCGACATGGGCAGCGCCTGCAGAACCTCTGGGCAGACA  
TTTGCCAACTCCAATTCCTCCCGCTCCACGCGTGTCTTCCAAATTATTCTTCGAGCTAAAGGGAGAATGC  
ATGGCAAGTTCTCTTTGGTAGATCTGGCAGGGAATGAGCGAGGCGCGGACACTTCCAGTGTGACCGGCA  
GACCCGCATGGAGGGCGCAGAAATCAACAAGAGTCTCTTAGCCCTGAAGGAGTGCATCAGGGCCCTGGGA  
CAGAACAAGGCTCACACCCGTTCCGTGAGAGCAAGCTGACACAGGTGCTGAGGGACTCCTTCATTGGGG  
AGAAGTCTAGGACTTGCATGATTGCCACGATCTCACCAGGCATAAGCTCCTGTGAATATACTTTAACAC  
CCTGAGATATGCAGACAGGTCAAGGAGCTGAGCCCCACAGTGGGCCAGTGGAGAGCAGTTGATTCAA  
ATGGAACAGAAGAGATGGAAGCCTGCTAACGGGGCGCTGATTCCAGGCAATTTATCCAAGGAAGAGG  
AGGAACTGTCTTCCAGATGTCCAGCTTAAACGAAGCCATGACTCAGATCAGGGAGCTGGAGGAGAAGGC  
TATGGAAGAGCTCAAGGAGATCATAAGCAAGGACCAGACTGGCTTGGAGCTCTCTGAGATGACCGAGCAG  
CCAGACTATGACCTGGAGACCTTTGTGAACAAAGCGGAATCTGCTCTGGCCAGCAAGCCAAGCATTTCT  
CAGCCCTGCGAGATGTATCAAGGCCTTGCCTGGCCATGCAGCTGGAAGAGCAGGCTAGCAGACAAAT  
AAGCAGCAAGAAACGGCCCCAG

AG**CGGACCC**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239300 representing NM\_001297655  
 Red=Cloning site Green=Tags(s)


MAMDSSLQARLFPGLAIKIQRSNGLIHSANVRTVNLEKSCVSEVAEGGATKGKEIDFDDVAAINPELLQ  
 LLPLHPKDNPLQENVTIQKQRRSVNSKIPAPKESLSRSTRMSTVSELRITAQENDMEVELPAAANSR  
 KQFSVPLRRKSLVKEVEKMKNKREEKKAQNSEMRMKRAQEYDSSFNWEFARMIKEFRATLECHPLTMT  
 DPIEEHRICVCVRKRPLNKQELAKKEIDVVISIPSKLLLVHEPKLVLDLTKYLENQAFCDFAFDEASN  
 EVVYRFTARPLVQTI FEGGKATCFAYGQTGSGKTHTMGGDLSGKAQNASKGIYAMASRDVFLLNQPCYR  
 KLGLEYVYVTFEYNGKLFDLLNKKAKLRVLEDGKQQVQVVGLEHLVNSADDVIKIDMGSA CRTSGQT  
 FANSNSRSHACFQIILRAKGRMHGKFSVLDLAGNERGADTSSADRQTRMEGAEINKSLALKECIRALG  
 QNKAHTPFRESKLTQVLRDSFIGENSRTCMIATISPGISSCEYTLNLTLYADRVKELSPHSGPSGEQLIQ  
 METEEMEACSNALIPGNLSKEEEEELSSQMSSFNEAMTQIRELEEKAMEELKEIQQGPDWLELSEMTEQ  
 PDYDLTFVNKAESALAQAKHFSALRDVIKALRLAMQLEEQASRQISSKKRPQ

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



## Plasmid Map:



ACCN: NM\_001297655

ORF Size: 2052 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001297655.2](#)

RefSeq Size: 2773 bp

RefSeq ORF: 2055 bp

Locus ID: 11004

UniProt ID: [Q99661](#)

Cytogenetics: 1p34.1

Protein Families: Druggable Genome

MW: 77.4 kDa

**Gene Summary:** This gene encodes a kinesin-like protein that functions as a microtubule-dependent molecular motor. The encoded protein can depolymerize microtubules at the plus end, thereby promoting mitotic chromosome segregation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]