

## Product datasheet for **RC239303**

### **IQCE (NM\_001287500) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	IQCE (NM_001287500) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	IQCE
Synonyms:	1700028P05Rik; PAPA7
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC239303 representing NM\_001287500  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGTTCTGGGCACCGGGGAGCCGGCCTTGGACACGAAAGCAAAAAGAAAGCTTCCACAAACCTCCAC  
CCACATCGCCAAAGTCACCTTATCTCTCTAAGCCGAGAAAAGTGGCCTCCTGGAGGTCCCTCAGGACGGC  
AGGGAGCATGCCTCTGGGCGGCCGAGCGTCCCTGACCCCGCAGAAGCTGTGGCTGGGAACCGCAAAGCCA  
GGAAGTCTGACCCAGGCCTGAACTCACCCCTCACCTGGGAGCATGCGTGGACTGGCGTCCCGGCGGCA  
CTCCTGACTGTCTGACAGACACCTTACAGTGAAGAGGCCACATCTCAGGCGCTCTGCCAGCAACGGTCA  
TGTCCCTGGGACTCCTGTCTACAGAGAAAAGAAGATATGTATGACGAGATTATTGAGTTAAGAAGTCA  
TTGCACGTGCAGAAGAGCGACGTGGACCTGATGAGAACGAAGCTCCGGCGCCTGGAGGAGAAAACAGCA  
GGAAGGACCGGCAGATAGAGCAGCTCCTGGATCCAGCCGCGGCACGGATTTTGTTCGGACTCTGGCAGA  
GAAAAGGCCCGATGCCAGTTGGGTCAATTAACGGGCTGAAGCAGAGGATCCTGAAGCTGGAACAGCAGTGC  
AAGGAGAAGGACGGCACCATCAGCAAATCCAGACCGATATGAAGACTACCAACCTGGAAGAGATGCGGA  
TCGCCATGGAGACATACTACGAGGAGGTGCATCGTCTCCAGACCTCTTGGAAGTTCTGAAACCCCGG  
AAGAAGCCCCTGGGGGAGAAGAAGACGGGCGCCAAAAGGCAGAAGAAGATGGGCAGTGCCTCCTGAGC  
TTGTCCCGAGGTGTCAGGAGCTCACGGAAGAGAACCAGAGCCTGAAGGAGGACCTGGACCCGCTGCTGA  
GCACCTCCCAACCATCTCCAAGACACAGGGTTATGTGGAGTGGAGCAAGCCCGGCTGTGAGGCGCAT  
TGTGGAGCTGGAGAAGAACTAAGTGTGATGGAGAGCTCAAAATCACACGCCGACAGCCAGTGCAGATCA  
CACCCGCCAGCCTGCCTTGATCCAGCTCTGCGCTGCACAGACAGCCACGAGGGGACCGCAACAAGGACC  
ACGAGCGTCTCCGAGGGGCTGTGAGAGACCTGAAGGAAGAGCGGACCCGCGCTGCAGGAGCAGTCTGTCA  
GAGAGATTTGGAGGTGAAGCAGCTCCTGCAGGCGAAGGCCGACCTGGAGAAGGAGCTGGAGTGCAGGAGG  
GAGGGCGAGGAGGAGGAGAGAGAGCGAGAGGAGTTTTGAGAGAGGAGATTACAGACACTACCAGCAAGC  
TCCAAGAATTGCAAGAAATGAAGAAAAGAAGAAAAGAGGATTGCCCGGAAGTTCTCATAAGGCCCAAGA  
GCTCCAGCTCCCACTCCAGCAGCAGGCACTGCGAGCAAGACTGGCCGCGGATTCCAGCGAGGAGGGG  
CTCCCGCGGCCCGCTCCCTGCTCTGATGGGAGAAGAGACGCCGCGGCCAGAGTCTGCAGGCCAGT  
GGAAGGTGTACAAGCACAAGAAAAAAGGCTGTTCTGGATGAGGCGGCTGTGGTCTTCAGGCAGCTTT  
CAGGGGACATCTACGCGGACAAAGCTTTAGCAAGCAAAGCACATGGCTCAGAGCCACCCAGCGTGCCA  
GGCCTCCAGACCAGAGCTCTCCTGTGCCCGGTTCCGAGCCCATCGCCAGGCCACGGGACAGCCCTG  
TGCAGGAGGAGGCCATCGTCATCATCCAGTCCGCTCTGCGGGCACACCTGGCCCGGGCCAGGCACAGTGC  
TACCGGTAAAAGAACCACCACCGCAGCTTCTACCAGGAGGAGATCGGCTTCAGCCACACAGGGGACGCC  
TCTCCCAACCTTCTCGCAGCTTCTCCTGGTAAATTTCAATTTGGATTCTGGCACCAGGGCAGCTG  
GTCTGCTGAGGTCTCAGCCAAAGAGTGGCCACCTCCAGGAAGCCCGGGCTGTGTGCGGACGGAAGGGAGG  
AGTGTCCCATCTGGAGCTGCCTCTGAAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC239303 representing NM\_001287500  
 Red=Cloning site Green=Tags(s)


MFLGTGEPALDTKAKRKAFHKPPPTSPKSPYLSKPRKVASWRSRLRTAGSMPLGGRASLTPQKLWLTAKP  
 GSLTQALNSPLTWEHAWTGVPGGTPDCLTDTFRVKRPHLRRSASNGHVPGTPVYREKEDMYDEIIELKKS  
 LHVQKSDVDLMRTKLRLEEENSRKDRQIEQLLDPSRGTDFVRTLAEKRPDASWVINGLKQRILKLEQQC  
 KEKDGTI SKLQTD MKTTNLEEMRIAMETYYEEVHRLQTL LASSETTGKKPLGEKKTGAKRQKKMGSALLS  
 LSRSVQELTEENQSLKEDLDRVLSTSPITISKTQGYVEWSKPRLLRRIVELEKKLSVMESSKSHAAEPVRS  
 HPPACLASSALHRQPRGDRNKDHERL R GAVRDLKEERTALQEQLLQRDLEVKQLLQAKADLEKELECAR  
 EGEEERREREVLR E E I Q T L T S K L Q E L Q E M K K E E K E D C P E V P H K A Q E L P A P T P S S R H C E Q D W P P D S S E E G  
 L P R P R S P C S D G R R D A A R V L Q A Q W K V Y K H K K K A V L D E A A V V L Q A A F R G H L T R T K L L A S K A H G S E P P S V P  
 G L P D Q S S P V R P V S P I A Q A T G S P V Q E E A I V I I Q S A L R A H L A R A R H S A T G K R T T T A A S T R R R S A S A T H G D A  
 S S P P F L A A L P G N F I H L D S G T R A A G L L R S Q P K S G H L Q E A R A V S G R K G G V S H L E L P L N


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001287500

**ORF Size:** 2058 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001287500.2](#)

**RefSeq Size:** 2466 bp

RefSeq ORF: 2061 bp

Locus ID: 23288

Cytogenetics: 7p22.3

MW: 77 kDa

**Gene Summary:** Component of the EvC complex that positively regulates ciliary Hedgehog (Hh) signaling (By similarity). Required for proper limb morphogenesis (PubMed:28488682).[UniProtKB/Swiss-Prot Function]