

Product datasheet for **RC239313**

CDT2 (DTL) (NM_001286230) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDT2 (DTL) (NM_001286230) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DTL
Synonyms:	CDT2; DCAF2; L2DTL; RAMP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239313 representing NM_001286230
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTCTTCAATTCGGTGTCCGCCAGCCAGCTTGGCGTCTGAGAAATGCTCCCAATATGGAACATG
TACTAGCAGTTGCCAATGAAGAAGGCTTTGTTCGATTGTATAACACAGAATCACAAAGTTTCAGAAAGAA
GTGCTTCAAAGAATGGATGGCTCACTGGAATGCCGTCTTTGACCTGGCCTGGGTTCTGGTGAACCTTAA
CTTGTTACAGCAGCAGGTGATCAAACAGCCAAATTTGGGACGTAAGAGCTGGTGAAGCTGATTGGAACAT
GCAAAGGTCATCAATGCAGCCTCAAGTCAGTTGCCTTTCTAAGTTTGAGAAAGCTGATTCTGTACGGG
TGGAAGAGATGGCAACATTATGGTCTGGGATACCAGGTGCAACAAAAAGATGGGTTTTATAGGCAAGTG
AATCAAATCAGTGGAGCTACAATACCTCAGACAAGCAAACCCCTTCAAACCCAAGAAGAACAGAAATT
CAAAGGACTTGCTCCTCTGTGGATTTCCAGCAAAGTGTACTGTGGTCTCTTCAAGACGAGAATAC
CTTAGTCTCAGCAGGAGCTGTGGATGGGATAATCAAAGTATGGGATTTACGTAAGAATTACTGCTTAT
CGACAAGAACCCATAGCATCCAAGTCTTTCCTGTACCCAGGTAGCAGCACTCGAAAACTTGATATTCAA
GTCTGATTTTGGATTCCACTGGCTCTACTTTATTTGCTAATTGCACAGACGATAACATCTACATGTTTAA
TATGACTGGGTTGAAGACTTCTCCAGTGGCTATTTTCAATGGACACCAGAAGCTACCTTTTATGTA
TCCAGCCTTAGTCCAGATGACCAGTTTTTATGTCAGTGGCTCAAGTATGAAGCTGCCTACATATGGAAGG
TCTCCACACCTGGCAACCTCCTACTGTCTCTGGGTCAATTCTCAAGAGGTCACGTCTGTGTGCTGGT
TCCATCTGACTTCACAAAGATTGCTACCTGTTCTGATGACAATACTAAAACTGGCGCTTGAATAGA
GGCTTAGAGGAGAAACCAGGAGGTGATAAACTTCCACGGTGGTGGGCTCTCAGAAGAAAAAGAGT
CAAGACTGGCTAGTAACAGTAACGAGTAGCCAGAGTACTCTGCCAAAGCCCCAGGGTAAAGTGCAA
TCCATCCAATTCTTCCCGTCATCCGCAGCTTGTGCCCAAGCTGTGCTGGAGACCTCCCTCTTCTTCA
AATACTCTACGTTCTCTATTAACCTCTCTGCAAGGCCGGTCTCCATCAACAGAAGAGGCTCTG
TCTCTCCGTCTCTCCAAGCCACCTTCTCTTCAAGATGTCGATTAGAAAAGGGTGACCCGAACACC
TTCCTCATCACCACCCATCACTCCACCTGCTTCGGAGACCAAGATCATGTCTCCGAGAAAAGCCCTTATT
CCTGTGAGCCAGAAGTCATCCAAGCAGAGGCTTGCTCTGAGTCTAGAAATAGAGTAAAGAGGAGGCTAG
ACTCAAGCTGTCTGGAGAGTGTGAAACAAAAGTGTGTGAAGAGTTGTAAGTGTGACTGAGCTTGTATGG
CCAAGTTGAAAACTTTCATTTGGATCTGTGCTGCCTTGTGTTAACCAGGAAGACCTTAGTAAGGACTCT
CTAGGTCCTACCAATCAAGCAAAATGAAGGAGCTGGTACCAGTATCTCAGAGCCTCCGCTCCTATCA
GTCCGTATGCTTCAAGAAAGCTGTGGAACGCTACCTCTTCTTTGAGACCTTGTGGAGAAGGGTCTGAAAT
GGTAGGCAAAGAGAATAGTTCCCCAGAGAATAAAAAGTGGTGTGGCCATGGCAGCCAAACGGAAGGCT
GAGAATCCATCTCCACGAAGTCCGTATCCCAGACACCAATTCCAGGAGACAGAGCGGAAAGACATTGC
CAAGCCCGGTACCATCACGCCAGCTCCATGAGGAAAACTGCACATACTTCCATAGAAAGTCCAGGA
GGACTTCTGTGGTCTGAACACTCAACAGAATTA


ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239313 representing NM_001286230
 Red=Cloning site Green=Tags(s)

MLFN^{SVLRQPQLGVL}RNAPNMEHVLAVANEEGFVRL YNTESQSF^{RKKCFKEWMAH}WNAVFDLAWVPGELK
 LVTAAGDQTAKFWDVKAGELIGTCKGHQCSLKSVAFSKFEKAVFTGGRDGNIMVWDTRCNKKDGFYRQV
 NQISGAHNTSDKQTPSKPKKKQNSKGLAPSVDFQQSVTVVLFQDENTLVSAGAVDGIK^{VWDLRKN}YNTAY
 RQEP^{IASKSF}LYPGSS^{TRKL}GYSSLILDSTGSTL^{FANCTDDNIYMFNMTGLK}TSPVAIFNGHQNSTFYVK
 SSLSPDDQFLVSGSSDEAAYIWKVSTPWQPPTVLLGHSQEVTSVCWCPSDFTKIATCSDDN^{LTKIWRLNR}
 GLEEKPGGDKLSTVGWASQKKESR^{PGLVTVTSSQSTPAKAPRVKCNPSN}SSPSSAACAPSCAGDLPLPS
 NTPTFSIKTSPAKARSPINRRGSVSSVSPKPPSSFKMSIRNWVTRTPSSSP^{ITPPASET}KIMSPRKALI
 PVSQKSSQAEACSESRNRV^{KRRLDSSCLESVKQKCVKSCNCVTELDGQVENLHDL}CCLAGNQEDLSKDS
 LGPTKSSKIEGAGT^{ISEPPSPISPYASESCGTLPLR}PCGEGSEMVGKENS^{SPENKNWLLAMA}AKRKA
 ENSPSRPSSQTPNSRRQSGKTLPSV^{ITPSSMRKICTYFHRKSQEDFCGPEH}STEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:
Cloning Scheme:



Plasmid Map:


ACCN: NM_001286230

ORF Size: 2064 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286230.1](#), [NP_001273159.1](#)

RefSeq Size: 4485 bp

RefSeq ORF: 2067 bp

Locus ID: 51514

Cytogenetics: 1q32.3

Protein Families: Druggable Genome

MW: 75.5 kDa

Gene Summary: Substrate-specific adapter of a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex required for cell cycle control, DNA damage response and translesion DNA synthesis. The DCX(DTL) complex, also named CRL4(CDT2) complex, mediates the polyubiquitination and subsequent degradation of CDT1, CDKN1A/p21(CIP1), FBH1, KMT5A and SDE2 (PubMed:16861906, PubMed:16949367, PubMed:16964240, PubMed:17085480, PubMed:18703516, PubMed:18794347, PubMed:18794348, PubMed:19332548, PubMed:20129063, PubMed:23478441, PubMed:23478445, PubMed:23677613, PubMed:27906959). CDT1 degradation in response to DNA damage is necessary to ensure proper cell cycle regulation of DNA replication (PubMed:16861906, PubMed:16949367, PubMed:17085480). CDKN1A/p21(CIP1) degradation during S phase or following UV irradiation is essential to control replication licensing (PubMed:18794348, PubMed:19332548). KMT5A degradation is also important for a proper regulation of mechanisms such as TGF-beta signaling, cell cycle progression, DNA repair and cell migration (PubMed:23478445). Most substrates require their interaction with PCNA for their polyubiquitination: substrates interact with PCNA via their PIP-box, and those containing the 'K+4' motif in the PIP box, recruit the DCX(DTL) complex, leading to their degradation. In undamaged proliferating cells, the DCX(DTL) complex also promotes the 'Lys-164' monoubiquitination of PCNA, thereby being involved in PCNA-dependent translesion DNA synthesis (PubMed:20129063, PubMed:23478441, PubMed:23478445, PubMed:23677613). The DDB1-CUL4A-DTL E3 ligase complex regulates the circadian clock function by mediating the ubiquitination and degradation of CRY1 (PubMed:26431207).[UniProtKB/Swiss-Prot Function]