

Product datasheet for **RC239325**

ITGB6 (NM_001282354) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ITGB6 (NM_001282354) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ITGB6
Synonyms:	A11H
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239325 representing NM_001282354
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGGATTGAACTGCTTTGCCCTGTTCTTTCTATTCTAGGAAGGAATGATCACGTACAAGTGGTGCGC
AGACTCTGCAGGTGCATGTCCGCCAGACTGAGGACTACCCGGTGGATTTGTATTACCTCATGGACCTCTC
CGCTCCATGGATGACGACCTCAACACAATAAAGGAGCTGGGCTCCCGGCTTTCCAAAGAGATGTCTAAA
TTAACCAGCAACTTTAGACTGGGCTTCGGATCTTTTGTGGAAAAACCTGTATCCCCTTTTGTGAAAAACA
CACCAGAAGAAATGCCAACCTTGCACTAGTATTCCATACTTCTGTTTACCTACATTTGGATTCAAGCA
CATTTTGGCATTGACAAATGATGCTGAAAGATTCAATGAAATTGTGAAGAATCAGAAAAATTTCTGCTAAT
ATTGACACACCCGAAGGTGGATTTGATGCAATTATGCAAGCTGCTGTGTGAAGGAAAAATTTGGCTGGC
GGAATGACTCCCTCCACCTCCTGGTCTTTGTGAGTGATGCTGATTCTCATTTTGGATGGACAGCAAACT
AGCAGGCATCGTATTCCCTAATGACGGGCTCTGCACTTGGACAGCAAGAATGAATACTCCATGTCAACT
GTCTTGGAAATATCCAACAATTGGACAACCTATTGATAAACTGGTACAAAACAACGTGTTATTGATCTTCG
CTGTAACCCAAGAACAAGTTCATTTATATGAGAATTACGCAAACTTATTCTGGAGCTACAGTAGGTCT
ACTTCAGAAGGACTCCGGAAACATTCTCCAGCTGATCATCTCAGCTTATGAAGAACTGCGGTCTGAGGTG
GAACTGGAAGTATTAGGAGACACTGAAGGACTCAACTTGTCAATTTACAGCCATCTGTAACAACGGTACCC
TCTTCCAACACAAAAGAAATGCTCTCACATGAAAGTGGGAGACACAGCTTCTTCAGCGTGACTGTGAA
TATCCCACTGCGAGAGAAGAAGCAGGCACATTATCATAAAGCCTGTGGGGCTGGGGATGCCCTGGAA
TTACTTGTGAGCCAGAATGCAACTGCGACTGTGCAAGAAAGTGAAGTGAACAGCTCCAAATGTCACC
ACGGGAACGGCTCTTCCAGTGTGGGTGTGTGCCTGCCACCCTGGCCACATGGGGCTCGCTGTGAGTG
TGGCGAGGACATGCTGAGCACAGATTCTGCAAGGAGGCCCCAGATCATCCCTCCTGCAGCGGAAGGGGT
GACTGCTACTGTGGCAGTGTATCTGCCACTTGTCTCCCTATGGAAACATTTATGGGCCTTATTGCCAGT
GTGACAATTTCTCCTGCGTGAGACACAAAAGGGCTGCTCTGCGGAGGTAACGGCGACTGTGACTGTGGTGA
ATGTGTGTGCAGGAGCGGCTGGACTGGCGAGTACTGCAACTGCACCACCAGCACGGACTCCTGCGTCTCT
GAAGATGGAGTGTCTGCAGCGGGCGGGGACTGTGTTTGTGGCAAGTGTGTTTGCACAAACCCTGGAG
CCTCAGGACCAACCTGTGAACGATGTCTACCTGTGGTGACCCCTGTAACCTAAACGGAGCTGCATTGA
GTGCCACCTGTCAGCAGCTGGCCAAGCCCAGAAAGATGTGTGGACAAGTGCAACTAGCTGGTGGCACC
ATCAGTGAAGAAGAAGATTTCTCAAAGGATGGTTCTGTTTCTGCTCTCTGCAAGGAGAAAAATGAATGTC
TTATTACATTCCTAATAACTACAGATAATGAGGGGAAAACCATCATTACAGCATCAATGAAAAAGATTG
TCCGAAGCCTCCAAACATTTCCATGATCATGTTAGGGGTTTCCCTGGCTATTCTTCTCATCGGGTTGTG
CTACTGTGCATCTGGAAGCTACTGGTGTCAATTCATGATCGTAAAGAAGTTGCCAAATTTGAAGCAGAAC
GATCAAAGCCAAGTGGCAAACGGGAACCAATCCACTCTACAGAGGATCCACAAGTACTTTTAAAAATGT
AACTTATAAACACAGGAAAAACAAAAGGTAGACCTTTCCACAGATTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239325 representing NM_001282354
 Red=Cloning site Green=Tags(s)

MGIELLCLFFLFLGRNDHVQGAQTLQVHVRQTEDYPVDLYYMLDL SASMDDLNTIKELGSRLSKEMSK
 LTSNFRFLGFGSFVEKPVSPFVKTTPEEIANPCSSIPYFCLPTFGFKHILPLTND AERFNEIVKNQKISAN
 IDTPEGGFDAIMQAAVCKEIKIGWRNDSLHLLVFVSDADSHFGMDSKLAGIVIPNDGLCHLDSKNEYSMST
 VLEYPTIGQLIDKLVQNNVLLIFAVTQEQQVHLYENYAKLIPGATVGLLQKDSGNILQLIISAYEELRSEV
 ELEVLGDTEGLNLSFTAICNNGTLFQHQQKCSHMKVGDTSF SVTVNIPH CERRSRHII IKPVGLGDALE
 LLVSPENCDCQKEVEVNSSKCHHGNGSFQCGVCAHPGHMGRPRCEGEDMLSTDSCKEAPDHPSCSRG
 DCYCGQCI CHLSPYGN IYGPYQC DNFCVVRHKGL LCGNGDCDCGECVCRSGWTGEYCNTTSTDCSVS
 EDGVLCSGRGDCVCGKCVCTNPGASGPTCERCPTCGDPCNSKRSCIECHLSAAGQAREECVDKCKLAGAT
 ISEEDFSKDGSVSCSLQGENECLITFLITTDNEGKTI IHSINEKDCPKPPNIPMIMLGVSLAILLIGV
 LLCIWKLLV SFHDRKEVAKFEAERSKAKWQTGTNPLYRGSTSTFKNVTYKHREKQKVDLSTDC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001282354

ORF Size: 2079 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282354.1](#), [NP_001269283.1](#)

RefSeq Size: 4408 bp

RefSeq ORF:	2082 bp
Locus ID:	3694
Cytogenetics:	2q24.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Regulation of actin cytoskeleton
MW:	76.3 kDa
Gene Summary:	<p>This gene encodes a protein that is a member of the integrin superfamily. Members of this family are adhesion receptors that function in signaling from the extracellular matrix to the cell. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. The encoded protein forms a dimer with an alpha v chain and this heterodimer can bind to ligands like fibronectin and transforming growth factor beta 1. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]</p>