

Product datasheet for **RC239330**

MS4A14 (NM_001261827) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MS4A14 (NM_001261827) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MS4A14
Synonyms:	MS4A16; NYD-SP21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239330 representing NM_001261827
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGTCAACATCCCAGGACAGAAGGGCAACTCAGTCATCACTATAAAACCAAACGAACTGTATTGA
 CTGCAATTTCCCTACAGACCTCATAGCTCTCTGCTGGATTTTCTGAAGGGAGAGCCAAGAGTCTTGGGGGC
 TACCCAGATCCTGCTTGCTCTAATCATTGTGGGCTTTGGAATATATTTGCACCTAATTACATCGGTTTC
 TCCCAAAGACTTCCCCTTGTTGCTCCTCACAGGATATCCATTCTGGGGAGCACTTATTGGTCAAGGTGTCA
 CGGGCATGAATGTTATCAGCTCCTTGGTTGCGATAACTGGGATTACTTTCACCATCTCAGCTACAGACA
 TCAAGACAAGTACTGCCAGATGCCATCCTTTGAAGAAATATGTGTTTTCACTAGAACTCTTTTCATTGGA
 ATTTTGTAACTTACTGATCATCAGCATAGCAGAGCTCAGCATCTCTGTGACTATTGCATCCTTTAGAA
 GCAAGTGTGGACACAGTCAGATGAGGTTCTGTTTTCTTGCCTTCGGATGTTACTCAAAATAGTGAACA
 ACCTGCCCCAGAAGAAATGATCAATTACAATTTGTGCTTCAAGAAGAGTTTTCCAGTGATGATTCAACA
 ACAAAATGCACAACTCTGTTATCTTTGGAGGCTATGCTTTCTTCAAGTTAACTCTCTAGGAGTCCCTTAG
 TCTCCCAACCAAGTAATAAAGGTAGAGAATTTGTGCCAGATGAACAAAAGCAAAGTATCCTTCCATCTCC
 CAAATTTTCAGAGGAAGAAATTTGAACCTTTGCCTCCCACTAGAGAAAAAGCCCTCAGAAAAATATGTCC
 ATTCAGCTAGACTCTACATTTAAACAAATGAAAGATGAAGATCTACAATCTGCTATTGTACAACCTTCTC
 AAATGCAAACCAAGCTTCTGCAGGACCAAGCTGCGTCACTCCAAGTTTTTCCATCCCATTCTGCACTAAA
 ACTCGAAGATATATCACCTGAAGACTTGCCATCCCAAGCTCTACCAGTAGAAGGCCTGTGAGAACAACCC
 ATGCCATCTAAGTCTACATCATCCCATGTCAAACAGTCTTCTAATCTGACAGCTAATGACCTGCCCCCTC
 AAGGCATACTATCCCAAGACACATCATCTCAAGATATGCTGTTTCATGACATGACATCCCAAGATATGCA
 ATCCCTAGATATGCTATCTCAAGACACACCATCCACGCCATGCCACCTCAAGACATACCTTCCCAAGAT
 ATGCTATCCCAAGCTCTATCAGCGCATGCCATATTACCTGAAGCCTCAACATCCCATATTGTGCAGTTCC
 CTGAAATACAACACCTACTTCAGCAGCCCCAGATCTTCAACCAGAAAACACTGAACCTCAAAACCAGCA
 AATTTTACAAATGTCATATCAAGATATTAGATCAGAAGTTATGGAAGAGACCAAAGAATGGAAATCTGAG
 GAGGAATCCATAGAAGAAAATCCTCAAGACGGCATTCTTAAACCAGCAAACCAAGCCTTGCAATACT
 TAAGGAGACATTCTTTAGACGTGCAAGCCAAAGGCCAGAAATCCTCAAAGAGGCATTCTTAGATCAGCA
 AAGCAAAGGCTGGCAATCTCAAAGCAGAAATCCTTAGACCAGCAAATCAAAGACTGGCTATCCCCAAAG
 AGGCACTCCGTAGATAAGCAAGCTCAACTTAATCAAACCTAAAGAGCAACTCCAGATCAGCAAGCTGAAG
 ATCAGCAAGCCAAAGGGGAACAATACCCAGAAGGACAATCTAAAGATGGACAAGTTAAAGACCAGCAGAC
 TGATAAGGAGCAAACTCAAAGAAGCAAACCCAGGATCAGCAAACCTGAAGACCAGCCGGCCCAAGAGAAG
 AAATCCCCGAAAGGACAATTCCAAAATGTTCAAGCCGAAGGACAGCAAGCTCAGGTGGAGAAAGTGCCAA
 AACTGTTATGCCAAGATTCAGAAATCCCAATACAGCAATACCAATTCTGGCAATTCACAAAGGCAATCT
 CCAGGCTGGACAACCCAGGACTGTCAATCTTTGGCCAAGAAATCCCCTGACTGGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239330 representing NM_001261827
 Red=Cloning site Green=Tags(s)


MESTSQDRRATHVITIKPNETVLTAFPYRPHSSLLDFLKGEPRVLGATQILLALIIVGFGTIFALNYIGF
 SQRLPLVLTGYPFWGALIGQVGTGMNVISSLVAITGITFTILSYRHQDKYQMPSFEEICVFSTRFIG
 ILLILLIISIAELSSIVTIAFRSKCWTQSDVLFLLPSDVTQNSEQPAPEENDQLQFVLQEEFSSDDST
 TNAQSVIFGGYAFFKLTLSRSPVLSQPGNKGREFVPDEQKQSILSPKFSEEEIEPLPPTLEKKPSNM
 IQLDSTFKQMKDEDLQSAIVQPSQMQLKLLQDQAASLVFSPHSALKLEDISPDLPSQALPVEGLSEQT
 MPKSTSSHVQSSNLTANDLPPQGILSQDTSSQDMLFDMTSQDMQSLDMLSQDTPSHAMPQDIPSD
 MLSQALSAHAILPEASTSHIVQFPEIQHLLQPPDLQPEPTEPQNQILQMSYQDIRSEVMEETKEWKSE
 EELHRRKSSRRHSLNQQTKALQYLRRHSLDVQAKGQKSSKRHSLDQSKGWQSPKQSLDQKIDWLSPK
 RHSVDKQAQLNQTKQLPDQAEDQAKGEQYPEGQSKDGQVKDQQTQDKEQNSKKQTQDQQTEDQPAQEK
 KSPKGQFQNVQAEGQQAQVEKVPKLLCQDSESQIQYQFWQFHKGNLQAGQPRVTNLLAKNPLTG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001261827

ORF Size: 2085 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001261827.1, NP_001248756.1

RefSeq Size: 3401 bp

RefSeq ORF: 2088 bp

Locus ID: 84689


Cytogenetics: 11q12.2

Protein Families: Transmembrane

MW: 78.7 kDa

Gene Summary: May be involved in signal transduction as a component of a multimeric receptor complex. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239330