

Product datasheet for **RC239337**

ZNF434 (ZSCAN32) (NM_001284527) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF434 (ZSCAN32) (NM_001284527) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF434
Synonyms:	HCCS-5; ZNF434
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239337 representing NM_001284527
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGGCTGCAGTGAAGAGTACCGAGGCACACCCATCCTCAAACAAGGATCCACACAGGGCCAGAAAT
 CAGCCCTCCAGGGTAACAGCCCTGACTCCGAGGCCTCCCGTCAGCGCTTCAGGCAGTTTTGCTACCAGGA
 GGTAAGTGGCCACATGAAGCTTTAGCAAACCTGCGGAACCTGTTGTAGTGGCTGAGGCCGAAGACC
 CACTCAAAGAGGAAATCCTGGAGCTGCTGGTTTTGGAGCAGTTCTGACTATCTTGCCAGAGGAGATCC
 AGACCTGGGTGAGGGAGCAGCATCCAGAAAACGGCGAGGAAGCTGTGGCTCTGGTTGAGGATGTACAGAG
 AGCTCCTGGACAACAGGTTCTAGATTCTGAGAAGGACTTGAAAGTACTCATGAAGGAGATGGCCCTTTG
 GGAGCAACCAGAGAATCACTGAGATCCCAATGGAAACAGGAGGTTGAGCCAGAGGAACCGACTTTTAAGG
 GATCACAGAGCTCACACAAAGACCAGGGGAACAGTCAGAAGCCTGGCTTGCTCCTCAGGCTCCCAGGAA
 CCTGCCTCAAACACAGGTCTCCACGACCAGGAGACAGGTGCTGTGGTCTGGACAGCTGGGTCCCAGGGA
 CCAGCCATGCGTGACAACAGAGCTGTATCCCTCTGTGACGAAGAATGGATGTGCCAGGCCTGCACAAA
 GGGCCCTCTACAGGGGTGCCACCCAGAGGAAGGACAGTCACGTCTCGCTGGCAACAGGTGTGCCCTGGGG
 CTATGAAGAGACCAAGACGCTCCTGGCTATTCTTAGTAGTTCTCAATTTTATGGAAGTCCAGACCTGT
 CAGCAGAACAGCCAGATCTACAGGGCCATGGCGGAAGGACTCTGGGAGCAGGGTTTTCTGCGGACCCAG
 AACAGTGTGCGACCAAGTTCAAAGCCTACAGTTGAGTTACCGCAAAGTGAGGAGAGGCCGTGTGCTGA
 GCCTTGATCTTTTATGAGGAAATGAATGCTCTTTCAGGCTCCTGGGCCTCTGCACCTCCTATGGCAAGC
 GATGCTGTTCTTGCCAAAGAAGGAAGTGATATTGAGGCTGGAGAGCTGAATCACCAGAATGGGGAACCCA
 CGGAGGTAGAAGATGGCACTGTGGATGGTGCAGACAGGGATGAAAAGGACTTCAGGAATCCTGGCCAGGA
 AGTCAGGAAACTAGACCTGCCAGTGTGTTCCCAAACAGACTTGTTTTGAGTTCAAGAACGAGATTAAA
 AAAGAAAATCTAAATGGGATGATTAGAGGAAGTAGAAATAACAAGGCTTTACAGAGAAAGTCCAGAG
 GAGTTTATTGGCACTCTGAGCTACAAAAGGCTTGAGAGTGAGCCAACATCAAGAAGGCAATGTAGAAA
 TTCTCCAGGGGAGAGTGAGGAGAAAACCCATCCAGGAGAAGATGAGTCACCAGAGTTTTTGTGCCAGG
 GACAAAGCCTGTACACATATCCTCTGTGGGAAAACTGCTCTCAGAGTGTGCACTCTCCCCACAAGCCAG
 CGCTCAAACCTGAAAAAGTATCTCAATGTCCTGAATGTGGGAAAACCTTTAGCCGAAGTTCTTATCTTGT
 TCGGCATCAAAGAATCCACACAGGCGAGAAGCCTCACAAGTGCAGTGAGTGCGGGAAGGGCTTTAGTGAG
 CGCTCAAACCTCACTGCCACCTACGAACCTCACACAGGGGAGAGGCCCTATCAGTGTGGCAATGTGGGA
 AAAGCTTCAACCAGAGTTCCAGCCTCATTGTCCACCAGAGGACCCATACCGGGGAAAAGCCTTACCAGTG
 CATTGTCTGTGGAAGAGATTCAACAACAGTTCACAGTTCAAGTGTCAACGGCGCATCCACACTGGGGAG
 AGCCCATACAAGTGTGCAAGTGTGTGGGAAAATCTTCAACAATAGCTCCCCTTCAAGTGTGCCACCGAAAA
 CCCCACTGGTGAAAAGCCTTACAGGTGTTCTCACTGTGAGAGAGGCTTCACTAAGAACTTGCCCTCAC
 CCGTCATCAGACAGTACACATGAAAGCAGTACTCTCATCACAGGAAGGAAGAGATGCGTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239337 representing NM_001284527
 Red=Cloning site Green=Tags(s)

MMAAAVKSTEHPSSNKPDTQQQKSALQGNPDSEASRQRFQFCYQEVTPGHEAFSKLWELCCQWLRPKT
HSKEEILELLVLEQFLTILPEEIQTWVREQHPENGEEAVALVEDVQRAPGQQVLDSEKDLKVLKMEMAPL
GATRESLSQWKQEVQPEEPTFKGSQSSHQRPGEQSEAWLAPQAPRNLPQNTGLHDQETGAVVWTAGSQG
PAMRDNRAVSLCQEQWMCPPAQRALYRGATQRKDSHVSLATGVPWGYEETKTLLAILSSSQFYGKLQTC
QQNSQIYRAMAEGLEWQGLRTPEQCRTKFKSLQLSYRKVRRGRVPEPCIFYEEMNALSGSWASAPPMAS
DAVPQGEGSDIEAGELNHQNGEPTVEVDGTVDGADRDEKDFRNPQGEVRKLDLPVLPFNRLGFEFKNEIK
KENLKWDDSEEVEINKLQRKSRGVYWHSELQKGLESEPTSRRCQRNSPGESEKTPSQEKMSHSFCAR
DKACTHILCGKNCSQSVHSPHKPALKLEKVSQCPECGKTFSSSYLVRHQRIHTGEKPHKCECGKFSE
RSNLTAHLRTHTGERPYQCGQCGKSFNQSSSLIVHQRTHTGEKPYQCIVCGKFRNNSQFSAHRRHTGE
SPYKCAVCGKIFNNSSHFSAHRKTHTGEKPYRCSHCERGFTKNSALTRHOTVHMKAVLSSOEGRDAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001284527

ORF Size: 2091 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001284527.2](#)

RefSeq Size: 3178 bp

RefSeq ORF: 2094 bp

Locus ID: 54925

UniProt ID: [Q9NX65](#)


Cytogenetics: 16p13.3

Protein Families: Transcription Factors

MW: 79.2 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239337