

## Product datasheet for **RC239348**

### NFAT2 (NFATC1) (NM\_001278672) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NFAT2 (NFATC1) (NM_001278672) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NFATC1
Synonyms:	NF-ATC; NF-ATc1.2; NFAT2; NFATc
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC239348 representing NM\_001278672  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGACGGGGCTGGAGGACCAGGAGTTCGACTTCGAGTTCCTCTTCGAGTTTAACACGCGCAGCAGGGCG  
 CCGCCGCGGCCGCCAGAACACTATGGCTATGCATCCTCCAAGTTCAGCCCCGCCCTGCCGCTCCAC  
 GGCGACTCCACCCTGCCGCGCCCGTGCACAACTTCAGACCTCCACACCGGGCATCATCCCGCCGGC  
 GATCACCCCTCGGGTACGGAGCAGCTTTGGACGGTGGGCCCGCGGGCTACTTCTCTCTCCGCCACA  
 CCAGGCCTGATGGGGCCCTGCCCTGGAGATCCTCGCATCGAGATAACCTCGTCTGGCCTGTACCA  
 CAACAATAACAGTTTTTCCACGATGTGGAGGTGAAGACGTCTCCCTAGTCCAAACGGTCCCCTCC  
 ACGGCCACGTGAGTCTGCCAGCCTGGAGGCTACAGAGACCCCTCGTCCCTGAGCCCGGCCAGCAGCC  
 TGTCTCCCGGAGCTGCAACTCAGAGGCTCCTCTACGAGTCCAACACTCGTACCCGTACGCGTCCCC  
 CCAGAGCTCGCATGGCAGTCTCCCTGCGTGTCTCCAAGACCACGGACCCCGAGGAGGGCTTTCCCGC  
 GGGCTGGGGCCCTGCACACTGCTGGGTTCCCGCGGCACTCCCCTCCACCTCGCCCCGCCAGCGTCA  
 CTGAGGAGAGCTGGCTGGGTGCCGCTCCTCCAGACCCGCTCCCCTTGAACAAGAGGAAGTACAGCT  
 CAACGGCCGGCAGCCGCCCTACTCACCCACCACTCGCCACGCGCTCCCAGCAGGCTCCCCGCGGGT  
 AGCGTGACCGAGACTCGTGGTGGGCAACACCACCACTACACAGCTCGGCCATCGTGGCCGCATCA  
 ACGCGCTGACCACCGACAGCAGCTGGACCTGGGAGATGGCGTCCCTGTCAAGTCCCAGAACACCCT  
 GGAGCAGCCGCCCTCAGTGGCGCTCAAGGTGGAGCCCTCGGGGAGGACTGGGCAGCCCCCGCCCCG  
 GCCGACTTCGCGCCCGAAGACTACTCTCTTTCCAGCACATCAGGAAGGGCGGCTTCTGCGACAGTACC  
 TGGCGGTGCCGACACCCCTACCAGTGGCGAAGCCCAAGCCCTGTCCCCTACGTCTACATGAGCCC  
 GACCTGCCCGCCCTGGACTGGCAGCTGCCGTCCCACTCAGGCCGATGAGCTTCGGATTGAGGTGCAG  
 CCCAAGTCCCACCACCGAGCCACTACGAGACGGAGGCAAGCCGGGGGGCGTGAAGGCGTGGCCGGAG  
 GACACCCATCGTGCAGCTGCATGGCTACTTGGAGAATGAGCCGCTGATGCTGCAGCTTTTCATTGGGAC  
 GGCGGACGACCGCTGCTGCGCCGCACGCCTTCTACCAGGTGCACCGCATCACAGGGAAGACCGTGTCC  
 ACCACCAGCCACGAGGCCATCCTCTCAACACCAAAAGTCTGGAGATCCCACTCCTGCCGGAACAGCA  
 TGCGAGCCGTCATTGACTGTGCCGAATCCTGAAACTCAGAACTCCGACATTGAACTTCGAAAGGAGA  
 GACGGACATCGGAGGAAGAACACACGGGTACGGCTGGTGTCCGCGTTCAGTCCCAGAACCCAGCGGC  
 CGCAGCTGTCCCTGCAGGTGGCTCCAACCCATCGAATGCTCCAGCGCTCAGCTCAGGAGCTGCCTC  
 TGGTGGAGAAGCAGACAGCAGGACAGCTATCCGGTCTGGGCGGGAAGAAGATGGTCTGTCTGGCCACA  
 CTTCTGCAGGACTCCAAGTCAATTTTCGTGGAGAAAGCCCAAGATGGCCACCATGTCTGGGAGATGGAA  
 GCGAAAAGTACCGGGACCTGTGCAAGCCGAATTTCTGGTGGTTGAGATCCCAGCATTTTCGGAATCAGA  
 GGATAACCGCCCGTTCAGTCAAGTTTCTACGTCTGCAACGGGAAGAGAAAGCGAAGCCAGTACCAGCG  
 TTTACCTACCTTCCCGCAACGTAATGAAATAATACGAAATGACCTCTCCAGCACGAGCACCCACTCC

AG**CGGACCC**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC239348 representing NM\_001278672  
 Red=Cloning site Green=Tags(s)


MTGLEDQEFDFEFLFEFNQRDEGAAAAAPEHYGYASSNVSPALPLPTAHSTLPAPCHNLQTSTPGIIPPA  
 DHPSGYGAALDGGPAGYFLSSGHTRPDGALESPRIEITSCLGLYHNNNQFFHDVEVEDVLPSSKRSPS  
 TATLSLPSLEAYRDPSCLSPASSLSSRSCNSEASSYESNYSYPYASPTSPWQSPCVSPKTTDPEEGFPR  
 GLGACTLLGSPRHSPSTSPRASVTEESWLGARSSRPASPCNKRYSLNGRQPPYSPHHSPTSPHGS  
 SVTDDSWLGNNTQYTSSAIVAAINALTTDSSLDLGDGVPVKSRKTTLEQPPSVALKVEPVGEDLGSPPPP  
 ADFAPEDYSSFHIRKGGFCQYLAVPQHPYQWAKPKPLSPTSYSMPTLPALDWQLPSHSGPYELRIEVQ  
 PKSHHRAHYETEGSRGAVKASAGGHPVQLHGYLENEPLMLQLFIGTADDRLLRPHAFYQVHRITGKTVS  
 TTSHEAILSNTKVLEIPLLPENSMRAVIDCAGILKLRNSDIELRKGETDIGRKNTRVRLVFRVHVQPQSG  
 RTLSLQVASNPIECSQRSQELPLVEKQSTDSYPVVGKKMVLSGHNFLQDSKIVFEKAPDGHVWEME  
 AKTDRDLCKPNSLVVEIPFRNQRITSPVHVSFVVCNGKRRKSQYQRFYLPANVNEIIRNDLSSTSTHS

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:**


Sgfl-RsrII

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

## Plasmid Map:



ACCN: NM\_001278672

ORF Size: 2100 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001278672.2](#)

RefSeq Size: 3967 bp

RefSeq ORF:	2103 bp
Locus ID:	4772
UniProt ID:	<a href="#">O95644</a>
Cytogenetics:	18q23
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Axon guidance, B cell receptor signaling pathway, Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway
MW:	76.9 kDa
Gene Summary:	<p>The product of this gene is a component of the nuclear factor of activated T cells DNA-binding transcription complex. This complex consists of at least two components: a preexisting cytosolic component that translocates to the nucleus upon T cell receptor (TCR) stimulation, and an inducible nuclear component. Proteins belonging to this family of transcription factors play a central role in inducible gene transcription during immune response. The product of this gene is an inducible nuclear component. It functions as a major molecular target for the immunosuppressive drugs such as cyclosporin A. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. Different isoforms of this protein may regulate inducible expression of different cytokine genes. [provided by RefSeq, Jul 2013]</p>