

Product datasheet for **RC239350**

Elastin (ELN) (NM_001278914) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Elastin (ELN) (NM_001278914) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ELN
Synonyms:	ADCL1; SVAS; WBS; WS
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239350 representing NM_001278914
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCGGGTCTGACGGCGCGGCCCGCGGCCGGAGTCTCCTGCTCCTGCTGTCCATCCTCCACCCT
 CTGCGCCTGGAGGGTCCCTGGGGCCATTCTGGTGGAGTTCCTGGAGGAGTCTTTTATCCAGCGCTGGG
 GCCTGGAGGCAAACCTCTTAAGCCAGTTCCTGGAGGGCTTGGCGGTGCTGGCCTTGGGGCAGGGCTCGGC
 GCCTTCCCGCAGTTACCTTTCCGGGGCTCTGGTGCCTGGTGGAGTGGCTGACGCTGCTGCAGCTATA
 AAGCTGCTAAGGCTGGCCTGGGCTTGGTGGTGTCCAGGAGTGGTGGCTTAGGAGTGTCTGCAGCCCC
 TTCTGTGCCAGGTGCGGTGGTTCCTCAGCCTGGAGCCGGAGTGAAGCCTGGGAAAGTGCCGGGTGTGGGG
 CTGCCAGGTGTATACCCAGGTGGCGTGTCCAGGAGCTCGTTCCCGGTGTGGGGGTGTCCCTGGAG
 TTCCACTGGAGCAGGAGTTAAGCCAAGGCTCCAGGTGTAGGTGGAGCTTTTGTGGAATCCAGGAGT
 TGGACCCTTTGGGGGACCGAACCTGGAGTCCACTGGGGTATCCCATCAAGCCCCCAAGCTGCCTGGT
 GGTATGGACTGCCCTACACCACAGGAAACTGCCCTATGGCTATGGGCCCGGAGGAGTGGCTGGTGCAG
 CGGGCAAGGCTGGTTACCAACAGGGACAGGGGTTGGCCCCAGGCAGCAGCAGCAGCGGAGCTAAAGC
 AGCAGCAAAGTTCGGTGTGGAGCAGCCGGAGTCTCCCTGGTGTGGAGGGGCTGGTGTTCCTGGCGTG
 CCTGGGGCAATTCCTGGAATGGAGGATCGCAGGCGTGGGACTCCAGCTGCAGCTGCAGCTGCAGCAG
 CAGCCGCTAAGGCAGCCAAGTATGGAGCTGTGCAGGCTTAGTGCCTGGTGGGCCAGGCTTTGGCCCCGG
 AGTAGTTGGTGTCCAGGAGCTGGCGTCCAGGTGTGGTGTCCAGGAGCTGGGATTCCAGTTGTCCCA
 GGTGCTGGGATCCAGGTGTGCGGTCCAGGGTGTGTACCAGAAGCAGCTGCTAAGGCAGCTGCCAA
 AGGCAGCCAAATACGGGGCCAGGCCGGAGTCCGAGTTGGAGGCATTCTACTACGGGGTTGGAGCTGG
 GGGCTTTCCCGCTTTGGTGTGCGAGTCCGAGGTATCCCTGGAGTCCGAGGTGTCCCTGGTGTGCGAGGT
 GTTCCCGAGTCCGAGGTGTCCCGGAGTTGGCATTTCGCCCAAGCTCAGGCAGCAGCTGCCGCCAAGG
 CTGCCAAGTACGGGTTAGTTCCTGGTGTGCGCGTGGCTCCTGGAGTTGGCGTGGCTCCTGGTGTGCGGTG
 GGCTCCTGGAGTTGGCTTGGCTCCTGGAGTTGGCGTGGCTCCTGGAGTTGGTGTGGCTCCTGGCGTTGGC
 GTGGCTCCCGCATTGGCCCTGGTGGAGTTGCAGCTGCAGCAAAATCCGCTGCCAAGGTGGCTGCCAAAG
 CCCAGCTCCGAGCTGCAGCTGGGCTTGGTGTGGCATCCCTGGACTTGGAGTTGGTGTGCGCGTCCCTGG
 ACTTGGAGTTGGTGTGTTCTGGACTTGGAGTTGGTGTGGTGTTCCTGGCTTCGGGGCAGTACCT
 GGAGCCCTGGCTGCCGCTAAAGCAGCCAAATATGGAGCAGCAGTGCCTGGGGTCTTGGAGGGCTCGGGG
 CTCTCGGTGGAGTAGGCATCCAGGCGGTGGTGGGAGCCGGACCCGCCCGCTGCCGCAGCCAA
 AGCTGTGCCAAAGCCGCCAGTTTGGCCTAGTGGGAGCCGCTGGGCTCGGAGGACTCGGAGTCCGAGGG
 CTTGGAGTTCCAGGTGTTGGGGGCTTGGAGGTATACCTCCAGCTGCAGCCGCTAAAGCAGCTAAATACG
 GTGCTGTGGCCTTGGAGGTGTCTAGGGGGTCCCGGGCAGTTCCCACTTGGAGGAGTGGCAGCAAGACC
 TGGCTTCGGATTGTCTCCATTTCCAGGTGGGGCTGCCTGGGAAAGCTTGTGGCCGGAAGAGAAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239350 representing NM_001278914
 Red=Cloning site Green=Tags(s)

MAGLTAAAPRPGVLLLLLSILHPSRPGGVPGAIPGGVPGGVFYPALGPGGKPLKPVPGLAGAGLGAGLG
 AFPAVTFPGALVPGGVADAAAAYKAAKAGAGLGGVPGVGGGLVSAAPSVPGAVPQPAGAVKPKVPGVG
 LPGVYPGGVLPGARFPVGVLPVPTGAGVKPKAPGVGGAFAGIPGVGPFGGPQPGVPLGYPIKAPKLP
 GYGLPYTTGKLPYGYGPGGVAGAAGKAGYPTGTGVGPQAAAAAAAAAKAAKFGAGAAGVLPVGGAGVPGV
 PGAIPGIGGIAGVGTAAAAAAAAAAKAAKYGAAAGLVPGGPGFPGVVGVPAGVPGVGPAGIPVVP
 GAGIPGAAVPGVVSPEAAAKAAKAAKYGARPGVGGIPTYGVGAGGFPFVGVGIPVAGVPGVGG
 VPGVGGVPGVVISPEAQAAAAKAAKYGLVPGVGVAPGVGVAPGVGVAPGVGLAPGVGVAPGVGVAPGVG
 VAPGIGPGGVAAAAKSAKVAAKAQLRAAAGL GAGIPGLGVGVGPGLGVGAGVPLGVGAGVPGF GAVP
 GALAAAKAAKYGAAVPGVLGGLGALGGVGPVGGVAGPAAAAAAAAAKAAKAAQFGLVGAAGLGGLVGG
 LGVPGVGLGGIPPAAAKAAKYGAAGLGGVLGGAGQFPLGGVAARPGFGLSPIFPGGACL GKACGRKRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001278914

ORF Size: 2100 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278914.2](#)

RefSeq Size: 3717 bp

RefSeq ORF: 2103 bp

Locus ID: 2006

Cytogenetics: 7q11.23

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

MW: 61.2 kDa

Gene Summary: This gene encodes a protein that is one of the two components of elastic fibers. Elastic fibers comprise part of the extracellular matrix and confer elasticity to organs and tissues including the heart, skin, lungs, ligaments, and blood vessels. The encoded protein is rich in hydrophobic amino acids such as glycine and proline, which form mobile hydrophobic regions bounded by crosslinks between lysine residues. Degradation products of the encoded protein, known as elastin-derived peptides or elastokines, bind the elastin receptor complex and other receptors and stimulate migration and proliferation of monocytes and skin fibroblasts. Elastokines can also contribute to cancer progression. Deletions and mutations in this gene are associated with supravalvular aortic stenosis (SVAS) and autosomal dominant cutis laxa. [provided by RefSeq, Aug 2017]