

Product datasheet for **RC239371**

IRAG1 (NM_001206880) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IRAG1 (NM_001206880) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	IRAG1
Synonyms:	IRAG; JAW1L; MRV11
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239371 representing NM_001206880
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTAAAAGCTCCCCAGAGTGAAGAGAGGCTGGCCAGAGGAGGAAAGGAGAATAACTCAGTTTTAGCCT
GTGGAGCCCAGGCCTCTTGGAGCATCTTTGGGGCTGACGCAGCGGAGTTCCGGGCACACGTGGCCACTC
CCAGCAGGAGGCTGCCATGCCCCACATTCGAGGACGAGGAGCCCCCGGAGAGCCACAGGCAGCCAG
AGCCCTGCCGGCCAAGTCTCTGCCGACAGGATATCTTGCAGTCCAACCTCCACGATTGTCTGACTG
GGGATGCCACTTCACCAGAAGGAGAAAACCGACAAAAACCTGGCCAACAGAGTTCACAGTCCCACAAGAG
GCTTTCTACCGACACTGAAGGTGTCCACTGCCTCCCTGACATCTGTGGACCCCGGGGCACATCATT
GACCTGGTGAATGACCAGCTGCCAGACATCAGCATCTCAGAGGAGGACAAGAAGAAAAACCTGGCGCTGC
TGAAGAAGCCAAGTTGGTGAGTGAGCGATTCTGACCCGCGTGGGAGGAAGTCCAGGAGCAGCCCGG
AGACTCCCACATCAGCTGTTTCCCGAACCTCAGCCCCAGCGCTTCTCCTACATCCTCTCGGAGCAACTCA
CTTACAGTCCCCACCCCGCAGGTTTGGATGTGTGCAGTGGCCCGCCATCCCCTCTGCCTGGAGCACCAC
CACAGAAGGGGGATGAGGCCGACGTCTTTCACCTCACCTGGCGAGCCTGAGCAGATCCTGATGAGAAA
TCAGAAGTGTGGGGCTCAAGCTTCCAGACCTTAGTGAAGCAGCTGAGCAGGAAAAAGGGCTTCTCTT
GAACTCTCCCAGCTATTGAGGAAGAAGAGTCAAAGAGTGGCTTAGATGTCATGCCTAATATTTCTGATG
TGCTGTGCGCAAACGCGGGTCCACAGGAGTCTCCCTGGAAGTGCCCTCCACTCACTGAAAAGGAAGT
TGAGAAGCTGTTGTGCAACTGTCTTGGCCTTAGAAATGACAGCTACACTCTGGAATCTAGAATTAAC
CAGGCTGAAAGGGAACGCAACCTGACAGAGGAGAACACTGAGAAAAGAACTGAAAACTTCAAAGCTTCCA
TTACGTCCTCAGTTCACTCTGGCACCACTGTGAGCACCAGGAAACCTACCAGAAGTTGCTGGAGGACAT
CGCTGTCTGCACCGCCTGGCTGCCCGCCTCTCCAGCCGAGCTGAGGTGGTAGGCCCGCTCCGCCAGGAA
AAGCGCATGTGCAAGCAACGGAAGTGTGATGTCAGTATGTGGAGAATCTAAAGAGGACGTATGAGAAGG
ACCATGCGGAGCTCATGGAGTTTAAAAAGCTTGCAAATCAGAATTCAGCCGACGCTGTGGCCCTCTGA
AGATGGGGTCCCTCGCACGGCACGGTCCATGTCCCTCACGCTGGGAAAGAATATGCCTCGCCGAGGGTC
AGCGTTGCTGTGGTTCTAAGTTAATGCCCTGAATCTGCCTGGCCAACTCCCAGCTCATCATCCATT
CCTCCTTACCAGCCTTGTGGAATCACCAATGGGAAAGGCAGCCTACCTGTCACCTCAGCACTGCCTGC
ACTTTTGGAAAAATGAAAGACAAATGGGACCCAGATTGTGAAGCCTCTGCTCTGCGCTGACCCTGAGC
TGCTTGAGGAGCTTAGTCAGGAGACCAAGGCCAGGATGGAGGAAGAAGCCTACAGCAAGGGATTCCAAG
AAGGTCTAAAGAAGACCAAGAAGTCAAGACCTGAAGGAGGAGGAGGAAGAACAGAAGAGTGAGAGTCC
TGAGGAACCTGAAGAGGTAGAAGAACTGAGGAAGAGGAAAAGGGCCCAAGAAGCAGCAAACTTGAAGAA
TTGGTCCATTTCTACAAGTCATGTATCCAACTGTGTGAGCACTGGCAAGTGTCTGGATGATGGCTG
CAGTGTGCTGGTCTTACTGTTGTGCTGGGGCTTACAATTCCTATAACTTTGTGAGAGCAGGCTGA
TGGGCCCTTGAAGATCCACTTGTCTGGCAGCCAGAGGGACTCCTGGTGGAGCTCAGGACTCCAGCAT
GAGCAGCTACAGAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239371 representing NM_001206880
 Red=Cloning site Green=Tags(s)


MVKAPQSEERLARGGKENNSVLACGAQASWSIFGADAAEVPGTRGHSQQEAAMPHIPEDEEPPGEPQAAQ
 SPAGQGPPAAGVSCSPTPTIVLTGDATSPGETDKNLANRVHSPHKRLSHRHLKVSTASLTSVDPAGHII
 DLVNDQLPDISISEEDKKKNLALLEEAKLVSERFLTRRGRKSRSSPGDPSAVSPNLSASPTSSRSNS
 LTVPTPPGLDVCSGPPSPLPGAPPQKGDEADVSSPHGPEPEHILMRNQNLVGLKLPDLSEAAEQEKGLPS
 ELSPAIEEEESKSGLDVMPNISDVLLRKLRVHRSLPGSAPPLTEKEVENVFVQLSLAFRNDYSYTLERIN
 QAERERNLTEENTEKELENFKASITSSASLWHHCEHRETYQKLLIEDI AVLHRLAARLSSRAEVGAVRQE
 KRMSKATEVMMQYVENLKRTYEKDHAELEMFKLANQNSSRSCGPESEDGVPRTARSMSLTGKNMPPRRRV
 SVAVVPKFNALNLPGQTPSSSSIPSLPALSEPNKGSLPVT SALPALLENGKTNGDPDCEASAPAL TLS
 CLEELSQETKARMEEEAYSKGFQEGLEKTKELQDLKEEEEQKSESPPEEVEEETEEEEKGRSSKLEE
 LVHFLQVMYPKLCQHWQVIWMMAAVMLVLTVVGLYNSYNCAEQADGPLGRSTCSAAQRDSWWSSGLQH
 EQPTEQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:


ACCN: NM_001206880

ORF Size: 2118 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001206880.2](#)

RefSeq Size: 5758 bp

RefSeq ORF: 2121 bp

Locus ID: 10335
UniProt ID: [Q9Y6F6](#)
Cytogenetics: 11p15.4
Protein Families: Transmembrane
Protein Pathways: Vascular smooth muscle contraction
MW: 77.3 kDa

Gene Summary: This gene is similar to a putative mouse tumor suppressor gene (Mrvi1) that is frequently disrupted by mouse AIDS-related virus (MRV). The encoded protein, which is found in the membrane of the endoplasmic reticulum, is similar to Jaw1, a lymphoid-restricted protein whose expression is down-regulated during lymphoid differentiation. This protein is a substrate of cGMP-dependent kinase-1 (PKG1) that can function as a regulator of IP3-induced calcium release. Studies in mouse suggest that MRV integration at Mrvi1 induces myeloid leukemia by altering the expression of a gene important for myeloid cell growth and/or differentiation, and thus this gene may function as a myeloid leukemia tumor suppressor gene. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene, and alternative translation start sites, including a non-AUG (CUG) start site, are used. [provided by RefSeq, May 2011]