

Product datasheet for **RC239386**

NFAT2 (NFATC1) (NM_001278670) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NFAT2 (NFATC1) (NM_001278670) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NFATC1
Synonyms:	NF-ATC; NF-ATc1.2; NFAT2; NFATc
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239386 representing NM_001278670
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCAAGCACCAGCTTCCAGTCCCTTCCAAGTTTCACTTGGCCCTGCGGCTGCGGTCTTCGGGAGAG
GAGAACTTTGGGGCCCGCGCCGCGCGCCGGCGCACCATGAAGTCAGCGGAGGAAGAACAATATGGCTA
TGCATCTCCAACGTCAGCCCCGCCCTGCCGCTCCCCACGGCGCACTCCACCCTGCCGGCCCCGTGCCAC
AACCTTCAGACCTCCACACCGGGCATCATCCCGCCGGCGGATCACCCCTCGGGGTACGGAGCAGCTTTGG
ACGGTGGGCCCGGGCTACTTCTCTCTCCGGCCACACCAGGCCTGATGGGGCCCTGCCCTGGAGAG
TCCTCGCATCGAGATAACCTCGTGTGGGCTGTACCACAACAATAACCAGTTTTTCCACGATGTGGAG
GTGAAGACGTCCTCCCTAGCTCAAACGGTCCCCTCCACGGCCACGCTGAGTCTGCCAGCCTGGAGG
CCTACAGAGACCCCTCGTGCCTGAGCCCGGCCAGCAGCCTGTCTCCGGAGCTGCAACTCAGAGGCCTC
CTCTACGAGTCCAACACTACTCGTACCCGTACGCGTCCCCCAGACGTCGCCATGGCAGTCTCCCTGCGTG
TCTCCAAGACCACGGACCCCGAGGAGGGCTTTCCCGCGGGCTGGGGGCTGCACACTGCTGGGTCCC
CGCGGCACTCCCCCTCCACCTCGCCCCGCGCCAGCGTCACTGAGGAGAGCTGGTGGGTGCCCGTCTCTC
CAGACCCCGCTCCCCTTGCAACAAGAGGAAGTACAGCCTCAACGGCCGGCAGCCGCCCTACTCACCCAC
CACTCGCCACGCCGTCCCCGCACGGCTCCCCGCGGGTCAAGCTGACCGACGACTCGTGGTTGGGCAACA
CCACCCAGTACACCAGCTCGGCCATCGTGGCCGCCATCAACGCGCTGACCACCGACAGCAGCCTGGACCT
GGGAGATGGCGTCCCTGTCAAGTCCCGCAAGACCACCTGGAGCAGCCGCCCTCAGTGGCGCTCAAGGTG
GAGCCCGTGGGGAGGACCTGGCAGCCCCCGCCCCGGCGGACTTCGCGCCGAAGACTACTCCTCTT
TCCAGCACATCAGGAAGGGCGGCTTCTGCGACCAGTACCTGGCGGTGCCGAGCACCCCTACCAGGCGC
GAAGCCCAAGCCCTGTCCCCTACGTCTACATGAGCCCGACCTGCCCGCCCTGGACTGGCAGCTGCCG
TCCCACTCAGGCCGTATGAGCTTCGGATTGAGGTGCAGCCCAAGTCCCACCACCGAGCCCACTACGAGA
CGGAGGGCAGCCGGGGGCGTGAAGGCGTCCGCCGGAGGACACCCCATCGTGCAGCTGCATGGTACTT
GGAGAATGAGCCGCTGATGCTGCAGTTTTTCAATTGGGACGGCGGACGACCGCCTGCTGCGCCGCACGCC
TTCTACCAGGTGCACCGCATCACAGGGAAGACCGTGTCCACCACCGCCAGGAGCCATCCTCTCCAACA
CCAAAGTCTGGAGATCCCACTCCTGCCGGAGAACAGCATGCGAGCCGTCATTGACTGTGCCGAATCCT
GAAACTCAGAACTCCGACATTGAACTTCGAAAGGAGAGACGGACATCGGGAGGAAGAACAACACGGGTA
CGGCTGGTGTCCGCTTACAGTCCCGCAACCCAGCGCCGCACGCTGTCCCTGCAGGTGGCCTCCAACC
CCATCGAATGCTCCAGCGCTCAGCTCAGGAGCTGCCTCTGGTGGAGAAGCAGAGCACGGACAGCTATCC
GGTCTGGGCGGGAAGAAGATGGTCTGTCTGGCCACAACCTCCTGCAGGACTCCAAGGTCATTTTCGTG
GAGAAAGCCCCAGATGGCCACCATGTCTGGGAGATGGAAGCGAAAAGTACCAGGGACCTGTGCAAGCCGA
ATTCTCTGGTGGTTGAGATCCCGCCATTTCCGAAATCAGAGGATAACCAGCCCCGTTACAGTCAAGTTCTA
CGTCTGCAACGGGAAGAGAAAGCGAAGCCAGTACCAGCGTTTACCTACCTTCCCGCCAACGTAATGAA
ATAATACGAAATGACCTCTCCAGCACGAGCACCCACTCC

AG**GCGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239386 representing NM_001278670
 Red=Cloning site Green=Tags(s)


MPSTSFVPVPSKFPLGPAAAVFRGETLGPAPRAGGTMKSAEEEHYGYASSNVSPALPLPTAHSTLPAPCH
 NLQTSTPGIIPPADHPSGYGAALDGGPAGYFLSSGHTRPDGAPALESPRIEITSCLGLYHNNNQFFHDVE
 VEDVLPSSKRSPSTATLSLPSLEAYRDPSCSPASSLSSRSCNSEASSYESNYSYPYASPQTSPWQSPCV
 SPKTTDPEEGFPRGLGACTLLGSPRHSPSTSPRASVTEESWLGARSSRPASPCNKRKYSLNQRQPPYSPH
 HSPTSPHGSPPRVSVTDDSWLGNTTQYTSSAIVAAINALTTDSSLDLGDGVPVKS RKTTLEQPPSVALKV
 EPVGEDLGSPPPADFPEDYSSFQHIRKGGFCQYLAVPQHPYQWAKPKPLSPTSYSMP TLPALDWQLP
 SHSGPYELRIEVQPKSHHRAHYETEGSRGAVKASAGGHPVQLHGYLENEPLMLQLFIGTADDRLLRPHA
 FYQVHRITGKTVSTTSHEAILSNTKVL EIPLLPENS MRAVIDCAGILKLRNSDIELRKGETDIGRKNTRV
 RLVFRVHVPQPSGRTL SLQVASNPIECSQRSAQELPLVEKQSTDSYPVVGKKMVL SGHNFLQDSKVIFV
 EKAPDGHVWEMEAKTRDLCKPNSLVVEIPFRNQRITSPVHVSFYVCNGKRRKSQYQRF TYLPANVNE
 IIRNDLSSTSTHS

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:


ACCN: NM_001278670

ORF Size: 2139 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278670.2](#)

RefSeq Size: 4341 bp

RefSeq ORF: 2142 bp

Locus ID:	4772
UniProt ID:	O95644
Cytogenetics:	18q23
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Axon guidance, B cell receptor signaling pathway, Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway
MW:	77.8 kDa
Gene Summary:	<p>The product of this gene is a component of the nuclear factor of activated T cells DNA-binding transcription complex. This complex consists of at least two components: a preexisting cytosolic component that translocates to the nucleus upon T cell receptor (TCR) stimulation, and an inducible nuclear component. Proteins belonging to this family of transcription factors play a central role in inducible gene transcription during immune response. The product of this gene is an inducible nuclear component. It functions as a major molecular target for the immunosuppressive drugs such as cyclosporin A. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. Different isoforms of this protein may regulate inducible expression of different cytokine genes. [provided by RefSeq, Jul 2013]</p>