

Product datasheet for **RC239394**

ADAM2 (NM_001278113) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADAM2 (NM_001278113) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ADAM2
Synonyms:	CRYN1; CRYN2; CT15; FTNB; PH-30b; PH30; PH30-beta
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide
Sequence:

>RC239394 representing NM_001278113
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTGGCGCTTGTGTTCTGCTCAGCGGGCTCGCGGGCTCGGATGGACAGTAATTTTATAGTTTAC
CTGTGCAAATTACAGTTCGGAGAAAATACGGTCAATAATAAAGGAAGGAATTGAATCGCAGGCATCCTA
CAAAATTGTAATTGAAGGAAAACCATATACTGTGAATTTAATGCAAAAAAATTTTTACCCATAATTTT
AGAGTTTACAGTTATAGTGGCACAGGAATTATGAAACCACTTGACCAAGATTTTCAGAATTTCTGCCACT
ACCAAGGGTATATTGAAGTTATCCAAAATCTGTGGTATGGTTAGCACATGACTGGACTCAGGGGCGT
ACTACAGTTTGAAAATGTTAGTTATGGAATAGAACCCTGGAGTCTTCAGTTGGCTTTGAACATGTAATT
TACCAAGTAAAACATAAGAAAGCAGATGTTTCTTATATAATGAGAAGGATATTGAATCAAGAGATCTGT
CCTTTAAATTCAAAGCGTAGAGTATAATCATATGGGGTCTGATACAACCTGTTGTCGCTCAAAAAGTTT
CCAGTTGATTGGATTGACGAATGCTATTTTTGTTTCATTAATATTACAATTATTCTGTCTTCATTGGAG
CTTTGGATAGATGAAAATAAAATTGCAACCACTGGAGAAGCTAATGAGTTATTACACACATTTTTAAGAT
GAAAAACATCTTATCTTGTTTTACGTCCTCATGATGTGGCATTTTTACTTGTTTACAGAGAAAAGTCAAA
TTATGTTGGTGCAACCTTCAAGGGAAGATGTGTGATGCAAACCTATGCAGGAGGTGTTGTTCTGCACCC
AGAACCATAAGTCTGGAATCACTTGCACTTATTTAGTCAATTAATTGAGCCTTAGTATGGGGATCACTT
ATGATGACATTAACAAATGCCAGTGCTCAGGAGCTGTCTGCATTATGAATCCAGAAGCAATTCATTTAG
TGGTGTGAAGATCTTTAGTAAGTGCAGCTTCAAGACTTTGCACATTTTATTTCAAAGCAGAAGTCCCAG
TGCTTTCACAATCAGCCTCGCTTAGATCCTTTTTCAAACAGCAAGCAGTGTGTGTAATGCAAAGCTGG
AAGCAGGAGAGGAGTGTGACTGTGGGACTGAACAGGATTGTGCCCTTATTGGAGAAAACATGCTGTGAT
TGCCACATGTAGATTTAAAGCCGGTTCAAACCTGTGCTGAAGGACCATGCTGCGAAAACATGCTATTATG
TCAAAAAGAAAGAATGTGTAGGCCTTCTTTGAAGAATGCGACCTCCCTGAATATTGCAATGGATCATCTG
CATCATGCCAGAAAACCACTATGTTTCAAGACTGGGCATCCGTGTGGACTGAATCAATGGATCTGTATAGA
TGGAGTTTGTATGAGTGGGGATAAAACAATGTACAGACACATTTGGCAAAGAAGTAGAGTTTGGCCCTTCA
GAATGTTATTCTCACCTTAATTCAAAGACTGATGTATCTGAAACTGTGGTATAAGTGAATTCAGGATACA
CACAGTGTGAAGCTGACAATCTGCAGTGCAGAAAATTAATATGTAATATGTAGGTAATTTTTATTACA
AATTCAGAGCCACTATTTTATGCCAACATAAGTGGACATCTCTGCATTGCTGTGGAATTTGCCAGT
GATCATGCAGACAGCCAAAAGATGTGGATAAAAGATGGAACCTTCTGTGGTTCAAATAAGGTTTGCAGGA
ATCAAAGATGTGTGAGTTCTTCACTTGGGTTATGATTGTAAGTACTGACAAATGCAATGATAGAGGTGT
ATGCAATAACAAAAGCACTGTCAGTGTAGTGCTTCAATTTACCTCCAGATTGCTCAGTTCAATCAGAT
CTATGGCCTGGTGGGAGTATTGACAGTGGCAATTTCCACCTGTAGTATACCAGCCAGACTCCCTGAAA
GGCGCTACATTGAGAACATTTACCATTCAAACCAATGAGATGGCCATTTTTCTTATTCTCTTTCTT
TATTATTTCTGTGACTGATTGCTATAATGGTGAAAGTTAATTTCAAAGGAAAAAATGGAGAAGTGAAG
GACTATTCAAGCGATGAGCAACCTGAAAGTGAAGTGAACCTAAAGGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001278113

ORF Size: 2148 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278113.2](#)

RefSeq Size: 2627 bp

RefSeq ORF: 2151 bp

Locus ID:	2515
UniProt ID:	Q99965
Cytogenetics:	8p11.22
Protein Families:	Druggable Genome, Protease, Transmembrane
MW:	80.6 kDa
Gene Summary:	<p>This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The encoded protein is a subunit of an integral sperm membrane glycoprotein called fertilin, which plays an important role in sperm-egg interactions. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2013]</p>