

Product datasheet for **RC239405**

KIAA1586 (NM_001286276) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1586 (NM_001286276) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KIAA1586
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239405 representing NM_001286276
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCAAAACGACAGGGTGATTTTTTGCATTTTTTAAATGTGAAGAAGGTGAAAACAGACACAGAAAATA
 ATGAAGTGAGCAAAAATCACTGCAGATTGTCTAAGGCAAGGAACACATTTTCGAGTATATTGAACAACC
 AATCATTTGAAGAAAAGCCATCACTTTCATCAAAGAAAAGAAATAGATAATCTTGTGCTTCCAGATTGTTGG
 AATGAAAAACAAGCATTTATGTTTACAGAACAAATACAAATGGCTTGAAATAAAAGAAGGTAAATTAGGAT
 GTAAGGATTGTTTCAGCAGTTTCGGCATTTGGGATCGAAAGCAGAAAAGCATGTCCATGTGTCCAAGGAATG
 GATTGCATATTTAGTAACCCCTAATGGCAGTAATAAACTACTAGGCAAGCTTCTCTACGAAAAAAATTT
 AGGGAACATGATGTTTCTAAAGCCCATGGTAAAATTCAGGATTTGTTAAAGGAATCAACTAATGATTCAA
 TTTGTAATTTAGTGCATAAACAAAATAATAAAATATTGATGCTACTGTAAAAGTTTTCAACTACTGTTTA
 CAGTTTAGTAAACATAACAGACCTTTATCTGATATTGAGGGGGCAGAGAATTACAGGAAAAAAATGGA
 GAGGTAAATTGTTTAAATACACGTTACAGTGCAACAAGAATAGCAGAACATATTGCAAAAGAAATGAAGA
 TGAAGATATTTAAGAATATTATAGAAGAGAATGCCAAAATCTGTATCATAATTGATGAGGCATCTACAGT
 TTCAAAGAAAACCACTAGTGATTATCTCCAGTGACAATTCAGTCAGCTCCTGCACCTGTTATGTTA
 TTTGTGGCTTTAAAGAATTGGTGTCAACTATAGCAGAGTGATTGTCAATACATTATTGACTACTTTAA
 ATGATTGTGGTTTTACAAATGAATATTTGAAAGCAAAATTAATTGCATTTTGTCTGATGGTGCTAATAC
 AATCCTGGGAAGAAAGTCTGGAGTAGCTACAAAATTTGTAGAAAATTTTCTGAAATCATCATTTGGAAC
 TGTTTAAATCATCGATTACAATTGTCACTTGATGATTCTATATCCGAAATAAAACAAATTAATCATTTAA
 AAATATTTATTGATAAAATTTATTCTATTTATCATCAACCTAATAAAAAATCAAACCAAGCTTCTAGGAAC
 GTAGCTAAAGAACTTGAAACTGAAATTATTTAAATTTGGTCGAGTAATGGGACCAAGATGGGCGGCATGT
 AGTTTACAAGCTGCTACTGCTGTATGGCATGCATATCCTATATTATATATGCATTTTTCTCATTCTTACT
 CTGGTTTGGCGAAGAGATTAGCTAACATTAATTTCTTACAAGACCTTGCTTTAATGATTGACATTCTTGA
 AGAATTTTCAGTACTTTCAACTGCATTACAGTCAAGATCAACTAATATTAAGAAAGCAGAAAAATGATC
 AAACGTACCATAAGAGCTTTGAAAAATTTAAAAATTTGGTACTGGAAAGTATGAATCTCAAATGAAGATT
 TGATCAAGTCAGATAAGTTTAAAGATATTCCATTTAATAAAAAACAATAAATTTAATGCTCTTCTAGGAG
 TATATTACTAGACAATATAATTCAGCACATGAACCTACGCCTTTTATCTGACAGAAACCATGAAGATATT
 TTTAATTACTTTGATTTGCTGGAACCTCCACATGGCCTTATGAAGAAATACTTCACCATGGATAGCTG
 GTGAAAAAACATTATTTCAATTTGTGTAATTTTAAAAATATGAAGTTGATTTGAATGATTTTCGGGAATT
 TGTAATAATAATATAAAATCAAACAATGTTTCAATTCCTACAACCTATATACAAAGCTAAAAAGATAGTT
 AGCACCATTGCAATCAATAGTGCTGAAGCTGAAAGGGGTTTCAATTTAATGAACATAATTTGTACAAGGG
 TGAGAAATAGTTTAAACAATAGATCATGTATCAGATTTAATGACAATAAATTTACTGGGGAAGAATTAGC
 AGATTGGGATGCAACACCGTTTGTAAAATCTTGGTCAAATTGCAACCACAGGTTGGCTACAGATACAAGA
 GTTCGGCAAAAGTCAACAAAAGCTTCCATGAGAATCAATTGGCTATATGGAACCTTAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239405 representing NM_001286276
 Red=Cloning site Green=Tags(s)


MPKRQGDFLHFLNVKKVKTDTENNEVSKNHCRLSKAKEPHFEYIEQPIIEEKPSLSSSKEIDNLVLPDCW
 NEKQAFMFTEQYKWLIEKEGKLCKDCSAVRHLGSKAEKHVHVSKEWIAVLVTPNGSNKTTQASLRKKI
 REHDVSKAHGKIQDLLKESTNDSICNLVHKQNNKIDATVKVFNTVYSLVKHNRPLSDIEGARELQEKNG
 EVNCLNTRYSATRIAEHIAKEMKMKIFKNIEENAKICIIIDEASTVSKKTTLVLYLQCTIQSAPAPVML
 FVALKELVSTIAECIVNTLLTTLNDCGFTNEYLKANLIAFCSDGANTILGRKSGVATKLLLENFPEIIWN
 CLNHRLQLSLDDSISEIKQINHLKIFIDKIYSIYHQPKNQTKLLGTVAKELETEIIGRVMGPRWAAC
 SLQAATAVWHAYPILYMHFSHYSGLAKRLANINFLQDLALMIDILEEFSVLSTALQSRSTNIKKAQKLI
 KRTIRALENLKIGTGKYESQIEDLIKSDKFKDIPFNKNKFNALPRSILLDNIIQHMLRLLSDRNHEDI
 FNYFDLLEPSTWPYEEITSPWIAGEKTLFHLCKILKYEVDLNDREFVNNNIKSNNVSIPTTIYKAKKIV
 STIAINSAAERGFNLMNIICTRVRNSLTIDHVSIDLMTINLLGKELADWDATPFVKSWSNCHRLATDTR
 VRQKSTKVFHENQLAIWNLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001286276

ORF Size: 2160 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001286276.1](#), [NP_001273205.1](#)

RefSeq Size: 2894 bp

RefSeq ORF: 2163 bp

Locus ID: 57691


UniProt ID: [Q9HCI6](#)

Cytogenetics: 6p12.1

MW: 82.9 kDa

Gene Summary: E3 SUMO-protein ligase; facilitates UBE2I/UBC9-mediated SUMO2 modification of target proteins (PubMed:26524493).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239405