

Product datasheet for **RC239433**

CCDC108 (CFAP65) (NM_001278296) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC108 (CFAP65) (NM_001278296) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CFAP65
Synonyms:	CCDC108; SPGF40
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239433 representing NM_001278296
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGCTCACCCAGGCTCCAAGCTCCGTCGTGAGGTCCAGGAACAGCAGGAACACACCGTGAACCTGTG
 GTGGATCCTGCCTGAGTGCCAGCACAGTGGCCATCCCTGCCATCAACGACAGCAGTGCAGCCATGAGTGC
 CTGCAGCACCATCAGCGCCAGCCGCAAGCTCCATGGACTCAGATGCACTCCCAAAGAAGCAGGAG
 AGAGTGAACAAGAGGGTTCATCTGGGGCATTGAGGTGGCTGAGGAGCTGCATTGAAAGGCTGGGAGCTAG
 GAAAGGAGACCACAAGGAATCTGGTTCTGAAAAATCGATCCTTAAACTCCAGAAGATGAAGTACAGGCC
 CCCCAGACCAAGTTCTTCTCACGGTCATCCCTCAGCCCATCTTCTGAGCCCAGGCATAACCCCTCACG
 CTCCCCATCGTCTCCGGCCTCTGGAGGCGAAGGAGTACATGGACCAGCTGTGGTTTGAGAAAGCGGAGG
 GGATGTTCTGTGTCGGCCTACGGGCCACCCTGCCACAGGCTGATCTGCCGCCACCATCCCTGCA
 GCTGCCCATGTGTGCTGTGGGAGATACGACTGAGGCCTTTTTCTGCCTGGATAATGTGGGGACCTGCC
 ACCTTCTCACCTGGGAGTCTCCAGCCATTCCAGATGCTGCCCGCCACGGGGCTCTGGAGCCAGGCC
 AGGCCTCTCAGATCAAGGTGACCTTTCAGCCCTTACAGCCGTCATCTACGAGGTGCAGGCCACGTGCTG
 GTACGGGGCGGGCAGCCGGCAGAGGAGCAGCATCCAGCTGCAGGCTGTGGCAAGTGCGCCCAGCTGCTG
 GTGAGCATAAAGCACAAGTGCCTGGAGGACCAGGATGCCGAGGGCTTCCAGAAGCTGTTGTACTTTGGCT
 CTGTTGCTGTGGGCTGCACCTCGGAGAGGCAGATCAGGCTACACAACCCGTCGGCGGTAATGCCCCCTT
 CAGGATTGAAATTTCCCGGATGAACTGGCCGAAGACCAGGCCTTCTCATGCCCCACGGCCCATGGCATC
 GTGCTTCCGGGAGAGAAGAAATGTGTGCGTGTCTTCCACCCCAAGACTCTGGACACCAGAAGTGTGG
 CCCTGCTGTGCCCTGCAGCACTACTGTGTCAACTTCAGCTGGTCAACCTTGGGAGCGCTCCGAGCAG
 CCCCTGTGGATTGAGAACCAATCGGACTGCACGGCCACTTCCAGTTTCCATCGACTGCTTGGAGAGTG
 TCTTTACCATCAGGCCTGCCTTTGGGACGCTGGTGGGCAAGGCCGATGACCCTGCACTGTGCCTTCCA
 GCCCACTACCCCATCATCTGCTTTCCGCGTGTGGCCTGTCTCATCCACCACCAGGACCCACTGTTCTG
 GACCTGATGGGGACCTGCCACTCGGACAGCACCAAGCCAGCCATCCTGAAGCCTCAGCACCTCACCTGGT
 ACCGCACACACCTGGCCCGGGCCTGACGCTTACCCCTGACATCCTGGATGCCATGCTGAAGGAGAA
 GAAGCTGGCACAGGACCAGAACGGGCTCTCATGATCCCATCCAGGATCTGGAGGACATGCCGGCCCCG
 CAGTACCCTTATATCCCCCATGACCGAGTTCTTCTTCGACGGCACCAGCGACATAACCATCTTCCCC
 CGCCCATCAGTGTAGAGCCTGTCGAGGTAGACTTCGGTGCCTGCCAGGGCCTGAGGCCCAACCTGT
 ACCCTGTGCCTGATGAACCACACCAAGGCAAGATCATGGTGGTCTGGACGCGAAGGTCTGACTGCCCC
 TTCTGGGTGACTCCAGAGAGCTGCGACGTGCCCCCACTCAAGTCCATGGCCATGCGCCTGCACTTCCAGC
 CGCCTACCCCAACTGCCTTTACACGGTGGAGCTCGAAGCCTTCGCCATCTATAAGGTGTGTGCACGCAA
 TGAGAGGGAGGAATGCGGGTCTCTGCTAGGAGCCTGAGTGGCTTGGTGGGGTGGCAGGAAGTGACCGAG
 GGCAGCTTCAAGCTCCATCCTCTGCGTGCCAGGCTTCTCTTGGCTGGACAGTGACCCCTATGAGTTGT
 CTCCTCAAAGCTCCTGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239433 representing NM_001278296
 Red=Cloning site Green=Tags(s)


MMLTQAPSSVVRSRNHTVNSGGSCLSASTVAIPAINSSAAMSACSTISAQPASSMDTQMHSPPKQE
 RVNKRVIWGIIEVAEELHWKGWELGKETTRNLVLKNRSLKQKMKYRPPKTKFFFTVIPQIFLSPGITLT
 LPIVFRPLEAKEYMDQLWFEKAEGMFCVGLRATLPCRHLICRPPSLQLPMCAVGDTEAFFCLDNVGDLP
 TFFTWEFSSPFQMLPATGLLEPGQASQIKVTFQPLTAVIYEVQATCWYGAGSRQRSSIQLQAVAKCAQLL
 VSIKHKCPEDQDAEGFQKLLYFGSVAVGCTSERQIRLHNPSAVNAPFRIEISPDELAEDQAFSCPTAHGI
 VLPGEKKCVSVFFHPKTLDRTRVDYCSIMPSGCASKTLLKVVGFRCRPAVSLQHYCVNFSWVNLGERSEQ
 PLWIENQSDCTAHFQFAIDCLESVFTIRPAFGTLVGKARMTLHCAFQPTHPIICFRRVACLIHHQDPLFL
 DLMGTCHSDSTKPAILKPQHLTWYRTHLARGLTLYPPDILDAMLKEKLAQDQNGALMIPIQDLEDMPAP
 QYPYIPPMTEFFFDGTSITIFPPPISVEPVEVDFGACPGPEAPNPVPLCLMNHTKGIKIMVWTRRSDCP
 FWVTPESCDVPLKSMAMRLHFQPPHPNCLYTVELEAFIYKVCARNERECEGVSARSLSGLVGWQEVTE
 GSFRLHPLRARLSLGTWVTPMSLSPKLLA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001278296

ORF Size: 2190 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278296.2](#)

RefSeq Size: 3287 bp

RefSeq ORF: 2193 bp

Locus ID: 255101

Cytogenetics:	2q35
Protein Families:	Transmembrane
MW:	81.8 kDa
Gene Summary:	The protein encoded by this gene has putative coiled-coil domains and may be a transmembrane protein. The chicken ortholog of this gene is involved in the Rose-comb mutation, which is a large chromosome inversion, resulting in altered comb morphology and defects in sperm motility. [provided by RefSeq, Aug 2016]