

Product datasheet for **RC239454**

ZNF84 (NM_001289972) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF84 (NM_001289972) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF84
Synonyms:	HPF2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC239454 representing NM_001289972
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGACCATGTTACAGGAGTCATTCTCATTGACGATTTATCTGTGGACTTCACCCAAAAGGAGTGGCAGC
 TACTGGATCCCTCTCAGAAGAATTTATACAAGGATGTGATGTTGGAGAAGTATAGCAGCCTAGTGTCAC
 GGGGTATGAAGTTATGAAACAGATGTCATCTTCAAATTGGAGCAAGGAGAAGAGCCGTGGGTAGGAGAT
 GGAGAAATCCAAGTTCAGATTCTCCAGAAGTCTGAAAGTAGATGGTAACATGATGTGGCACCAGGATA
 ACCAAGACAAGCTTAAAAATTAAAAAGAGGTCATGAATGTGATGCATTTGGAAGAAATTTCAATCTGAA
 CATGAACCTTTGTTCTTTAAGGAAATCAACAGTGAAGGTGACTTAGATGGATTGATTTTAAACATCAT
 TTAGATTTGCTTATTCAAAAGGAGATTATGGAAGCAGAATCAGATGACTTTAATGTGTTTGATAATT
 TTTTCTCCATTCCAAGCCTGAGGATACTGATACCTGGTTAAAAACTATGACTGTGATAAATAAAGA
 GAGCTATAAAAGTCACAGATTATCATATATCATAGAAATCGTTAGGGGAGAACTCTATGAATGCAGT
 GAATGTAGGAAGCGCTTCAGTAAGAAACCAAGTCTCATTAAACATCAGAGCAGACATATAAGAGACATAG
 CCTTTGGCTGTGGTAATTGTGGCAAAACCTTTCCCAAGAGTCTCAGTTTATTACACATCAGAACTCA
 TACAGGAGAAAAACCTTATAATTGTAGCCAGTGTGGGAAAGCCTTCTCCAAAAGTCACAGCTCACATCC
 CATCAGCGGACACATACAGGAGAGAAACCTTATGAGTGTGGTGAATGTGGGAAAGCCTTCTCCCGAAGT
 CACATCTCATATCGCATTGGAGAACACACAGGAGAGAAACCTTATGGATGCAATGAATGTGGGAGGGC
 CTTTAGTGAAGTCCAATCTCATTAAACATCAGAGAATTCATACCGGTGAGAAGCCTTTTGAATGCAGG
 GAATGTGGGAAAGCCTTCAGCAGGAAGTCACAACTCGTTACACATCAGAACTCACACAGGAACAAAAC
 CCTTTGGATGTAGTGATTGTAGAAAAGCATTCTTTGAGAAGTCAGAGCTTATTAGACATCAGACAATTCA
 TACTGGAGAGAAACCTATGAATGCAGCGAGTGTAGGAAAGCATTAGAGAGAGGTGAGTCTCATTAAAT
 CATCAGAGAACACATACAGGAGAGAAACCTCATGGATGCATTGAGTGTGGGAAAGCCTTCTCCAGAAAT
 CACATCTCATATCACATCAGATGACACACAGGAGAGAAACCTTTATATGCAGTAAATGTGGGAAAGC
 CTTTACAGGAGAAATCACAGCTCGTTAGACATCAGAGAAGTATACGGGAGAGAAACCGTATGAATGCAGT
 GAGTGTGGGAAAGCTTTCAGTGAAGAAATTAAGTCTCACTAATCATCAAAGAATTCATACAGGAGAAAAAC
 CATATGTATGCAGTGAATGTGGGAAAGCCTTTTGTGAGAAGTCACATCTCATATCACATCAGAGGACACA
 TACAGGGGAGAAACCTATGAATGCAGTGAATGTGGGAAAGCCTTTGGTGAAGTCAAGTCTTGCAACT
 CATCAGAGAAGTCACTGGAGAAAAACCGTATGAATGCAGGAGTGTGAAAAAGCTTTCTCCAGAAAT
 CACAGCTAAATACCCATCAGAGAATTCACACTGGAGAGAAACCTATGAATGCAGTCTTTGTAGGAAAGC
 TTTTTTGTAGAAGTCGGAGCTAATTAGACATCTGAGAAGTATACAGGAGAAAAACCTTATGAATGCAAT
 GAATGTAGAAAAGCCTTCAGGGAGAAGTCAAGTCTCATCAATCATCAGAGAATACATACAGGAGAGAAGC
 CTTTTGAATGCAGTGAATGTGGCAAGCCTTCTCTCGGAAGTCACACCTTATACCACATCAAAGGACACA
 TACGGGTGAGAAACCTATGGATGCAGTGAATGTAGGAAGGCCTTCTCTCAGAAGTCACAGCTGGTTAAT
 CATCAGAGAATTCATACAGGAGAGAAGCCTTATCGATGCATTGAATGTGGGAAAGCCTTTCTCACAGAAGT
 CACAGCTCATCAATCATCAGAGAAGTATACAGTAAAAAATCC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239454 representing NM_001289972
 Red=Cloning site Green=Tags(s)


MTMLQESFSFDDLSVDFTQKEWQLLDPSQKNLYKDVMLENYSSLVSLGYEVMKPDVIFKLEQGEEPWVG
 GEIPSSDSPEVWKVDGNMMWHDNDQDKLIIKRGHECDAFGKNFNLNMNFVPLRKSNSGDLGLILKHH
 LDLLIPKGDYGAESDDFNVDNFLLHSPEDTDWLKYYDCDKYKESYKKSQIIYHRNRLGEKLYECS
 ECRKRFSSKPSLIKHSRHIRDIAFGCGNCGKTFPQKSQFITHHRTHTGEKPYNCSCQCKAFSQQSQTLS
 HQRTHTEKPYECGECGKAFSRKSHLISHWRTHTEKPYGCNECGRAFSEKSNLINHQRTHTEKPFECR
 ECGKAFSRKSQLVTHHRTHTGTPFGCSDCRKAFFEKSELIRHQTHTGEKPYECSECRKAFFRSSLIN
 HQRTHTEKPHGCIQCGKAFSQQSHLISHQMTHTGEKPFICSKCGKAFSRKSQLVRHQRTHTEKPYECS
 ECGKAFSEKLSLTNHQRTHTEKPYVCSECGKAFQKSHLISHQRTHTEKPYECSECGKAFGEKSSLAT
 HQRTHTEKPYECRDCEKAFSQQSQLNTHQRTHTEKPYECSLCRKAFFEKSELIRHLRTHTEKPYECN
 ECRKAFREKSSLINHQRTHTEKPFECSECGKAFSRKSHLIPHQRTHTEKPYGCSECRKAFSQQSQLVN
 HQRTHTEKPYRCIECGKAFSQQSQLINHQRTHTVKKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001289972

ORF Size: 2214 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001289972.1](#), [NP_001276901.1](#)

RefSeq Size: 7175 bp

RefSeq ORF: 2217 bp

Locus ID: 7637

UniProt ID: [P51523](#)


Cytogenetics: map from Rosati ref via FISH [AFS]

Protein Families: Transcription Factors

MW: 85.9 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239454