

Product datasheet for **RC239477**

Periostin (POSTN) (NM_001286666) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Periostin (POSTN) (NM_001286666) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	POSTN
Synonyms:	OSF-2; OSF2; PDLPOSTN; PN
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239477 representing NM_001286666
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGATTCCCTTTTACCCATGTTTTCTCTACTATTGCTGCTTATTGTTAACCCATAAACGCCAACAAATC
 ATTATGACAAGATCTTGGCTCATAGTCGTATCAGGGTTCGGACCAAGGCCCAAATGTCTGTGCCCTTCA
 ACAGATTTGGGCACAAAAAGAAATACTTCAGCACTTGAAGAAGTGGTATAAAAAAGTCCATCTGTGGA
 CAGAAAACGACTGTGTTATATGAATGTTGCCCTGGTTATATGAGAATGGAAGGAATGAAAGGCTGCCAG
 CAGTTTGGCCATTGACCATGTTTATGGCACTCTGGGCATCGTGGGAGCCACCACAACGCAGCGCTATTC
 TGACGCCTCAAACTGAGGGAGGAGATCGAGGGAAAGGGATCCTTCACTTACTTTGCACCGAGTAATGAG
 GCTTGGGACAACCTGGATTCTGATATCCGTAGAGGTTTGGAGAGCAACGTGAATGTTGAATTACTGAATG
 CTTTACATAGTCACATGATTAATAAGAGAATGTTGACCAAGGACTTAAAAATGGCATGATTATTCCTTC
 AATGTATAACAATTTGGGGCTTTTCATTAACCATTATCCTAATGGGGTGTCACTGTTAATTGTGCTCGA
 ATCATCCATGGGAACAGATTGCAACAAATGGTGTGTCCATGTCAATTGACCGTGTGCTTACACAAATTG
 GTACCTCAATTCAGACTTCATTGAAGCAGAAGATGACCTTTCATCTTTAGAGCAGCTGCCATCACATC
 GGACATATTGGAGGCCCTTGAAGAGACGGTCACTTCACACTCTTTGCTCCCACCAATGAGGCTTTTGGAG
 AAACCTCCACGAGGTGCTCTAGAAAGGATCATGGGAGACAAAAGTGGCTTCCGAAGCTCTTATGAAGTACC
 ACATCTTAAACTCTCCAGTGTCTGAGTCTATTATGGGAGGAGCAGTCTTTGAGACGCTGGAAGGAAA
 TACAATTGAGATAGGATGTGACGGTGACAGTATAACAGTAAATGGAATCAAAATGGTGAACAAAAAGGAT
 ATTGTGACAAATAATGGTGTGATCCATTTGATTGATCAGGTCCTAATTCCTGATTCTGCCAAACAAGTTA
 TTGAGCTGGCTGGAAAACAGCAAACCACCTTACGGATCTTGTGGCCCAATTAGGCTTGGCATCTGCTCT
 GAGGCCAGATGGAGAATACACTTTGCTGGCACCTGTGAATAATGCATTTTCTGATGATACTCTCAGCATG
 GATCAGCGCTCCTTAAATTAATTCTGCAGAATCACATATTGAAAGTAAAAGTTGGCCTTATGAGCTTT
 ACAACGGGCAAATACTGGAAACCATCGGAGGCAACAGCTCAGAGTCTTCGTATATCGTACAGCTGTCTG
 CATTGAAAATTCATGCATGGAGAAAGGAGTAAGCAAGGGAGAAAACGGTGCATTACATATTCGCGGAG
 ATCATCAAGCCAGCAGAGAAATCCCTCCATGAAAAGTTAAAACAAGATAAGCGCTTAGCACCTTCTCA
 GCCTACTTGAAGCTGCAGACTTGAAGAGCTCCTGACACAACCTGGAGACTGGACATTATTTGTGCCAAC
 CAATGATGCTTTTAAAGGAATGACTAGTGAAGAAAAAGAAATTCGATACGGGACAAAAATGCTCTTCAA
 AACATCATCTTTATCACCTGACACCAGGAGTTTTTCATTGAAAAGGATTTGAACCTGGTGTACTAACA
 TTTTAAAGACCACACAAGGAAGCAAAATCTTTCTGAAAGAAGTAAATGATACACTTCTGGTGAATGAATT
 GAAATCAAAAGAACTGACATCATGACAACAAATGGTGTAAATTCATGTTGTAGATAAACTCCTCTATCCA
 GCAGACACACCTGTTGGAATGATCAACTGCTGGAAATACTTAATAAATTAATCAAAATACATCCAAATTA
 AGTTTGTTCGTGGTAGCACCTTCAAAGAAATCCCCGTGACTGTCTATAGTCCTGAAATAAAATACACTAG
 GATTTCTACTGGAGGTGGAGAAACAGAAGAACTCTGAAGAAATGTTACAAGAAGAGGTCACCAAGGTC
 ACCAAATTCATTGAAGGTGGTGTGATGGTCATTTATTTGAAGATGAAGAAATTTAAAGACTGCTTCAGGGAG
 ACACACCGTGAGGAAGTTGCAAGCCAACAAAAAGTTCAAGGATCTAGAAGACGATTAAGGGAAGGTCG
 TTCTCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239477 representing NM_001286666
 Red=Cloning site Green=Tags(s)


MIPFLPMFSLLLLLIVNPINANNHYDKILAHSRIRGRDQGNVCAQQILGTTKKYFSTCKNWKYSICG
 QKTTVL YECCPGYMRMEGMKGPVLPIDHVVYGLGIVGATTTQRYSDASKLREEIEGKGSFTYFAPSNE
 AWDNLDSDIRRGLESNVNVELLNALHSHMINKRMLTKDLKNGMIIPSMYNNLGLFINHYPNGVVTVNCAR
 IIHGNQIATNGVVHVIDRVLTQIGTSIQDFIEAEDDLSSFRAAAITSDILEALGRDGHFTLFAPTNEAFE
 KLPRGVLERIMGDKVASEALMKYHILNLTQCSESIMGGAVFETLEGNTIEIGCDGDSITVNGIKMVNKKD
 IVTNGV IHLIDQVLIPDSAKQVIELAGKQQTTFDLVAQLGLASALRPDGEYTLAPVNNAFSDDTLSM
 DQRLLKILQNHILKVKVGLNELYNGQILETIGGKQLRVFVYRTAVCIENSCMEKGSQGRNGAIHFRE
 I IKPAEKSLHEK LQDKRFSTFLSLEADLKELLTQPGDWTLFVPTNDAFKGMTSEEKELIRDKNALQ
 NIILYHLTPGVFIGKGFEPGVTNLIKTTQGSKIFLKEVNDTLLVNELKSKESDIMTNGVIHVVDKLLYP
 ADTPVGNQDQLEILNLIKYIQIKFVRG STFKEIPVTVYSPEIKYTRISTGGGETEETLKKLLQEEVTKV
 TKFIEGGDGHLEDEEIKRLLQGDTPVRK LQANKKVQGSRRRLREGRSQ

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:


ACCN: NM_001286666

ORF Size: 2247 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286666.1](#), [NP_001273595.1](#)

RefSeq Size: 3129 bp

RefSeq ORF: 2250 bp

Locus ID: 10631

UniProt ID: [Q15063](#)

Cytogenetics: 13q13.3

Protein Families: Druggable Genome, Secreted Protein

MW: 84 kDa

Gene Summary: This gene encodes a secreted extracellular matrix protein that functions in tissue development and regeneration, including wound healing, and ventricular remodeling following myocardial infarction. The encoded protein binds to integrins to support adhesion and migration of epithelial cells. This protein plays a role in cancer stem cell maintenance and metastasis. Mice lacking this gene exhibit cardiac valve disease, and skeletal and dental defects. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]