

Product datasheet for **RC239515**

CSGLCAT (CHPF2) (NM_001284295) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CSGLCAT (CHPF2) (NM_001284295) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CSGLCAT
Synonyms:	chPF-2; ChSy-3; CSGLCA-T; CSGLcAT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC239515 representing NM_001284295
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGTCTCTGGCGCGCCTCCGCTCCCGCCGACTGGCCTGAGAACGAGCCTGAGCCTCCTGCGGGTTT
 CCTGGATCCAGGGGAGGGAGAAGATCCCTGTGTCGAGGCTGTAGGGGAGCGAGGAGGGCCACAGAATCC
 AGATTCCAGAGCTCGGCTAGACCAAAGTATGAAGACTTCAAACCCCGGATTGTCCCTACTACAGGGAC
 CCCAACAGCCCTACAAGAAGGTGCTCAGGACTCGGTACATCCAGACAGAGCTGGGCTCCCGTGAGCGGT
 TGCTGGTGGTGTCTGACCTCCCGAGCTACACTGTCCACTTTGGCCGTGGCTGTGAACCGTACGGTGGC
 CCATCACTTCCCTCGTTACTCTACTTCACTGGGCAGCGGGGGCCGGGCTCCAGCAGGGATGCAGGTG
 GTGTCTCATGGGATGAGCGGCCCGCTGGCTCATGTGAGAGACCCTGCGCCACCTTACACACACTTTG
 GGGCCGACTACGACTGGTCTTTCATCATGCAGGATGACACATATGTGCAGGCCCGCCGCTGGCAGCCCT
 TGCTGGCCACCTCAGCATCAACCAAGACCTGTACTTAGGCCGGGCAGAGGAGTTCATTGGCCGAGGGCAG
 CAGGCCCGGACTGTGTCATGGGGGCTTTGGCTACCTGTTGTGACGGAGTCTCTGCTTCGCTGCGGCCAC
 ATCTGGATGGCTGCCGAGGAGACATTCTCAGTGCCCGTCTGACGAGTGGCTTGGACGCTGCCTCATTGA
 CTCTCTGGGCGTCGGCTGTGTCTCACAGCACCAGGGGAGCAGTATCGCTCATTTGAACTGGCCAAAAAT
 AGGGACCTGAGAAGGAAGGGAGCTCGGCTTTCCTGAGTGCCTTCGCGCTGCACCCTGTCTCCGAAGGTA
 CCCTCATGTACCGGCTCCACAACCGCTTTCAGCGCTCTGGAGTTGGAGCGGGCTTACAGTGAATAGAACAA
 ACTGCAGGCTCAGATCCGAACCTGACCGTGTGACCCCGAAGGGGAGGCAGGGCTGAGCTGGCCCGTT
 GGGCTCCCTGCTCCTTTCACACCACACTCTCGCTTTGAGGTGCTGGGCTGGGACTACTTACAGAGCAGC
 ACACCTTCTCTGTGAGATGGGGCTCCAAGTGCCCACTACAGGGGGCTAGCAGGGCGGACGTGGGTGA
 TCGGTTGGAGACTGCCCTGGAGCAGCTCAATCGGCGCTATCAGCCCGCCTGCGCTTCCAGAAGCAGCGA
 CTGCTCAACGGCTATCGGCGCTTCGACCCAGCACGGGGCATGGAGTACACCCTGGACCTGCTGTTGGAAT
 GTGTGACACAGCGTGGGCACCGCGGGCCCTGGCTCGCAGGGTACGCTGCTGCGGCCACTGAGCCGGGT
 GGAATCTACCTATGCCCTATGTCACTGAGGCCACCGAGTGCAGCTGGTGTGCCACTCCTGGTGGCT
 GAAGCTGCTGCAGCCCGGCTTTCCTCGAGGCCTTTCAGCCAATGTCTGGAGCCACGAGAATGCAT
 TGCTCACCTGTTGCTGGTCTACGGGCCACGAGAAGGTGGCCGTGGAGCTCCAGACCCATTTCTTGGGGT
 GAAGGCTGCAGCAGCGGAGTTAGAGCGACGGTACCCTGGGACGAGGCTGGCTGGCTCGCTGTGCGAGCA
 GAGGCCCTTCCAGGTGCGACTCATGGACGTGGTCTCGAAGAAGCACCTGTGGACACTCTTCTTCC
 TTACCACCGTGTGACAAGGCCTGGGCCGAAGTCCCAACCGCTGTGCGATGAATGCCATCTCTGGCTG
 GCAGGCCTTCTTCCAGTCCATTTCCAGGAGTTCAATCCTGCCCTGTACCCACAGAGATCACCCCGAGG
 CCCCCGGGGGCTGGCCCTGACCCCCCTCCCTCCTGGTGTGACCCCTCCCGGGGGGCTCCTATAGGGG
 GGAGATTTGACCGGCAGGCTTCTGCGGAGGGCTGCTTCTACAACGCTGACTACCTGGCGGCCGAGCCCG
 GCTGGCAGGTGAAGTGGCAGGCCAGGAAGAGGAGGAAGCCCTGGAGGGGCTGGAGGTGATGGATGTTTTT
 CTCCGGTTCTCAGGGCTCCACCTCTTTCGGGCCGTAGAGCCAGGGCTGGTGCAGAAGTTCTCCCTGCGAG
 ACTGCAGCCACGGCTCAGTGAAGAAGTCTACCACCGCTGCCCGCTCAGCAACCTGGAGGGGCTAGGGG
 CCGTGCCAGCTGGCTATGGCTCTTTTGGAGCAGGAGCAGGCCAATAGCACT

ACGCGTACGCGGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001284295.2
RefSeq Size:	2730 bp
RefSeq ORF:	2295 bp
Locus ID:	54480
Cytogenetics:	7q36.1
Protein Pathways:	Chondroitin sulfate biosynthesis, Metabolic pathways
MW:	85.6 kDa
Gene Summary:	Transfers glucuronic acid (GlcUA) from UDP-GlcUA to N-acetylgalactosamine residues on the non-reducing end of the elongating chondroitin polymer. Has no N-acetylgalactosaminyltransferase activity.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC239515