

## Product datasheet for **RC239534**

### HA1 (ARHGAP45) (NM\_001282334) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	HA1 (ARHGAP45) (NM_001282334) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ARHGAP45
Synonyms:	HA-1; HLA-HA1; HMHA1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC239534 representing NM\_001282334  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGGTGGGGAGGAAAGGAGGAGCTGGGGAGACTGAGTCCCATCCGAGGATAGGGTTGGAAGTGGCCCT  
CCTGGCTCCACACCCACAGCAAGAGGGGAGTCCAACCTGCGCAAGGCCAAGCAGGGTTACGTGCAGCG  
CTGCGAGGACACGACAAGGCTCGTTCTCGTGGCCAAGGCGGAGGAGGAGCAGGCTGGCAGCGCGCCG  
GGAGCAGGCAGCACGGCCACCAAGACCTGGACAAGCGCGCGGCTGGAGGAGGAGGCCAAGAACAAGG  
CGGAGGAAGCTATGGCCACCTACCGCACCTGCGTGGCCGACGCGAAGACGAGCAGGAGCTGGAGGA  
TACCAAGGTGACGGCGCTGCGGCAGATCCAGGAGGTATCCGGCAGAGCGACCAACCATCAAGTCGGCC  
ACGATCTCTACTACCAGATGATGCATATGCAGACGGCGCCGCTGCCCGTGCCTCCAGATGCTGTGTG  
AGAGCAGCAAGCTGTATGACCCAGGCCAGCAGTACGCCTCCACGTGCGCCAGCTGCAGCGGGACCAGGA  
GCCCGATGTGACTACGACTTTGAGCCACAGTCTCCGCCAACGCCTGGTCCCCGTGCCTGCGTGGCCGG  
AAGAGCAGCTTCAACGTGAGTGTATGGCGCGGCCGGAGGCTGCCGGGAGCCCCCAGAAGAAGGCGGGT  
GCACTGAGGGCACACCTGCCAAGGACCACAGGGCCGGGCGAGGACACCAGGTTCAACAAGTCATGGCCGCT  
CTCGATCTCAGACTCGGACAGTGGGCTGGACCCCGGCCCTGGCGCAGGGGACTTTAAGAAGTTCGAGCGG  
ACGTATCCAGTGGTACCATGTCGTCCACGGAGGAGCTGGTGGACCCAGACGGTGGAGCCGGGGCTTCAG  
CCTTTGAGCAGGCTGACCTCAACGGCATGACCCCGAGCTGCCGGTGGCCGTGCCAGTGGACCGTTCCG  
CCACGAGGGGCTGTCCAAGGCGGCCGACTACCGGCTCCGGAAGCTCCGCACGCCGCCAAGTGGCCG  
GAGTGAACAGCTACGTCTACTCCAGGGTGTGAGTGTGAAGAGTGTGCCTGGCTGCCACAAGAAAT  
GTCTGGAGAGCTGGCCATACAGTGGCGGCACAAGAAGCTGCAAGGCCGCTGCAGCTGTTCCGGCAGGA  
CTTCAGCCACGCGGCCCGCAGCGCCCCGACGGCGTGCCTTCATCGTCAAGAAGTGCCTGCGAGATC  
GAGCGGCGGGCGCTGCGCACCAAGGGCATCTACCGGTTCAATGGGGTAAAGACACGCGTGGAGAAGCTGT  
GCCAGGCTTCGAGAACGGCAAGGAGCTGGTCGAGCTGTGCGAGGCTCGCCCCACGACATCAGCAACGT  
CCTCAAGCTCTACCTGCGTCAGCTTCCCGAGCCGCTCATCTCCTTCCGCTCTACCACGAGCTCGTAGGG  
CTGGCCAAGGACAGCCTGAAGGCAGAGGCCGAGGCCAAGGCGGCTCCCGGGCCGGCAGGACGGCTCGG  
AGAGCGAGGCAGTGGCGGTGGCCCTGGCAGGTGGCTGCGGGAGCTCCTGCGGGACCTGCCGCTGAGAA  
CCGGCCCTCGCTGCAGTACCTGCTGCGTCACTACGCAGGATCGTGGAGGTGGAGCAGGACAACAAGATG  
ACCCCCGGGAACCTGGGCATCGTGTTCGGGCCACGCTGCTTCGGCCACGGCCACCGAGGCCACCGTGT  
CCCTCTCCTCCCTGGTGGATTATCCCCATCAGGCCCGGCTCATCGAGACTCTCATCGTCCACTACGGCCT  
GGTCTTCGAGGAGGAGCCGGAGGAGACCCCGGGGCCAGGACGAGTCATCCAACCAGCGAGCTGAGGTA  
GTCGTCCAGGTGCCGTACCTGGAGGCGGGCGAGGCGGTGGTCTACCCGCTGCAGGAGGCGGCGGGGACG  
GGTGCAGAGAATCCCGAGTTGTGTCCAACGATTCGGACTCGGACCTAGAGGAGGCTCCGAGCTGCTGTC  
CTCATCGGAGGCCAGTGCCTGGGCCACCTCAGCTTCTGGAGCAGCAGCAGAGCGAGGCCAGCCTAGAG  
GTGGCTTCTGGCAGCCACAGCGGCAAGTGGAGCAGCTGGAGGCCACAGCCGGGAGGACGGGGACGGGG  
ACGAGGACGGCCCGCCAGCAGCTCTCAGGATTCAACACCAACCAGTCCAACAACGTGCTGCAGGCCCC  
ACTGCCCCCATGAGGCTCCGTGGCGGGCGGATGACACTGGGCTCCTGCAGGAAAGGCAGCCGGAAATTC  
GTG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC239534 representing NM\_001282334  
 Red=Cloning site Green=Tags(s)

```


MGVGRKGGAGETESHPRIGLELASWLPHPQQE AESNLRKAKQGYVQRCEHDHKARFLVAKAE EEQAGSAP
GAGSTATKTLDKRRRLEEEAKNAEEAMATYRTC VADAKTQKQLEEDTKVTALRQIQEVI RQSDQTIKSA
TISYYQMMHMQTAPLPVHFQMLCESSKLYDPGQ QYASHV RQLQRDQEPDVHYDFEPHVSANAWSPVMRAR
KSSFNVSDVARPEAAGSPPEEGGCTEGTPAKDHR AGRGHQVHKSWPLSISDSDSGLDPGP GAGDFKKFER
TSSSGTMSSTEELVDPDGAGASAFEQADLNGMTPE LPAVAVPSGPF RHEGLSKAARTHRLRKL RTPAKCR
ECNSYVYFQGAECCECLACHKKCLETLAIQC GHKKLQGRLLQFGQDFSHAARSAPDGV PFIKKCVCEI
ERRALRTKGIYRVNGVKTRVEKLCQAFENGKELV ELSQASPHDISNVLKLYLRQLPEPLISFRLYHEL VG
LAKDSLKAEAEAKAASRGRQDGSESEAVAVALAG RLRELLRDLPPENRASLQYLLRHLRRI VEVEQDNKM
TPGNLGI VFGPTLLRPRPTEATVSSLVDYPHQARV IETLIVHYGLVFEEPEETPGGQDESSNQRAEV
VVQVPYLEAGEAVVYPLQEAAADGCRESRVVSND SDSL EEASELLSSSEASALGHL SFLEQQQSEASLE
VASGSHSGSEEQLEATAREDDGDEDPAAQQLSGF NTNQSNV LQAPLPPMRLRGRMTLGSCRERQPEF
V
  
```

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:**


SgfI-RsrII

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

## Plasmid Map:



ACCN: NM\_001282334

ORF Size: 2313 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001282334.1](#), [NP\\_001269263.1](#)

RefSeq Size: 5286 bp

RefSeq ORF: 2316 bp

Locus ID: 23526

UniProt ID: [Q92619](#)

Cytogenetics: 19p13.3

**MW:** 85.1 kDa

**Gene Summary:** Contains a GTPase activator for the Rho-type GTPases (RhoGAP) domain that would be able to negatively regulate the actin cytoskeleton as well as cell spreading. However, also contains N-terminally a BAR-domain which is able to play an autoinhibitory effect on this RhoGAP activity. [UniProtKB/Swiss-Prot Function]