

Product datasheet for **RC239576**

UHRF1 (NM_001290050) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UHRF1 (NM_001290050) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UHRF1
Synonyms:	hNP95; hUHRF1; huNp95; ICBP90; Np95; RNF106; TDRD22
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC239576 representing NM_001290050
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGGATCCAGTTTCGACCATGGACGGGAGGCAGACCCACACGGTGGACTCGTCTCCAGGCTGACCA
AGGTGGAGGAGCTGAGGCGGAAGATCCAGGAGCTGTTCCACGTGGAGCCAGGCCTGCAGAGGCTGTTCTA
CAGGGGCAAACAGATGGAGGACGCCATACCCTCTTCGACTACGAGGTCCGCCTGAATGACACCATCCAG
CTCCTGGTCCGCCAGAGCCTCGTCTCCCCACAGCACCAAGGAGCGGGACTCCGAGCTCTCCGACACCG
ACTCCGGTCTGCTGGCCAGAGTGAGTGCAGACAAGTCTCCACCCACGGTGGAGCGGCCGCCGAGAC
TGACAGCAGGCCAGCCGATGAGGACATGTGGGATGAGACGGAATTGGGGCTGTACAAGGTCAATGAGTAC
GTCGATGCTCGGACACGAACATGGGGCGTGGTTTGGGCGCAGGTGGTCAAGGTGACCGGAAGGCC
CCTCCCGGACGAGCCCTGCAGTCCACGTCCAGGCCGGCTGGAGGAGACGTCATTTACCACGTGAA
ATACGACGACTACCCGGAGAACGGCGTGGTCCAGATGAACTCCAGGGACGTCGAGCGCGCGCCCGCACC
ATCATCAAGTGGCAGGACCTGGAGGTGGCCAGGTGGTCAATGCTCAACTACAACCCCGACAACCCCAAGG
AGCGGGGCTTCTGGTACGACCGGAGATCTCCAGGAAGCGGAGACCAGGACGGCGCGGGAACCTCTACGC
CAACGTGGTGTGGGGATGATTCTCTGAACGACTGTCCGATCATCTTCGTGGACGAAGTCTTCAAGATT
GAGCGGCCGGGTGAAGGGAGCCCATGGTTGACAACCCATGAGACGGAAGAGCGGGCCGTCTGCAAGC
ACTGCAAGGACGACGTGAACAGACTCTGCCGGTCTGCGCTGCCACCTGTGCGGGGCCGCGAGGACCC
CGACAAGCAGCTCATGTGCGATGAGTGCACATGGCCTCCACATCTACTGCCTGGACCCGCCCTCAGC
AGTGTTCACAGCAGGACGAGTGGTACTGCCCTGAGTCCGGAATGATGCCAGCGAGGTGGTACTGGCGG
GAGAGCGGCTGAGAGAGAGCAAGAAGAAGGCGAAGATGGCTCGGCCACATCGTCTCACAGCGGGACTG
GGCAAGGGCATGGCCTGTGTGGCCGCAACCAAGGAATGTACCATCGTCCCGTCCAACCACTACGGACCC
ATCCCGGGATCCCGTGGGCACCATGTGGCGTTCCGAGTCCAGGTGAGCGAGTGGGTGTCCATCGGC
CCCACGTGGTGGCATAACGCGCCGAGCAACGACGGAGCGTACTCCCTAGTCTGGCGGGGGCTATGA
GGATGATGTGGACCATGGGAATTTTTTACATACACGGGTAGTGGTGGTGGAGATCTTCCGGCAACAAG
AGGACCGCGAACAGTCTTGTGATCAGAACTCACCAACCAACAGGGCGTGGCTCTCAACTGCTTTG
CTCCCATCAATGACCAAGAAGGGCCGAGGCCAAGGACTGGCGGTGGGGAAGCCGGTCAAGGTGGTGGC
CAATGTCAAGGGTGGCAAGAATAGCAAGTACGCCCGCTGAGGGCAACCGCTACGATGGCATCTACAAG
GTTGTGAAATACTGGCCCGAGAAGGGGAAGTCCGGTTTTCTCGTGTGGCGCTACCTTCTGCGGAGGGACG
ATGATGAGCCTGGCCCTTGGACGAAGGAGGGGAAGGACCGGATCAAGAAGCTGGGGCTGACCATGCAGTA
TCCAGAAGGCTACCTGGAAGCCCTGGCCAACCGAGAGCGAGAGAAGGAGAACAGCAAGAGGGAGGAGGAG
GAGCAGCAGGAGGGGGGCTTCGCGTCCCCAGGACGGGCAAGGGCAAGTGGAAAGCGGAAGTGGCAGGAG
GTGGCCCGAGCAGGGCCGGTCCCCGCGCCGACATCCAAGAAAACCAAGGTGGAGCCCTACAGTCTCAC
GGCCAGCAGAGCAGCCTCATCAGAGAGGACAAGAGCAACGCCAAGCTGTGGAATGAGGTCTGGCGTCA
CTAAGGACCGCCGCGGAGCGGACGCCGTTCCAGTTGTTCTGAGTAAAGTGGAGGAGACGTTCCAGT
GTATCTGCTGTGAGGAGTGGTGTCCGGCCATCACGACCGTGTGCCAGCACACGTGTGCAAGGACTG
CCTGGACAGATCCTTTCCGGCACAGGTGTTTCAGCTGCCCTGCCCTGCGCTACGACCTGGGCCGAGCTAT
GCCATGCAGGTGAACAGCCTCTGCAGACCGTCTCAACCAGCTTTCGCCGCTACGGCAATGGCCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC239576 representing NM_001290050
 Red=Cloning site Green=Tags(s)


MWIVRTMDGRQTHTVDSL SRL TKVEELRRKIQELFHVEPGLQRLFYRGKQMEDGHTLFDYEVRLNDTIQ
 LLVLRQSLVLPSTKERDSELSDTDSGCCLGQSESDKSSSTHGEAAAETDSRPAEDMWDTELGLYKVNEY
 VDARDTNMGAWFEAQVVRVTRKAPSRDEPCSSSRPALEEDVIYHVKYDDYPENGVVQMNSRDVRAART
 IIKWQDLEVGQVVMLNYPDNPKERGFWDYDAEISRKRETRTARELYANVVLGDDSLNDCRIIFVDEVFKI
 ERPGEGSPMVDNPMRRKSGPSCCHKKDDVNRLCRVCACHLCGGRQDPDKQLMCDECDMAFHICYLDPPLS
 SVPSEDEWYCPECRNDASEVVLAGERLRESKKKAKMASATSSSQRDWKGCMACVGRKTECTIVPSNHYGP
 IPGIPVGTMRFRVQVSESGVHRPHVAGIHGRSNDGAYSLVLAGGYEDDVHDGNFFTYTGSGRDLSGNK
 RTAEQSCDQKL TNTNRALALNCFAPINDQEGAEAKDWRSGKPVVRVVRNVKGGKNSKYAPAEGNRYDGIYK
 VVKYWPEKKGSGFLVWRVLLRRDDDEPGPWTKEGKDRIKLLGLTMQYPEGYLEALANREREKENSREEE
 EQQEGGFASPRTGKGKWRKSAGGGPSRAGSPRRTSKTKVEPYSHTAQSSLIREDKSNAKLWNEVLAS
 LKDRPASGSPFLFLSKVEETFQICCCQELVFRPITTVQCQHNVCCKDLDRSFRAQVFSACPACRYDLGRSY
 AMQVNVQLQTVLNQLFPGYGNR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001290050

ORF Size: 2379 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290050.1](#), [NP_001276979.1](#)

RefSeq Size: 4339 bp

RefSeq ORF: 2382 bp

Locus ID: 29128

UniProt ID: [Q96T88](#)

Cytogenetics: 19p13.3

Protein Families: Druggable Genome, Transcription Factors

MW: 90.3 kDa

Gene Summary: This gene encodes a member of a subfamily of RING-finger type E3 ubiquitin ligases. The protein binds to specific DNA sequences, and recruits a histone deacetylase to regulate gene expression. Its expression peaks at late G1 phase and continues during G2 and M phases of the cell cycle. It plays a major role in the G1/S transition by regulating topoisomerase IIalpha and retinoblastoma gene expression, and functions in the p53-dependent DNA damage checkpoint. It is regarded as a hub protein for the integration of epigenetic information. This gene is up-regulated in various cancers, and it is therefore considered to be a therapeutic target. Multiple transcript variants encoding different isoforms have been found for this gene. A related pseudogene exists on chromosome 12. [provided by RefSeq, Feb 2014]