

Product datasheet for **RC239582**

HERC4 (NM_001278186) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HERC4 (NM_001278186) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	HERC4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239582 representing NM_001278186
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGATGAAAATGGAAGGTGGAGTGTACTTTGGAGCTGGAGGGTATGGTCAGTTGGCCATAATTCTA
 CCAGTCATGAAATAAACCAAGGAAAGTTTTGAACTTATGGGAAGCATTGTCACCTGAGATTGCTTGTGG
 ACGGCAGCACACTTCTGCTTTTGTCTTCATCAGGACGAATTTACTCTTTGGGCTTGGTGAATGGG
 CAGCTGGGAACCGTTCAACAAGCAACAGGAAAAGCCCTTTACTGTAAGGAAATTTGGTACCCCTATA
 ATGGGCAGTGTCTACCAGATATTGATTCTGAAGAATATTTCTGTGTAAAAAGAAATTTCTCAGGGGGAGA
 TCAAAGCTTTTACATTACTCTAGTCCCCAGAAGTGGGCCACCCAGATGACTTCAGATGTCCCAATCCG
 ACAAGCAGATCTGGACAGTGAATGAAGCTCTAATTCAGAAATGGCTGAGCTATCCTTCTGGAAGGTTTC
 CTGTGGAGATAGCCAATGAGATAGATGGAACGTTTTCTCTCTGGTGCCTAAATGGAAGTTTTTAGC
 TGTTAGCAATGATGATCACTATAGAACAGTACCAGATTTTCAGGGGTTGATATGAATGCTGCTAGGCTT
 TTATCCACAAACTTATACAACCTGATCATCCGCAGATATCTCAGCAGGTGGCAGCTAGTTTGGAAAAGA
 ATCTTATTCCTAAACTGACTAGCTCCTTACCTGATGTTGAAGCATTGAGGTTTTATCTTACTCTACCAGA
 ATGTCCCCTGATGAGTGATTCCAACAATTTCAACAATAGCAATTCCTTTGGTACAGCTCTTGTGAAC
 CTAGAAAAGGCCACTGAAAGTACTTGAAAAGTGGTGGTCACTTGAACCTCCACTATTCTCAAGA
 TAGTAGAACTTTTTAAGGAAGTTGGTACATCTTTGAAACTCTACAAGATCGGATTCCCCCTCTGA
 AAGAAGATTTTCAACAGTTTTCTCATACTGCATTAAGGTTTTAGAAACTACATAGGGTAAATGAG
 AAAATGGGACAGATTATACAGTATGATAAATTTATATACATGAAGTACAAGAATTGATAGACATAAGAA
 ATGATTATCAACTGGGTCCAACAGCAGCCTATGGAATGTTGGCAGATATCCCTGTTACAATCTGTAC
 ATATCCATTTGTATTTGATGCCCAAGCAAAAAGTACTCTGTTACAGACCGATGCAGTCTTACAGATGCAG
 ATGGCTATTGATCAGGCCACAGGCAGAATGTCTCCTCTTTTTCTCCAGTGATTGAATCTGTGAATC
 CCTGCTTAATTCTAGTGGTGCAGAGAAAAATTTGTAGGAGATGCAATGGAAGTCTTAGGAAAACAAA
 GAACATAGATTACAAGAAGCCACTCAAGGTTATATTTGTTGGAGAAGATGCTGTGGATGCAGGAGGGGTG
 CGCAAAGAATTTTCTTGCTCATCATGAGGGAATTATTGGATCCTAAATACGGCATGTTTAGGTATTATG
 AAGATTCCAGGCTCATTGGTTTTCTGATAAGACATTTGAAGACAGTGATTTGTTCCATTTGATTGGTGT
 TATCTGTGGCTTAGCAATTTATAATTGTACCATTGTGGACCTCCATTTCTTTGGCTTTATATAAGAAA
 CTACTGAAAAGAAGCCATCCTTGGATGATTTGAAAGAATAATGCCTGATGTTGGGAGAAGCATGCAAC
 AGTTACTGGATTATCCAGAAGATGACATAGAGGAAACATTTTGTCTTAATTTTACGATCACAGTTGAAAA
 CTTTGGTGCAACAGAAGTGAAGAGCTGGTTCTAAATGGTGCAGACACAGCTGTTAACAAACAAAATCGG
 CAAGAGTTTGTGATGCTTATGTGGATTACATATCAATAAATCAGTGGCTTCTTATTTGATGCTTTTC
 ATGCGGGCTTTCATAAGGTCTGTGGAGGAAAAGTCTTCTGCTCTTTCAGCCTAATGAACTACAAGCAAT
 GGTCAATTGAAAATACAAATTATGATTGGAAGGAACTGGAAAAGAATACAGAATACAAAGGGGAATATTGG
 GCAGAATCCTACGATAAAAAATTTTTGGGAAGTATTTACGAATTACCATTGGAAAAGAAGAAACAGT
 TTCTGTTATTTTACAGGTAGTATCGCATTCTTCTTGGTATGAAGAGTCTGAAACTAGTCATCCA
 GTCCACAGGAGGTGGTGGAGGATATCTCCAGTTTCCCACTTGTGTTTAACTTCTGGATCTTCCAAAA
 TATACAGAAAAAGAACTCTACGCTCTAACTGATCCAAGCTATTGATCACAATGAAGGCTTCAGTTTAA
 TA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC239582 representing NM_001278186
Red=Cloning site Green=Tags(s)

```
MMKMEGGVFTFGAGGYGQLGHNSTSEINPRKVFELMGSIVTEIACGRQHTSAFVPSGRIYSFGLGGNG
QLGTGSTSNRKSPTVKGNWYPYNGQCLPDIDSEEYFCVKRIFSGGDQSFSHYSSPQNCGPPDDFRCPNP
TKQIWTVNEALIQKWLSPSGRFPVEIANEIDGTFSSSGCLNGSFLAVSNDHDYRTGTRFSGVDMNAARL
LFHKLIQPDHPQISQVVAASLEKNLIPKLTSSLPDVEALRFYLTLPCEPLMSDSNNFTTIAIPFGTALVN
LEKAPLKVLENWWSVLEPPLFLKIVELFKEVVVHLLKLYKIGIPPSERRIFNSFLHTALKVLEILHRVNE
KMGQIIQYDKFYIHEVQELIDIRNDYINWVQQQAYGMLADIPVTICTYPFVFDQAQAKTLLQTDVAVLQMQ
MAIDQAHQRQNVSSLFLPVIESVNPCLILVVRRENIVGDAMEVLRKTKNIDYKPLKIVIFVGEDAVIDAGGV
RKEFFLLIMRELLDPKYGMFRYYEDSRLIWFSDKTFEDSDLFHLIGVICGLAIYNCTIVDLHFPLALYKK
LLKKKPSLDDLKELMPDVGRSMQQLLDYPEDDIEETFCLNFTITVENFGATEVKELVLNGADTAVNKQNR
QEFVDAYVDYIFNKSVALSFDFAFHAGFHKVCGGKVLVLLFQPNELQAMVIGNTNYDWKELEKNTEYKGEYW
AEHPTIKIFWEVFHELPLEKKKQFLFLTGSDRIPILGMKSLKLVISTGGGEEYLPVSHTCFNLLDLPK
YTEKETLRSKLIQAIDHNEGFSLI
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001278186

ORF Size: 2382 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278186.2](#)

RefSeq Size: 4255 bp

RefSeq ORF: 2385 bp

Locus ID: 26091

UniProt ID: [Q5GLZ8](#)

Cytogenetics: 10q21.3

Protein Families: Druggable Genome

Protein Pathways: Ubiquitin mediated proteolysis

MW: 90.6 kDa

Gene Summary: HERC4 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).[supplied by OMIM, Mar 2008]