

## Product datasheet for **RC239706**

### **Dynamin 1 (DNM1) (NM\_001288739) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Dynamin 1 (DNM1) (NM_001288739) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DNM1
Synonyms:	DEE31; DNM; EIEE31
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide  
Sequence:**


>RC239706 representing NM\_001288739  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGCAACCGCGCATGGAAGATCTCATCCCGCTGGTCAACCGGCTGCAAGACGCCTTCTCTGCCATCG  
 GCCAGAACCGCGGACCTCGACCTGCCGACAGATCGCTGTGGTGGCGGCCAGAGCGCCGCAAGAGCTCGGT  
 GCTCGAGAATTTCTAGGCAGGGACTTCTTGCCTCGAGGATCTGGCATTGTCACCCGACGTCCTCCCTGGTC  
 TTGCAGCTGGTCAATGCAACCACAGAATATGCCGAGTTCCTGCACTGCAAGGAAAGAAATTCACCGACT  
 TCGAGGAGGTGCGCCTTGGATCGAGGCCGAGACCGAGGGTCAACCGCACCAACAAGGCATCTCGCC  
 GGTGCCTATCAACCTCCGCTCTACTCGCCGACGCTGCTGAACCTGACCTGGTGGACCTGCCCGGAATG  
 ACCAAGTCCCGTGGGGACCAACCTCCCGACATCGAGTTCAGATCCGAGACATGCTTATGCAGTTTG  
 TCACCAAGGAGAACTGCCTCATCTGGCCGTGTCCTCCCGCAACTCTGACCTGGCCAATTCTGACGCCCT  
 CAAGTTCGCCAAGGAGGTGGACCCCGAGGCCAGCGCACCATCGGGTTCATCACAAGCTGGACCTGATG  
 GACGAGGGCACAGATGCCCGTGTGTGCTGGAGAACAAGCTGCTCCCGTGCAGAGGGTACATTTGGAG  
 TGGTGAACCGGAGCCAGAAGGACATTGATGGCAAGAAGGACATTACCGCCGCTTGGCTGCTGAACGAAA  
 GTTCTTCTCTCCATCCATCTTATCGCCACTTGGCTGACCGTATGGGCACGCCCTACCTGCAGAAGGTC  
 CTCAATCAGCAACTGACGAACCACATCCGGGACACACTGCCGGGCTGCGGAACAAGCTGCAGAGCCAGC  
 TACTGTCCATTGAGAAGGAGGTGGAGGAATACAAGAACTCCGCCCTGATGACCCAGCTCGCAAGACCA  
 GGCCCTGCTGCAGATGGTCCAGCAGTTCGCCGTAGACTTTGAGAAGCGCATTGAGGGCTCAGGAGATCAG  
 ATCGACACCTACGAAGTGTGAGGGGAGCCCGCATTAAACCGAATCTTCCACGAGCGCTTCCCTTCGAGC  
 TGGTCAAGATGGAGTTTGTGAGAAGAACTCCGAAGGGAGATCAGCTATGCTATCAAGAATATCCATGG  
 CATTAGGACGGGCTCTTACACCTGACCTCGCTTTTGAAGCCACAGTGAAAAAGCAGGTGCAGAAGCTC  
 AAAGAGCCAGTATCAAGTGTGTGGATATGGTAGTCAGTGAAGTCAAGCCACCATCAGAAAAGTGTAGCG  
 AAAAGTCCAGCAGTACCGCGGCTACGGGAGGAGATGGAGCGCATCGTGACCACCCACATCCGGGAGCG  
 CGAGGGCCGCACTAAGGAGCAGGTATGCTTCTCATCGATATCGAGCTGGCTTACATGAACCAACCAT  
 GAGGACTTCATAGGCTTTGCCAATGCTCAGCAGAGGAGCAACAGATGAACAAGAAGAAGACTTCAGGGA  
 ACCAGGATGAGATTCTGGTATCCGCAAGGGCTGGCTGACTATCAATAATATTGGCATCATGAAAGGGGG  
 CTCCAAGGAGTACTGGTTTGTGCTGACTGCTGAGAATCTGTCTGGTACAAGGATGATGAGGAGAAAGAG  
 AAGAAATACATGTGTCTGTGGACAACCTCAAGCTGCGGGACGTGGAGAAGGGCTTTATGTGAGCAAGC  
 ATATCTTTGCCCTCTTTAACACGGAGCAGAGGAATGTCTACAAGGATTATCGGCAGCTGGAGCTAGCCTG  
 TGAGACACAGGAGGAGGTGGACAGCTGGAAGGCCTCCTTCTGAGGGCTGGCGTGTACCCTGAGCGTGT  
 GGGGACAAAGAGAAAGCCAGCGAGACCGAGGAGAATGGCTCCGACAGCTTCATGCATTCCATGGACCCAC  
 AGCTGGAACGGCAAGTGGAGACCATCCGGAATCTTGTGGACTCATACATGGCCATTGTCAACAAGACCGT  
 GAGGGACCTCATGCCAAGACCATATGCACCTCATGATTAACAATACCAAGGAGTTCATCTTCTCGGAG  
 CTGCTGGCCAACCTGTACTCGTGTGGGACCAGAACACGCTGATGGAGGAGTCCGGGAGCAGGCACAGC  
 GGCAGCAGGATGCTGCGCATGTACCACGCACTGAAGGAGGCGCTCAGCATCATCGCGACATCAACAC  
 GACCACCGTCAGCAGCCCATGCCCGCCCGTGGACGACTCCTGGCTGCAGGTGCAGAGGCTACCGGCC  
 GGACGCAGGTCCGCCACGTCCAGCCCCACGCCGAGCGCCGAGCCCCGCGTGCCTCCCGAGCCCCCG  
 GGTCCGCGGGCCCTGCTCCTGGGCTCCGCTGCTGGTCCGCCCTGGGGGGGGCCCCCGTGCCTC  
 CAGGCCGGGGCTTCCCCTGACCCTTTCGGCCCTCCCCCTCAGGTGCCCTCGCGCCCAACCGCGCCCCG  
 CCCGGGTCCCCAGCCGATCGGGTCAAGCAAGTCCATCCGTCCTGAGAGCCCCAGGCCCTTCGACC  
 TC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



**Plasmid Map:**


**ACCN:** NM\_001288739

**ORF Size:** 2592 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001288739.2](#)

**RefSeq Size:** 3260 bp

**RefSeq ORF:** 2595 bp

**Locus ID:** 1759

**UniProt ID:** [Q05193](#)

**Cytogenetics:** 9q34.11

<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Endocytosis, Fc gamma R-mediated phagocytosis
<b>MW:</b>	97.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the dynamamin subfamily of GTP-binding proteins. The encoded protein possesses unique mechanochemical properties used to tubulate and sever membranes, and is involved in clathrin-mediated endocytosis and other vesicular trafficking processes. Actin and other cytoskeletal proteins act as binding partners for the encoded protein, which can also self-assemble leading to stimulation of GTPase activity. More than sixty highly conserved copies of the 3' region of this gene are found elsewhere in the genome, particularly on chromosomes Y and 15. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]</p>