

Product datasheet for **RC400039**

p53 (TP53) (NM_000546) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	p53 (TP53) (NM_000546) Human Mutant ORF Clone
Mutation Description:	V173M
Affected Codon#:	173
Affected NT#:	c.517
Nucleotide Mutation:	TP53 mutant (V173M), Myc-DDK-tagged ORF clone of Homo sapiens tumor protein p53 (TP53), transcript variant 1 as transfection-ready DNA
Effect:	Missense
Symbol:	p53
Synonyms:	BCC7; BMFS5; LFS1; P53; TRP53
E. coli Selection:	Ampicillin (100 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-Myc-DDK (PS100007)
Tag:	Myc-DDK
ACCN:	NM_000546
ORF Size:	1183 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC400039 representing NM_000546
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGGAGCCGAGTCAGATCCTAGCGTCGAGCCCCCTCTGAGTCAGGAAACATTTTCAGACCTATGGA
 AACTACTTCCTGAAAACAACGTTCTGTCCCTTGCCGTCCTCAAGCAATGGATGATTTGATGCTGTCCCC
 GGACGATATTGAACAATGGTTCCTGAAGACCCAGGTCCAGATGAAGCTCCAGAAATGCCAGAGGCTGT
 CCCCCGTGGCCCTGCACCAGCAGCTCTACCCGGCGGCCCTGCACCAGCCCCCTCCTGGCCCTGT
 CATCTTCTGTCCCTCCAGAAAACCTACCAGGGCAGCTACGGTTTCCGTCTGGGCTTCTTGCACTTCTGG
 GACAGCCAAGTCTGTGACTTGCACGTAATCCCCTGCCCTCAACAAGATGTTTTGCCAACTGGCCAAGACC
 TGCCCTGTGCAGCTGTGGTTGATTCCACACCCCGCCCGCACCCGCTCCGCGCCATGGCCATCTACA
 AGCAGTCACAGCACATGACGGAGGTTATGAGGCGCTGCCCCACCATGAGCGCTGCTCAGATAGCGATGG
 TCTGGCCCTCCTCAGCATCTTATCCGAGTGAAAGGAAATTTGCGTGTGGAGTATTTGGATGACAGAAA
 CACTTTTCGACATAGTGTGGTGGTGCCTATGAGCCGCTGAGGTTGGCTCTGACTGTACCACCATCCAC
 TACAACACTACATGTGTAACAGTTCCTGCATGGGCGGCATGAACCGGAGGCCATCCTCACCATCATCACAC
 TGGAAAGACTCCAGTGGTAATCTACTGGGACGGAACAGCTTTGAGGTGCGTGTGTTGTGCCTGTCTGGGAG
 AGACCGGCGCACAGAGGAAGAGAATCTCCGCAAGAAAGGGGAGCCTCACCACGAGCTGCCCCAGGGAGC
 ACTAAGCGAGCACTGCCCAACAACAGCTCCTCTCCCAGCCAAGAAGAAACCCTGGATGGAGAAT
 ATTTACCCTTCAGATCCGTGGGCGTGAGCGCTTCGAGATGTTCCGAGAGCTGAATGAGGCTTGGAACT
 CAAGGATGCCAGGCTGGGAAGGAGCCAGGGGGAGCAGGGCTCACTCCAGCCACCTGAAGTCCAAAAAG
 GTTCAGTCTACCTCCCGCCATAAAAACTCATGTTCAAGACAGAAGGGCCTGACTCAGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC400039 representing NM_000546
 Red=Cloning site Green=Tags(s)

MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIEQWFTEDPGPDEAPRMPEAA
 PPVAPAAPAPTPAAPAPAPSWPLSSSVPSQKTYQGSYGFRLGFLHSGTAKSVTCTYSPALNKMFCQLAKT
 CPVQLWVDSTPPPGRVTRAMAIYKQSQHMTEVMRRCPHHERCSDSDGLAPPQHLIRVERKFACGVFG*QK
 HFST*CGGAL*AA*GWL*LYHHPLQLHV*QFLHGRHEPEAHPHHHHTGRLQW*STGTEQL*GACCLCSWE
 RPAHRGRESPQERGASPRAPREH*ASTAQHQQLLSPAKEETGWRIHFHSPDPA*ALRDVPAE*GLGT
 QGCPGWEGARGEQGLQPPEVQKGSVYLPP*KTHVQDRRA*LR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

/chromatograms/ja1116_h12.zip

Restriction Sites:

Sgfl-Mlul

Protein Pathways:	Amyotrophic lateral sclerosis (ALS), Apoptosis, Basal cell carcinoma, Bladder cancer, Cell cycle, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, Glioma, Huntington's disease, MAPK signaling pathway, Melanoma, Neurotrophin signaling pathway, Non-small cell lung cancer, p53 signaling pathway, Pancreatic cancer, Pathways in cancer, Prostate cancer, Small cell lung cancer, Thyroid cancer, Wnt signaling pathway
MW:	43.6 kDa
Gene Summary:	This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons from identical transcript variants (PMIDs: 12032546, 20937277). [provided by RefSeq, Dec 2016]