

Product datasheet for **RC400068**

p53 (TP53) (NM_000546) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	p53 (TP53) (NM_000546) Human Mutant ORF Clone
Mutation Description:	L257P
Affected Codon#:	257
Affected NT#:	c.770
Nucleotide Mutation:	TP53 mutant (L257P), Myc-DDK-tagged ORF clone of Homo sapiens tumor protein p53 (TP53), transcript variant 1 as transfection-ready DNA
Effect:	Missense
Symbol:	p53
Synonyms:	BCC7; BMFS5; LFS1; P53; TRP53
E. coli Selection:	Ampicillin (100 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-Myc-DDK (PS100007)
Tag:	Myc-DDK
ACCN:	NM_000546
ORF Size:	1182 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC400068 representing NM_000546
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGGAGCCGAGTCAGATCCTAGCGTCGAGCCCCCTCTGAGTCAGGAAACATTTTCAGACCTATGGA
 AACTACTTCTGAAAACAACGTTCTGTCCCTTGCCGTCCCAAGCAATGGATGATTTGATGCTGTCCCC
 GGACGATATTGAACAATGGTTCCTGAAGACCCAGGTCCAGATGAAGCTCCAGAATGCCAGAGGCTGT
 CCCCCGTGGCCCTGCACCAGCAGCTCTACCCGGCGGCCCTGCACCAGCCCCCTCCTGGCCCTGT
 CATCTTCTGTCCCTCCAGAAAACCTACCAGGGCAGCTACGGTTTCCGTCTGGGCTTCTTGCACTTCTGG
 GACAGCCAAGTCTGTGACTTGCACGTACTCCCCTGCCCTCAACAAGATGTTTTGCCAACTGGCCAAGACC
 TGCCCTGTGCAGCTGTGGTTGATTCCACACCCCGCCCGCACCCGCTCCGCGCCATGGCCATCTACA
 AGCAGTCACAGCACATGACGGAGGTTGTGAGGCGCTGCCCCACCATGAGCGCTGCTCAGATAGCGATGG
 TCTGGCCCTCCTCAGCATCTTATCCGAGTGGAAGGAAATTTGCGTGTGGAGATTTGGATGACAGAAAC
 ACTTTTCGACATAGTGTGGTGGTCCCTATGAGCCGCTGAGGTTGGCTCTGACTGTACCACCATCCACT
 ACAACTACATGTGTAACAGTTCCTGCATGGCGGCATGAACCGGAGGCCATCCTCACCATCATCACACC
 GGAAGACTCCAGTGGTAATCTACTGGGACGGAACAGCTTTGAGGTGCGTGTTTGTGCCTGTCTGGGAGA
 GACCGGCGCACAGAGGAAGAAATCTCCGCAAGAAAGGGGAGCCTCACCAGAGCTGCCCCAGGGAGCA
 CTAAGCGAGCACTGCCAAACAACACCAGCTCCTCTCCCCAGCCAAAGAAGAAACCACTGGATGGAGAATA
 TTTACCCCTTCAGATCCGTGGGCGTGAGCGCTTCGAGATGTTCCGAGAGCTGAATGAGGCCTTGGAACTC
 AAGGATGCCAGGCTGGGAAGGAGCCAGGGGGAGCAGGGCTCACTCCAGCCACCTGAAGTCCAAAAAGG
 GTCAGTCTACCTCCCGCCATAAAAACTCATGTTCAAGACAGAAGGGCCTGACTCAGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC400068 representing NM_000546
 Red=Cloning site Green=Tags(s)

MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIEQWFTEDPGPDEAPRMPEAA
 PPVAPAPAAPTPAAPAPAPSWPLSSSVPSQKTYQGSYGFRLGFLHSGTAKSVTCTYSPALNKMFCQLAKT
 CPVQLWVDSTPPPGRVTRAMAIYKQSQHMTEVVRRCPPHHERCSDSDGLAPPQHLIRVEGNLRVEYLDDRN
 TFRHSVVVPYEPPEVGSDCCTIIHYNMCSNCSMGMNRRPILTIITPEDSSGNLLGRNSFEVRVCACPGR
 DRRTEENLRKKGEPHELPPGSTKRALPNNTSSSPQPKKPLDGEYFTLQIRGRERFEMFRELNEALEL
 KDAQAGKEPGGSAHSSHLKSKKGQTSRHKKLMFKTEGPDS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

/chromatograms/ja1119_e12.zip

Restriction Sites:

Sgfl-Mlul

Protein Pathways: Amyotrophic lateral sclerosis (ALS), Apoptosis, Basal cell carcinoma, Bladder cancer, Cell cycle, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, Glioma, Huntington's disease, MAPK signaling pathway, Melanoma, Neurotrophin signaling pathway, Non-small cell lung cancer, p53 signaling pathway, Pancreatic cancer, Pathways in cancer, Prostate cancer, Small cell lung cancer, Thyroid cancer, Wnt signaling pathway

MW: 43.6 kDa

Gene Summary: This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons from identical transcript variants (PMIDs: 12032546, 20937277). [provided by RefSeq, Dec 2016]