

Product datasheet for **RC400171**

TPOR (MPL) (NM_005373) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	TPOR (MPL) (NM_005373) Human Mutant ORF Clone
Mutation Description:	S505N
Affected Codon#:	505
Affected NT#:	c.1514
Nucleotide Mutation:	MPL Mutant (S505N), Myc-DDK-tagged ORF clone of Homo sapiens myeloproliferative leukemia virus oncogene (MPL) as transfection-ready DNA
Effect:	Missense
Symbol:	MPL
Synonyms:	C-MPL; CD110; MPLV; THCYT2; THPOR; TPOR
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_005373
ORF Size:	1905 bp
Restriction Sites:	Sgfi-MluI



[View online »](#)

ORF Nucleotide
Sequence:

>RC400171 representing NM_005373
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCCCTCTGGGCCCTTTCATGGTCACCTCCTGCCTCCTCGCCCTCAAACCTGGCCCAAGTCA
GCAGCCAAGATGTCTCCTTGCTGGCATCAGACTCAGAGCCCTGAAGTGTTCCTCCGAACATTTGAGGA
CCTCACTTGCTTCTGGGATGAGGAAGAGGCAGGCCCCAGTGGGACATACCAGCTGCTGTATGCCTACCCG
CGGGAGAAGCCCCGTGCTTGCCCTGAGTTCCAGAGCATGCCCACTTTGGAACCCGATACGTGTGCC
AGTTTCCAGACCAGGAGGAAGTGCCTCTCTTTCCGCTGCACCTCTGGGTGAAGAATGTGTTCTTAA
CCAGACTCGGACTCAGCGAGTCTCTTTGTGGACAGTGTAGGCCTGCCGGCTCCCCCAGTATCATCAAG
GCCATGGTGGGAGCCAGCCAGGGAACTTCAGATCAGCTGGGAGGAGCCAGTCCAGAAATCAGTGATT
TCCTGAGGTACGAACTCCGCTATGGCCCAAGATCCCAAGAACTCCACTGGTCCACGGTCATACAGCT
GATTGCCACAGAAACCTGCTGCCCTGCTCTGCAGAGCCCTCACTCAGCCTCTGCTCTGGACCAGTCTCCA
TGTGCTCAGCCACAATGCCCTGGCAAGATGGACCAAAAGCAGACCTCCCAAGTAGAGAACTTCAAGTCT
TGACAGCAGAGGGTGAAGTGCCTCATCTCAGGACTCCAGCCTGGCAACTCCTACTGGCTGCAGCTGCG
CAGCGAACCTGATGGGATCTCCCTCGGTGGTCTCTGGGATCCTGGTCCCTCCCTGTGACTGTGGACCTG
CCTGGAGATGCAGTGGCACTTGGACTGCAATGCTTTACCTTGGACCTGAAGAATGTTACCTGTCAATGGC
AGCAACAGGACCATGCTAGTCCCAAGGCTTCTTACACAGCAGGGCAGCGGTGCTGCCCCAGAGACAG
GTACCCCATCTGGGAGAACTGCGAAGAGGAAGAGAAAACAAATCCAGGACTACAGACCCACAGTTCTCT
CGTGCCACTTCAAGTCAGAAATGACAGCATTATTACATCCTTGTGGAGGTGACCACAGCCCCGGGTA
CTGTTACAGCTACCTGGGCTCCCTTTCTGGATCCACCAGGCTGTGCGCTCCCAACCCAACTTGA
CTGGAGGAGATCTCCAGTGGGATCTGGAATTGGAGTGGCAGCACCCATCGTCTGGGAGCCCAAGAG
ACCTGTTATCAACTCCGATACACAGGAGAAGGCCATCAGGACTGGAAGGTGCTGGAGCCGCTCTCGGGG
CCCAGGAGGGACCCTGGAGCTGCGCCCGCATCTCGCTACCGTTTACAGCTGCGCGCCAGGCTCAACGG
CCCCACCTACCAAGGTCCTGGAGCTCGTGGTCCGACCCAACTAGGGTGGAGACCCGACCGAGACCGCC
TGGATCTCCTTGGTGACCCTCTGCATCTAGTGTGGGCTCAACGCCGTCTGGGCTGCTGCTGCTGA
GGTGGCAGTTTCTGCACACTACAGGAGACTGAGGCATGCCCTGTGGCCCTCACTTCCAGACCTGCACCG
GGTCTAGGCCAGTACCTTAGGGACTGCAGCCCTGAGCCCGCCCAAGGCCACAGTCTCAGATACTGT
GAAGAAGTGAACCCAGCCTCCTTGAATCCTCCCAAGTCTCAGAGAGGACTCCTTTGCCCTGTGTT
CCTCCAGGCCAGATGGACTACCGAAGATTGCAGCCTTCTTGCCTGGGACCATGCCCTGTCTGTGTG
CCCACCATGGCTGAGTCAGGGTCTGCTGTACCACCCACATTGCCAACCATTCTACTACCTACCACTAAGC
TATTGGCAGCAGCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC400171 representing NM_005373
 Red=Cloning site Green=Tags(s)

MPSWALFMVTSCLLLAPQNLAQVSSQDVSLASDSEPLKCFSRTEFDLTCFWDEEEAAPSGTYQLLYAYP
 REKPRACPLSSQSMPHFGTRYVCQFPDQEEVRLFFPLHLWVKNVFLNQTRTQRVLFVDSVGLPAPPSIIK
 AMGGSQPGELQISWEEPAPPEISDFLRYELRYGPRDPKNSTGPTVIQLIATETCCPALQRPHSASALDQSP
 CAQPTMPWQDGPQKQTSREASALTAEGGSCLISGLQPGNSYWLQLRSEPDGISLGGSWGWSLPTVDL
 PGDAVALGLQCFITLTLKNTVCWQQQDHASSQGFYHSRARCPRDRYPWENECHEEETNPGLQTPQFS
 RCHFCSRNDISIIHILVEVTTAPGTVHSYLGSFVIHQAVRLPTPNLHWREISSGHLELEWQHPSSWAAQE
 TCYQLRYTGEHQDWKVLPEPLGARGGTLELRPRSRYRLQLRRLNGPTYQGPWSSWSDPTRVETATETA
 WISLVTAHLVLGLNAVLGLLLRWQFPAHYRRLRHAWPSLPDLHRVLGQYLRDTAALSPPKATVSDTC
 EEVEPSLLEILPKSSERTPLPLCSSQAQMDYRRLQPSCLGTMPLSVCPMAESGSCCTTHIANHSYLPLS
 YWQQP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
RefSeq:	NP_005364
RefSeq Size:	3645 bp
RefSeq ORF:	1908 bp
Locus ID:	4352
Cytogenetics:	1p34.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway
MW:	71 kDa
Gene Summary:	<p>In 1990 an oncogene, v-mpl, was identified from the murine myeloproliferative leukemia virus that was capable of immortalizing bone marrow hematopoietic cells from different lineages. In 1992 the human homologue, named, c-mpl, was cloned. Sequence data revealed that c-mpl encoded a protein that was homologous with members of the hematopoietic receptor superfamily. Presence of anti-sense oligodeoxynucleotides of c-mpl inhibited megakaryocyte colony formation. The ligand for c-mpl, thrombopoietin, was cloned in 1994. Thrombopoietin was shown to be the major regulator of megakaryocytopoiesis and platelet formation. The protein encoded by the c-mpl gene, CD110, is a 635 amino acid transmembrane domain, with two extracellular cytokine receptor domains and two intracellular cytokine receptor box motifs. TPO-R deficient mice were severely thrombocytopenic, emphasizing the important role of CD110 and thrombopoietin in megakaryocyte and platelet formation. Upon binding of thrombopoietin CD110 is dimerized and the JAK family of non-receptor tyrosine kinases, as well as the STAT family, the MAPK family, the adaptor protein Shc and the receptors themselves become tyrosine phosphorylated. [provided by RefSeq, Jul 2008]</p>