

Product datasheet for **RC400235**

beta Catenin (CTNNB1) (NM_001904) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	beta Catenin (CTNNB1) (NM_001904) Human Mutant ORF Clone
Mutation Description:	A5_A80>D
Affected Codon#:	5
Affected NT#:	c.14_241
Nucleotide Mutation:	CTNNB1 Mutant (A5_A80>D), Myc-DDK-tagged ORF clone of Homo sapiens catenin (cadherin-associated protein), beta 1, 88kDa (CTNNB1), transcript variant 1 as transfection-ready DNA
Effect:	deletion
Symbol:	beta Catenin
Synonyms:	armadillo; CTNNB; EVR7; MRD19; NEDSDV
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_001904
ORF Size:	2115 bp
Restriction Sites:	Sgfl-MluI

[View online »](#)

**ORF Nucleotide
Sequence:**

>RC400235 representing NM_001904
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTACTCAAGATATTGATGGACAGTATGCAATGACTCGAGCTCAGAGGGTACGAGCTGCTATGTTCC
 CTGAGACATTAGATGAGGGCATGCAGATCCCATCTACACAGTTTGATGCTGCTCATCCCACTAATGTCCA
 GCGTTTGGCTGAACCATCACAGATGCTGAAACATGCAGTTGTAAACTTGATTAACATCAAGATGATGCA
 GAACTTGCCACACGTGCAATCCCTGAAGTACAAAACTGCTAAATGACGAGGACCAGGTGGTGGTTAATA
 AGGCTGCAGTTATGGTCCATCAGCTTTCTAAAAAGGAAGCTTCCAGACACGCTATCATGCGTTCTCCTCA
 GATGGTGTCTGCTATTGTACGTACCATGCAGAATACAAATGATGTAGAAACAGCTCGTTGTACCGCTGGG
 ACCTTGCATAACCTTTCCCATCATCGTGAGGGCTTACTGGCCATCTTTAAGTCTGGAGGCATTCTGCCC
 TGGTGAAGATGCTTGGTTCACCACTGGATTCTGTGTTGTTTTATGCCATTACAACTCTCCACAACCTTTT
 ATTACATCAAGAAGGAGCTAAAAATGGCAGTGCCTTTAGCTGGTGGGCTGCAGAAAATGGTTGCCTTGCTC
 AACAAAACAAAATGTTAAATCTTGGCTATTACGACAGACTGCCTTCAAATTTTAGCTTATGGCAACCAAG
 AAAGCAAGCTCATCACTACTGGCTAGTGGTGGACCCCAAGCTTTAGTAAATATAATGAGGACCTATACTTA
 CGAAAACTACTGTGGACCACAAGCAGAGTGTGAAGGTGCTATCTGTCTGCTCTAGTAATAAGCCGGCT
 ATTGTAGAAGCTGGTGAATGCAAGCTTTAGGACTTCACCTGACAGATCCAAGTCAACGCTTGTTCAGA
 ACTGTCTTTGGACTCTCAGGAATCTTTGAGATGCTGCACTAAACAGGAAGGGATGGAAGGTCTCCTTGG
 GACTCTTGTTCAGCTTCTGGGTTGAGATGATATAAATGTGGTCACCTGTGCAGCTGGAATCTTTCTAAC
 CTCCTTGAATAATTATAAGAACAAGATGATGGTCTGCCAAGTGGGTGGTATAGAGGCTCTTGTGCGTA
 CTGCTCTCGGGCTGGTGACAGGAAGACATCACTGAGCCTGCCATCTGTGCTCTTCGTATCTGACCA
 CCGACACCAAGAAGCAGAGATGGCCAGAATGCAGTTTCGCTTCACTATGGACTACCAAGTTGTGGTTAAG
 CTCTTACCCACCATCCCACTGGCCTCTGATAAAGGCTACTGTTGGATTGATTCGAAATCTTGCCCTTT
 GTCCCGCAAAATCATGCACCTTTGCGTGAGCAGGGTGCCATTCCACGACTAGTTCAGTTGCTTGTTCGTGC
 ACATCAGGATACCCAGCGCCGTACGTCCATGGGTGGGACACAGCAGCAATTTGTGGAGGGGGTCCGCATG
 GAAGAAATAGTTGAAGGTTGTACCGGAGCCCTTACATCCTAGCTCGGGATGTTCAACCGAATTGTTA
 TCAGAGGACTAAATACCATTCCATTGTTTGTGAGCTGCTTTATTCTCCATTGAAAACATCCAAAGAGT
 AGCTGCAGGGGTCTCTGTGAATTTGCTCAGGACAAGGAAGCTGCAGAAGCTATTGAAGCTGAGGGAGCC
 ACAGCTCCTCTGACAGAGTTACTTCACTCTAGGAATGAAGGTGTGGCGACATATGCAGCTGCTGTTTTGT
 TCCGAATGTCTGAGGACAAGCCACAAGATTACAAGAAACGGCTTTCAGTTGAGCTGACCAGCTCTCTCTT
 CAGAACAGAGCCAATGGCTTGAATGAGACTGCTGATCTTGGACTTGATATTGGTGCCCGAGGAGAACCC
 CTTGGATATCGCCAGGATGATCCTAGCTATCGTTCTTTTCACTCTGGTGGATATGGCCAGGATGCCTTGG
 GTATGGACCCCATGATGGAACATGAGATGGGTGGCCACCACCTGGTGTGCTGACTATCCAGTTGATGGGCT
 GCCAGATCTGGGCGATGCCAGGACCTCATGGATGGGCTGCCTCCAGGTGACAGCAATCAGCTGGCCTGG
 TTTGATACTGACCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC400235 representing NM_001904
Red=Cloning site Green=Tags(s)

MATQDIDGQYAMTRAQVRVRAAMFPETLDEGMQIPSTQFDAAHPTNVQRLAEPQMLKHAVVNLINYQDDA
ELATRAIPELTKLLNDEDQVVVNKAAMVHQLSKKEASRHAIMRSPQMVSIVRTMQNTNDVETARCTAG
TLHNLSSHREGLLAIFKSGGIPALVKMLGSPVDSVLFYAITTLHNLHLLHQEGAKMAVRLAGGLQKMVALL
NKTNVKFLAITTDCLQILAYGNQESKLIILASGGPQALVNIMRTYTYEKLWTTSRVLKVL SVCSSNKPA
IVEAGGMQALGLHLTDPSQRLVQNCWTLRNLSDAATKQEGMEGLLGTLVQLLGSDDINVVTCAAGILSN
LTCNNYKNKMMVCQVGGIEALVRTVLRAGDREDITEPAICALRHLLSRHQEAEMAQNAVRLHYGLPVVVK
LLHPPSHWPLIKATVGLIRNLALCPANHAPLREQGAIPRLVQLLVRAHQDTQRRTSMGGTQQQFVEGV
RMEEIVEGCTGALHILARDVHNRI VIRGLNTIPLFVQLLYSPIENIQRVAAGVLC ELAQDKEAAEAIEAEGA
TAPL TELLHSRNEG VATYAAAVLFRMSDKPQDYKKRLSVELTSSLFRTEPMAWNETADLGLDIGAQGE
P LGYRQDDPSYRSFHSGGYGQDALGMDPMMEHEMGGHHPGADYPVDGLPDLGHAQDLMDGLPPGDSNQLAW
FDTDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NP_001895
RefSeq Size:	3720 bp
RefSeq ORF:	2346 bp
Locus ID:	1499
Cytogenetics:	3p22.1
Domains:	Armadillo_seg
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors
Protein Pathways:	Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Basal cell carcinoma, Colorectal cancer, Endometrial cancer, Focal adhesion, Leukocyte transendothelial migration, Melanogenesis, Pathogenic Escherichia coli infection, Pathways in cancer, Prostate cancer, Thyroid cancer, Tight junction, Wnt signaling pathway
MW:	77 kDa

Gene Summary:

The protein encoded by this gene is part of a complex of proteins that constitute adherens junctions (AJs). AJs are necessary for the creation and maintenance of epithelial cell layers by regulating cell growth and adhesion between cells. The encoded protein also anchors the actin cytoskeleton and may be responsible for transmitting the contact inhibition signal that causes cells to stop dividing once the epithelial sheet is complete. Finally, this protein binds to the product of the APC gene, which is mutated in adenomatous polyposis of the colon. Mutations in this gene are a cause of colorectal cancer (CRC), pilomatrixoma (PTR), medulloblastoma (MDB), and ovarian cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2016]