

Product datasheet for **RC401613**

MLH1 (NM_000249) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	MLH1 (NM_000249) Human Mutant ORF Clone
Mutation Description:	E512X
Affected Codon#:	512
Affected NT#:	1534
Nucleotide Mutation:	MLH1 Mutant (E512X), Myc-DDK-tagged ORF clone of Homo sapiens mutL homolog 1, colon cancer, nonpolyposis type 2 (E. coli) (MLH1), transcript variant 1 as transfection-ready DNA
Effect:	Colorectal cancer, non-polyposis
Symbol:	MLH1
Synonyms:	COCA2; FCC2; hMLH1; HNPCC; HNPCC2; MMRCS1
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_000249
ORF Size:	1533 bp
Restriction Sites:	SgfI-MluI



[View online »](#)

ORF Nucleotide Sequence:

>RC401613 representing NM_000249
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCGTTCTGGCAGGGGTTATTCGGCGGCTGGACGAGACAGTGGTGAACCGCATCGCGCGGGGGAAG
 TTATCCAGCGGCCAGCTAATGCTATCAAAGAGATGATTGAGAAGTGTGTTAGATGCAAAATCCACAAGTAT
 TCAAGTGATTGTTAAAGAGGGAGGCCTGAAGTTGATTCAGATCCAAGACAATGGCACCGGGATCAGGAAA
 GAAGATCTGGATATTGTATGTGAAAGGTTCACTACTAGTAAACTGCAGTCCTTTGAGGATTTAGCCAGTA
 TTTCTACCTATGGCTTCGAGGTGAGGCTTTGGCCAGCATAAGCCATGTGGCTCATGTTACTATTACAAC
 GAAAACAGCTGATGAAAGTGTGCATACAGAGCAAGTACTCAGATGGAAAAGTAAAGCCCTCTCTAAA
 CCATGTGCTGGCAATCAAGGGACCCAGATCACGGTGGAGGACCTTTTTTACAACATAGCCACGAGGAGAA
 AAGCTTTAAAAATCCAAGTGAAGAATATGGGAAAATTTGGAAGTTGTTGGCAGGTATTCAGTACACAA
 TGCAGGCATTAGTTTCTAGTTAAAAACAAGGAGAGACAGTAGCTGATGTTAGGACACTACCCAATGCC
 TCAACCGTGGACAATATTCGCTCCATCTTTGAAATGCTGTTAGTCGAGAAGTATAGAAATGGATGTG
 AGGATAAAACCTAGCCTTCAAATGAATGGTTACATATCCAATGCAAACACTCAGTGAAGAAGTGCAT
 CTTCTTACTCTTCATCAACCATCGTCTGGTAGAATCAACTTCTTGAGAAAAGCCATAGAAAACAGTGTAT
 GCAGCCTATTTGCCAAAAACACACACCATTCTGTACCTCAGTTTAGAAAATCAGTCCCAGAAATGTGG
 ATGTTAATGTGCACCCCAAAAGCATGAAGTTCACCTCCTGCACGAGGAGAGCATCCTGGAGCGGGTGCA
 GCAGCAGATCGAGAGCAAGCTCCTGGGCTCAAATTCCTCCAGGATGTACTTACCCAGACTTTGCTACCA
 GGACTTGCTGGCCCTCTGGGGAGATGGTTAAATCCACAACAAGTCTGACCTCGTCTTCTACTTCTGGAA
 GTAGTGATAAGGTCTATGCCACCAGATGGTTCGTACAGATTCCCGGGAACAGAAGCTTGATGCATTTCT
 GCAGCCTCTGAGCAAACCCCTGTCCAGTCAGCCCGAGCCATTGTACAGAGGATAAGACAGATATTTCT
 AGTGGCAGGGCTAGGCAGCAAGATGAGGAGATGCTTGAACCTCCAGCCCTGCTGAAGTGCTGCCAAAA
 ATCAGAGCTTGGAGGGGATACAACAAAGGGGACTTCAAAAATGTGAGAGAAGAGAGGACCTACTTCCAG
 CAACCCAGAAAGAGACATCGGGAAGATTCTGATGTGGAATGGTGAAGATGATTTCCGAAAAGGAAATG
 ACTGCAGCTTGTACCCCGGAGAAGGATCATTAACCTCACTAGTGTGTTTGTGCTCCAGGAA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:

>RC401613 representing NM_000249
 Red=Cloning site Green=Tags(s)

MSFVAGVIRRLDETVVNRIAAGEVIQRPANAIEKIENCLDAKSTSIQVIVKEGGLKLIQIQDNGTGIRK
 EDLDIVCERFTTSKLFQFEDLASISTYGRGEALASISHVAHVITTKTADGKAYRASYSKGKLPKPPK
 PCAGNQGTQITVEDLFYNIATRRLKKNPSEYKILEVVGRYSVHNAGISFSVKKQGETVADVRLPNA
 STVDNIRSIIFGNVSRLEIEIGCEDKTLAFKMNGYISNANYSVKKCIFLLFINHRLVESTSLRKAIEVY
 AAYLPKNTHPFLYLSLEISPNVDVNVHPTKHEVHFLHEESILERVQQHIESKLLGNSSSRMFTQTLPL
 GLAGPSGEMVKSTTSLTSSSTSGSSDKVYAHQMVRTDSREQLDAFLQPLSKPLSSQPQAIPTEDKTDIS
 SGRARQQDEMLELPAPAEVAANKQSLEGDTTKGTSEMSEKRGPTSSNPKRHRESDVEMVEDDSRKEM
 TAACTPRRRIINLTSVLSLQE

SGPTRTRRLEQKLI**SEEDLAANDILDYKDDDDKV**

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_000240](#)

RefSeq Size:

1533 bp

RefSeq ORF:

2271 bp

Locus ID:

4292

Cytogenetics:

3p22.2

Domains:

DNA_mis_repair, HATPase_c

Protein Families:

Druggable Genome

Protein Pathways:

Colorectal cancer, Endometrial cancer, Mismatch repair, Pathways in cancer

MW:

56.2 kDa

Gene Summary:

The protein encoded by this gene can heterodimerize with mismatch repair endonuclease PMS2 to form MutL alpha, part of the DNA mismatch repair system. When MutL alpha is bound by MutS beta and some accessory proteins, the PMS2 subunit of MutL alpha introduces a single-strand break near DNA mismatches, providing an entry point for exonuclease degradation. The encoded protein is also involved in DNA damage signaling and can heterodimerize with DNA mismatch repair protein MLH3 to form MutL gamma, which is involved in meiosis. This gene was identified as a locus frequently mutated in hereditary nonpolyposis colon cancer (HNPCC). [provided by RefSeq, Aug 2017]