

Product datasheet for **RC401614**

MLH1 (NM_000249) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	MLH1 (NM_000249) Human Mutant ORF Clone
Mutation Description:	G517X
Affected Codon#:	517
Affected NT#:	1549
Nucleotide Mutation:	MLH1 Mutant (G517X), Myc-DDK-tagged ORF clone of Homo sapiens mutL homolog 1, colon cancer, nonpolyposis type 2 (E. coli) (MLH1), transcript variant 1 as transfection-ready DNA
Effect:	Colorectal cancer, non-polyposis
Symbol:	MLH1
Synonyms:	COCA2; FCC2; hMLH1; HNPCC; HNPCC2; MMRCS1
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_000249
ORF Size:	1548 bp
Restriction Sites:	SgfI-MluI

[View online »](#)

ORF Nucleotide Sequence:

>RC401614 representing NM_000249
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGTTCTGTCGAGGGGTTATTCGGCGGCTGGACGAGACAGTGGTGAACCGCATCGCGCGGGGGAAG
 TTATCCAGCGGCCAGCTAATGCTATCAAGAGAGTATTGAGAAGTGTGTTAGATGCAAAATCCACAAGTAT
 TCAAGTATTGTTAAAGAGGGAGGCCTGAAGTTGATTGATCCAAAGACAATGGCACCGGGATCAGGAAA
 GAAGATCTGGATATTGTATGTGAAAGGTTCACTACTAGTAAACTGCAGTCCTTTGAGGATTTAGCCAGTA
 TTTCTACCTATGGCTTTCGAGGTGAGGCTTTGGCCAGCATAAGCCATGTGGCTCATGTTACTATTACAAC
 GAAAACAGCTGATGAAAGTGTGCATACAGAGCAAGTTACTCAGATGGAAAAGTAAAGCCCTCCTAAA
 CCATGTGCTGGCAATCAAGGGACCCAGATCACGGTGGAGGACCTTTTTTACAACATAGCCACGAGGAGAA
 AAGCTTTAAAAAATCAAGTGAAGAATATGGGAAAATTTTGAAGTTGTTGGCAGGTATTCAGTACACAA
 TGCAGGCATTAGTTTCTAGTTAAAAACAAGGAGAGACAGTAGCTGATGTTAGGACACTACCCAATGCC
 TCAACCGTGGACAATATTCGCTCCATCTTTGAAATGCTGTTAGTCGAGAACTGATAGAAATTGGATGTG
 AGGATAAAACCTAGCCTTCAAATGAATGGTTACATATCCAATGCAAACTACTCAGTGAAGAAGTGCAT
 CTTCTTACTCTTCATCAACCATCGTCTGGTAGAATCAACTTCCTTGAGAAAAGCCATAGAAAACAGTGTAT
 GCAGCCTATTTGCCCAAAACACACACCCATTCTGTACCTCAGTTTAGAAATCAGTCCCCAGAATGTGG
 ATGTTAATGTGCACCCCAAAAGCATGAAGTTCACCTCCTGCACGAGGAGAGCATCCTGGAGCGGGTGCA
 GCAGCACATCGAGAGCAAGCTCCTGGGCTCCAATTCTCCAGGATGTACTTCACCCAGACTTTGCTACCA
 GGACTTGCTGGCCCTCTGGGGAGATGGTTAAATCCACAACAAGTCTGACCTCGTCTTCTACTTCTGGAA
 GTAGTGATAAGGTCTATGCCACCCAGATGGTTCGTACAGATTCCCGGGAACAGAAGCTTGATGCATTTCT
 GCAGCCTCTGAGCAAACCCCTGTCCAGTCAGCCCCAGGCCATTGTCACAGAGGATAAGACAGATATTTCT
 AGTGGCAGGGCTAGGCAGCAAGATGAGGAGATGCTTGAATCCCAGCCCCTGCTGAAGTGCTGCCAAAA
 ATCAGAGCTTGGAGGGGGATACAACAAGGGGACTTCAGAAATGTCAGAGAAGAGAGGACCTACTTCCAG
 CAACCCAGAAAGAGACATCGGGAAGATTCTGATGTGGAATGGTGAAGATGATTCCCGAAAGGAAATG
 ACTGCAGCTTGTAACCCCGGAGAAGGATCATTAACCTCACTAGTGTGTTTGTAGTCTCCAGGAAGAAATTA
 ATGAGCAG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:


>RC401614 representing NM_000249
 Red=Cloning site Green=Tags(s)

MSFVAGVIRRLDETVVNRIAAGEVIQRPANAIEKMIENCLDAKSTSIQVIVKEGGLKLIQIQDNGTGIRK
 EDLDIVCERFTTSKLQSFEDLASISTYGRGEALASISHVAHVITTKTADGKCAYRASYSYGKLGKAPPK
 PCAGNQGTQITVEDLFYNIATRRKALKNPSEEGKILEVVGRYSVHNAGISFSVKKQGETVADVRLPNA
 STVDNIRSIIFGNAVSRELIEIGCEDKTLAFKMNGYISNANYSVKKCIFLLFINHRLVESTSLRKAETVY
 AAYLPKNTHPFLYLLEISPNQVDVNVHPTKHEVHFLHEESILERVQQHIESKLLGSNSSRMFTQTLPL
 GLAGPSGEMVKSTSLTSSSTSGSSDKVYAHQMVRTDSREQKLDAFLQPLSKPLSSQPQAIVTEDKTDIS
 SGRARQQDEMLELPAPAEVAAKNQSLGDTTKGTSEMSEKRGPTSSNPRKRHRESDVEMVEDDSRKEM
 TAACTPRRRIINLTSVLSLQEEINEQ

SGPTRRRLE**QKLI**SEED**LA**AND**ILDYKDDDDKV**

Restriction Sites:

SgfI-MluI

Cloning Scheme:

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_000240](#)

RefSeq Size:

1548 bp

RefSeq ORF:

2271 bp

Locus ID:

4292

Cytogenetics:

3p22.2

Domains:

DNA_mis_repair, HATPase_c

Protein Families:

Druggable Genome

Protein Pathways:

Colorectal cancer, Endometrial cancer, Mismatch repair, Pathways in cancer

MW:

56.8 kDa

Gene Summary:

The protein encoded by this gene can heterodimerize with mismatch repair endonuclease PMS2 to form MutL alpha, part of the DNA mismatch repair system. When MutL alpha is bound by MutS beta and some accessory proteins, the PMS2 subunit of MutL alpha introduces a single-strand break near DNA mismatches, providing an entry point for exonuclease degradation. The encoded protein is also involved in DNA damage signaling and can heterodimerize with DNA mismatch repair protein MLH3 to form MutL gamma, which is involved in meiosis. This gene was identified as a locus frequently mutated in hereditary nonpolyposis colon cancer (HNPCC). [provided by RefSeq, Aug 2017]